


# AXIOM improv

INDOOR / OUTDOOR MOBILITY PUSHCHAIR



www.adaptivestar.com 1-800-278-9626



# Family First We believe everyone using our products becomes family

Our Mission is to fulfill the dreams of families with products and services that benefit their daily lives. Our commitment grows with every member that joins our family, to listen, to learn and to improve our service, our innovation and our products.

Family is LOVE
Family is TRUST
Family is RESPECT
Family is ACCEPTANCE

Our Family grows with every phone call, every shipment, every giddy success and discouraging failure.

Thank you for joining our FAMILY





# Frequently Asked Questions

## Can I roller blade with my Axiom Push Chair?

We do not recommend this, it is just too easy to lose control and cause possible harm to your passenger or yourself

## What if I break my Axiom Pushchair?

No problem! Just call us for help at 800-278-9626. We'll get you back up and on the road in no time.

## Is the fabric on my Axiom Pushchair washable?

Yes! The seat, canopy and basket are removable and machine washable in cold water with mild detergent. After washing, place the seat back onto the Pushchair for drying. Machine drying in not recommended. Remember to remove the seat panels before washing and to replace them when you're finished.

## Reminders

ALWAYS use your safety wrist strap any time a passenger is in the Pushchair.

ALWAYS engage rear food brake when loading and unloading a passenger.


NEVER leave a passenger unattended in or around the Pushchair.

NEVER overload the Pushchair. Overloading the Pushchair may cause it to become unstable and tip over.

DO NOT allow your passenger to stand in the Pushchair. This may cause the Pushchair to tip and cause injury to the passenger, yourself or others.



Your AXIOM Pushchair is almost completely assembled when you receive it. The following short series of instructions will help you through the final few steps so you can get out and start having fun!



- A Axis swivel wheel kit (boxed)
- B AXIOM pushchair (folded)
- C Wheel set 2 Rear 1 Front
- D Rain Canopy
- E Assembly instructions and Pump

Remove all the contents of the box and confirm that they match the above list. If you determine that you are missing something, please call 1-800-278-9626 and talk with one or our helpful technicians.

First unfold your Axiom and place it into the upright position as shown here by pulling up on the handle until it locks in place.

## Lets get started!



From this position...



To upright and locked.



Make sure the units is unfolded and lock in place before loading the passenger and using the pushchair.



## Sun Canopy Installation





The sun canopy is preinstalled for your convenience. The following instructions outline how to reinstall the sun canopy. The next step is to install the sun canopy onto your Axiom Pushchair. Orient the canopy so that the Velcro flap is toward the push handle and the mesh visor is toward the front of the pushchair. Slip the plastic canopy clip through the fabric openings on each side of the seat and while applying a slight twisting pressure onto the canopy clip snap them individually onto the pushchair frame.





Once you have the canopy clips installed, connect the Velcro flap to the on the canopy to the seat back and pull the mesh front out until it looks like the photo above.



## Rear Wheel Installation

First unlock the rear foot brake by pulling the black loop upward





Next, push the black button and while holding it in... ...insert the wheel axle all the way into the pushchair frame and release the button.

Now repeat to install the other wheel.





Check to make sure your wheels are locked in before each use.

#### Front Wheel Installation

Start by tipping the push chair backward so the fork assembly is pointing upward. This will make it easier to install the wheel.



The front wheel is the same as a bicycles. First, find the brake caliper release (see pg 11) and open it. Place the tire between the brake pads and orient your wheel so that the flat sections of the wheel axle aligns with the slotted metal

fork tips and push the wheel so that it seats all the way in the back of the fork tip.



Position the silver lever straight out from the frame as shown. Now tighten the black knob on



the opposite side of the lever until it just starts to become snug.

Be sure to dispose of the packaging from the wheels, they can be a choking hazard to young children.

Now swing the lever closed as shown here. If you

are unable to swing the lever completely into the locked position, open it back straight out, loosen the black knob



slightly, and close the lever again.



## Axis Swivel Wheel Adaptor

The Axis swivel wheel adaptor kit comes

standard with each Axiom Pushchair. In order to install the swivel wheels first you must remove the



standard front wheel by unfolding the quick release lever and sliding the wheel out of the fork tips. Next loosen the black knob on the top of the Axis swivel wheel adaptor and unfold the hinged lid.



Now slide both forks onto the Axis swivel wheel adaptor until the fork tips slide over the round cross bar.



Marin Ing

Check to make sure your Axis Swivel Wheel Adaptor is properly installed before each use. Finally, fold the hinge flap closed and turn the black knob until snug.

Give the Axis swivel wheel kit a tug and make sure the wheels are on tightly.





## Folding your AXIOM Improv...

Locate the Primary folding lockout release lever located on the left side upright of the AXIOM.

Primary Release Latch Lever



Squeeze and hold the lever, while pushing down slightly on the AXIOM handle. Once the AXIOM has started to fold, you may release the lever.



Your AXIOM improv folds to a convenient size that will allow you to transport it in many automobile trunks. If you require a more compact folded size, simply remove the wheels.

When restoring your AXIOM improve to its upright and locked use position, be certain that there are no obstructions where the latch mechanism comes together.

Make sure to confirm you have the Axiom Improv locked into the upright position by pushing down on the handle





## AXIOM Improv Front wheel splash guard, installation guide

Start by locating the front brake caliper retaining nut, located behind the front wheel. Use a wrench to remove the top nut as shown:



Now set the removed nut aside while you slip the splash guard onto the caliper mounting bolt. Hold the splash guard in place and reinstall the mounting nut.

Tighten snuggly.



Please note this is a splash guard only and is not to be used as a foot rest



# Things you need to know about....



## All Weather Canopy

If you get stuck in the rain or heavy wind, the all weather canopy will help keep the pushchair passenger dry and comfortable. Just drape the cover over your bubble style sun canopy and seal the back of the Velcro strips together around the handle like a sleeve. Next, buckle the straps around the foot plate, which attach underneath between the brake caliper and the pushchair frame.

## Adjustable Height Handle

Your pushchair is equipped with a height adjustable handle. To adjust the handle to best suit your comfort, simply press the buttons located on the inside of the handle near the end of the foam handle grip. While holding the buttons in, adjust the handle up or down to the height that is correct for you. Once there, release the buttons and the handle will lock into place.

Braking systems



Your pushchair is equipped with a brake for the front wheel and a brake for the rear wheels.

The front wheel brake ONLY works with the standard front wheel and



WILL NOT

work with the Axis swivel wheel adaptor. pressed position when it is locked.

The rear wheel brake works by pushing down on the black bar with the loop that runs between the rear wheels. It will lock into the black toothed hubs on the rear wheels.



Always use the rear footbrake when loading and unloading your passenger. Always use the rear footbrake anytime a passenger is in the pushchair and you are not pushing or within reach of the pushchair.



# More things you need to know about....

## Seat Anchor Straps

The seat harness system requires that the seat anchor straps be attached at all times. These straps are removable so that you can take off the seat for laundering etc. Please add this to your visual inspection prior to using the pushchair.

## Wrist Strap


Your Axiom pushchair is equipped with a safety wrist strap which you should always use when there is a passenger loaded. The wrist strap is the link between you and the pushchair and keeps the pushchair from rolling away from your control by accident. Never leave a passenger unattended in the pushchair. The wrist strap is constructed with loops on each end, one large and one small. The wrist strap should be attached to the pushchair at the frame member that runs between the rear wheels, by looping around it, and inserting the large looped end through the small loop and pulling it tight. The large loop may now be wrapped around your wrist attaching you to the pushchair.





Never leave a passenger in the pushchair unattended

Always use the wrist strap anytime a passenger is in the pushchair



# Adjustments

#### Seat Recline

To recline the seat, simply slide the two straps behind the seat. To put the seat back to the upright position, pull the webbing taut until in the upright position.



## Front Brake Adjustment

The front brake on your Axiom pushchair may require periodic adjustment. By turning the silver ring, indicated by the red arrow, outward it will cause the brake pads to move closer to the front wheel. There is a similar adjuster located underneath the pushchair, near the front wheel, on the brake caliper itself. This caliper style brake system is the same as used on a bicycle, and if you require more detailed assistance, a local bicycle shop will be able to provide it for you.



## Tracking

If your pushchair is pulling to the left or right while using the fix front wheel, check to see if it is one of the following:

Tire pressure, check to see if all 3 tires are fully inflated.

Be sure that the front wheel hasn't been installed **too** tight. Release the QR lever, push the front wheel all the way back into the forktips and when closing the lever, make sure that it is only snug, not *He-man* tight. If it is, loosen the thumb nut on the other side and close the lever. You may have to make several small adjustments until the lever closes snuggly.



Inspect the wheels for bent spokes or axles.



## Harness Adjustment

You will want to do a one time general positioning of the harness before placing the passenger in the seat. Once the passenger is in the seat then do the final harness adjustments as follows:

Undo both the buckles from the padded center strap.



Push down on the silver center button located between the passengers knees and pull the webbing out to the end.



Pull up the webbing through the seat and up on the shoulder harness area.

Place the passenger in the seat.

Adjust the padded center strap forward or backward to fit your passenger.

Adjust the padded shoulder harness up or down on the Velcro strips to fit your passenger.

Buckle the two buckles to the padded center strap.

Pull the adjustment strap until the padded center strap and the padded shoulder harness are both snug. Adjust the chest strap to a comfortable length and connect the buckles.

If the harness is too small for the passenger harness extensions are available for purchase.

<u>Leminder</u>

Always check that the seat anchor straps are attached properly before each use.



Avoid serious injury from falling or sliding out; Always us the seat harness system



## Maintenance

We recommend that you have an experienced dealer perform these maintenance check at least once each year. Lubrication of these following items is necessary:

- · Pivot point and cable of hand brake
- Pivot point and cable of brake caliper
- Front wheel axle
- Front wheel quick release lever

- Insert in frame for rear wheels
- · Wheel hub and bearings of rear wheels
- Locking bars of folding mechanism
- All folding joints of pushchair

Check to make sure all of the bolts are tight and that springs are secure on the rear foot brake.

Check fabric seams and straps for fraying and tearing.

Keep pushchair our of the sun and in a dry place when not in use, to prevent mold and mildew.

Check the frame for loose rivets, cracks and wear.

Check to make sure all nuts and bolts are tight on the frame.

We recommend that you put a rust preventative type coating on your wheels, nuts, bolts, brake springs and brake caliper on your pushchair. These parts can also be maintained by cleaning with a damp clothe and keeping these parts dry and free from debris.

The tires on your pushchair are a standard bicycle size and have an inner tube. In the event that one of your tires will not hold air, the inner tube can be replaced and both tires and inner tubes can be found at professional bicycle stores and most major retailers.

The front brake uses a standard bicycle style caliper and brake pad. If replacement brake pads are needed, they can be found at professional bicycle stores and most major retailers.



#### **WARRANTY CARD**

Adaptive Star, LLC

Po Box 369

Thank you for purchasing an Adaptive Star Indoor/Outdoor Mobility Push Chair. We have created a warranty that will keep your unit in top condition and headache free. The following warranty is on workmanship and materials from the date of **original consumer** purchase:

- A 5-year warranty on the frame.
- A 1-year warranty on parts and soft goods. Exclusions include:
- Parts that require replacement under normal use.
- Damage or malfunction caused by negligence or abuse.
- Foam grips, inner tubes, tires and brake pads are specifically excluded.

Adaptive Star will, at its option, repair or replace defective parts upon delivery to our service center accompanied by proof of date of first consumer purchase, such as a copy of a sales receipt. Please refer to the store where you made your purchase for returns and refunds as Adaptive Star does not give direct refunds. Purchaser is responsible for initial shipping charges required to ship the product for warranty service. Purchaser is also responsible for all applicable duties and/or custom charges. The return charges will be at Adaptive Star's expense. Please include your name, address, telephone number, email address and a brief description of the problem, Adaptive Star does not accept delivery of any product without a Return Authorization # clearly marked on the return shipment box and shipping label.

All implied warranties, including warranties of merchantability and fitness for a particular purpose are limited in duration to the length of this warranty. Adaptive Star will not be liable for any incidental, consequential, or other damages, including, without limitation, damages resulting from loss of use. Some states do not allow limitations on how long an implied warranty last and/or do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation may not apply to you. This limited warranty gives you specific legal rights, and you may also have other rights, which vary from state to state.

Your complete satisfaction is our goal. Please let us know if you are dissatisfied with any of our products. We will do our very best to correct the situation and make you happy. Feel free to contact us and send photos. We love to hear from you!

Selah, WA 98902 cus Telephone: 800-	stomerservice@adaptivestar.com			
FAX: 509-	453-8326			
Name				
Address				
	State	Zip	Phone	
E-Mail				
Model	odelDate of Purchase			
Where did you pu	ırchase your unit?			
Ο,	gestions for improvements to p			
	, ,		send you some information?	

