

atHand™

OVERBED TABLE SYSTEM

ASSEMBLY INSTRUCTIONS:

Accessory Kit Adding an Additional Open Shelf Adding an Additional Drawer

Warnings:

Max power rating for power outlets: 125V

Max weight capacity of storage compartments: 5 pounds each

Max weight capacity of tray table: 100 pounds

This product is not intended for use as support or bed rail

For indoor use only

Electrical components are safety agency listed: UL, CSA or equivalent.

Warranty:

At Integrant, we take great pride in the quality and craftsmanship of our products. Attention to design, materials, construction and safety is our priority. We carefully inspect your order prior to shipment.

Upon receipt, please inspect your purchase and notify us of any damage or defect immediately. We will arrange for a prompt replacement.

If, within 30 days of receipt, you are dissatisfied for any reason, you may return your purchase for a full refund of the merchandise value. Simply return the unused item in its original packaging to the address below. You are responsible for return shipping costs.

The product is warranted against defects for three years.

Visit integrantllc.com for details.



INTEGRANT™

Manufactured by:
Integrant LLC
12315 York Delta Drive
North Royalton, OH 44133
440-628-9550
www.integrantllc.com


#16780

ACCESSORY KIT


Kit includes the accessories shown below.



Accessories can be placed in any accessory hole. An example of placement is shown below.



To position the cupholder, open bin or trash bin, select a location for each. Line up the accessory hooks with mounting holes. Push the accessory hooks into the holes, then slide down to lock into place.



Shelf basket fits on the top of the upper shelf or inside either open shelf.

Bags for the trash bin can be secured in place with the bag retainer, as shown below.



ADDING AN ADDITIONAL OPEN SHELF

1. Unplug the ambient light from the USB port.
2. Slide the ambient light towards you to remove from the shelf and set aside.
3. Align the rails on the top of the new open shelf to the rails on the existing assembly and slide into place. Repeat with additional shelves, if necessary.
4. Replace the ambient light by aligning the rails on the top of the light to the rails on the bottom of the shelf.
5. Plug in the USB cord of the ambient light into any USB port.



ADDING AN ADDITIONAL DRAWER

The additional drawer can be placed in any open shelf.

Simply tilt the drawer slightly and slide into place.

Drawer can be easily removed for cleaning.

OPTIONAL ACCESSORIES AVAILABLE SEPARATELY

- An accessory pack includes one each of:
 - a. Cup/Remote holder
 - b. Trash bin with lid
 - c. Shelf basket
 - d. Small bin
- Additional drawer
- Additional open shelf

Visit www.integrantllc.com for ordering details.