



The Aviator Series is the only complete line of aviation-specific transfer chairs available. These chairs meet the Air Carrier Access Act and 14 CFR Part 382 guidelines requiring airlines to provide safe and secure assistance to travelers with disabilities. The broad range of products in the Aviator Series allows the airline to select the chair offering the right blend of features to meet their aircraft, budget and customer requirements.

## AisleMaster™ Boarding Transfer Chair



The AisleMaster is the premier aisle boarding chair in the North American market. The redesign AisleMaster now boasts weight capacity up to 400 lb. With an innovative, compact, mid-wheel layout designed for the tightest turning radius, 'deadman' brake and lifting handles, the AisleMaster combines ease of use with safety and functionality. The narrow model (AWA-8010) offers all the same great features and an extra narrow wheel track for use in smaller commercial/regional jets with narrower aisles.

### FEATURES

- Two wheel track widths: 16" standard and 13" for smaller jets
- Heavy-duty frame construction boasts weight capacity up to 400 lb
- Compact, mid-wheel turning design for the tightest turning radius
- 'Deadman' brake
- Retractable lifting handles
- Flip-up armrests for transfer
- Padded back and seat supports
- Height adjustable padded headrest
- Color-coded crossbody chest belts
- Leg and ankle restraint belts
- Meet Air Carrier Access Act and 14 CFR Part 382

## SkyMaster™ Onboard Wheelchair

### FEATURES

- 300 lb weight capacity lightweight frame
- Compact, folding design for onboard storage
- Front locking wheels
- Rear locking swivel casters allow for easy maneuvering
- Chest and lap safety belts
- ABS seat and back supports meet FAR requirements
- Meet Air Carrier Access Act and 14 CFR Part 382



## ORDER FORM **FAX TO 310-305-1718 OR SUBMIT ONLINE**

Dealer: \_\_\_\_\_ Marked For: \_\_\_\_\_ Account #: \_\_\_\_\_

Date: \_\_\_\_\_ PO#: \_\_\_\_\_ Ordered By: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_ E-mail: \_\_\_\_\_

Ship to address: \_\_\_\_\_

| BASE UNIT                                  | AisleMaster Boarding Chairs |         |                   |       | SkyMaster        |         |
|--|-----------------------------|---------|-------------------|-------|------------------|---------|
|  | AWA-8000                    |         | AWA-8010          |       | AWS-9000         |         |
|  | \$3,548                     |         | \$3,548           |       | \$3,531          |         |
| Wheel Track Width                          | 16"                         |         | 13"               |       | 15.5"            |         |
| Weight Capacity                            | 400 lb                      |         | 400 lb            |       | 300 lb           |         |
| Seat Width                                 | 12.5"                       |         | 12.5"             |       | 12.5"            |         |
| Seat Depth                                 | 15"                         |         | 15"               |       | 14"              |         |
| Seat-to-Floor Height (front)               | 19"                         |         | 19"               |       | 18.5"            |         |
| Back Height                                | 17"                         |         | 17"               |       | 29"              |         |
| Armrest Height                             | 26"                         |         | 26"               |       | 25.5"            |         |
| Push Handle Height                         | 37"                         |         | 37"               |       | 32"              |         |
| Overall Dimensions: L x W x H              | 32" x 16" x 48"             |         | 32" x 16" x 48"   |       | 30" x 16" x 32"  |         |
| Folded Dimensions: L x W x H               | N/A                         |         | N/A               |       | 34" x 17" x 8"   |         |
| Product Weight                             | 40.5 lb (18.4 kg)           |         | 39.5 lb (17.9 kg) |       | 16.7 lb (7.6 kg) |         |
| <b>PARTS</b>                               |                             |         |                   |       |                  |         |
| 8" Main Wheel w/Bearings (pr)              | AWA-R8002                   | \$132   | AWA-R8002         | \$132 |                  |         |
| 3" Front Wheel (pr)                        | AWA-R8003                   | \$66    | AWA-R8003         | \$66  |                  |         |
| Lower Brake Assembly (ea)                  | AWA-R8004                   | \$226   | AWA-R8027         | \$116 |                  |         |
| Upper Brake Grab Bar (ea)                  | AWA-R8005                   | \$127   | AWA-R8005         | \$127 |                  |         |
| Brake Cable & Cable End Lock (ea)          | AWA-R8006                   | \$143   | AWA-R8006         | \$143 |                  |         |
| Brake Springs (pr)                         | AWA-R8007                   | \$39    | AWA-R8007         | \$39  |                  |         |
| Brake Shoe Pad (pr)                        | C0800-80-23                 | \$10    | C0800-80-23       | \$10  |                  |         |
| Seat Cushion                               | AWA-R8008                   | \$358   | AWA-R8008         | \$358 | 9000P-01         | \$284   |
| Back Cushion                               | AWA-R8009                   | \$358   | AWA-R8009         | \$358 | 9000P-02         | \$147   |
| Headrest Complete Kit                      | AWA-R8010                   | \$358   | AWA-R8010         | \$358 |                  |         |
| Rear Flip-Up Handle SPECIFY SIDE (ea)      | AWA-R8012                   | \$182   | AWA-R8012         | \$182 |                  |         |
| Front Telescoping Handle SPECIFY SIDE (ea) | AWA-R8013                   | \$149   | AWA-R8013         | \$149 |                  |         |
| Plastic Grip Handles (pr)                  | AWA-R8014                   | \$50    | AWA-R8014         | \$50  |                  |         |
| Crisscross Body Belt w/Buckle (ea)         | AWA-R8015                   | \$275   | AWA-R8015         | \$275 | 9000P-07         | \$275   |
| Leg Restraint Belt w/Buckle (ea)           | AWA-R8016                   | \$171   | AWA-R8016         | \$171 | 9000P-06         | \$171   |
| Ankle Restraint Belt w/Buckle (ea)         | AWA-R8017                   | \$55    | AWA-R8017         | \$55  |                  |         |
| Armrest, RIGHT SIDE                        | AWA-R8018                   | \$127   | AWA-R8018         | \$127 | 9000P-13         | \$39    |
| Armrest, LEFT SIDE                         | AWA-R8019                   | \$121   | AWA-R8019         | \$121 | 9000P-12         | \$39    |
| Replacement Frame                          | AWA-R8021                   | \$2,266 | N/A               |       | 9000P-17         | \$2,294 |
| Push Handle (ea)                           | AWA-R8022                   | \$154   | AWA-R8022         | \$154 | 9000P-10         | \$165   |
| Footplate (ea)                             | AWA-R8023                   | \$154   | AWA-R8024         | \$154 |                  |         |