

## Nitro DLX Rollator

HCPCS:



FIGURE A




FIGURE B



FIGURE C

- Attractive, Euro-style design
- Handle height easily adjusts with unique push button
- Back support height easily adjusts with tool-free thumb screw
- Caster fork design enhances turning radius
- 10" front casters allow excellent steering and rolling comfort
- Easily folds with one hand to ultra-compact size for storage
- Cross-brace design allows for side-to-side folding and added stability
- Durable and comfortable firm seat separates allowing unit to fold
- Handsome and removable deluxe zippered storage bag with unique attachment to keep bag securely in place when Nitro is open or folded
- Lightweight frame

### PRODUCTS

### Description

RTL10266BL-HS	Nitro,FirmFoam,BlueFrame,RTL 1/ea
RTL10266CH-HS	Nitro,FirmFoam,ChampagneFrame,RTL 1/ea

### SPECIFICATIONS

### Description

Handle Height	33.5" - 38.5"
Seat Dimensions	18" (W) x 10.75" (D)
Seat Height	20"
Unit Dimensions	28" (L) x 24" (W)
Warranty	Frame: Limited Lifetime   Brake Handle and Housing: 5 Years
Weight	22 lbs
Weight Capacity	300 lbs

### ACCESSORIES

### Description

892	Cane Holder, Nitro Rollator, 1/ea, EA
-----	---------------------------------------