

Escape™

THE NEXT STEP IN MOBILITY



Spring-assist Folding
Folds with a single lift of the seat handle.



Large adjustable padded backrest



Anatomical handles and
Cable-free braking system



Easy to adjust handle height



Easy to use frame lock for
safe storage and transport

NEW!

Freedom to go anywhere you want to go...

Escape is available in 3 seat heights to fit almost any user and has a compact folding system that is easy to fold with a simple lift of the release handle. This design makes it easier to lift into the trunk of a car but, as an added benefit, it also stands when folded for easy storage when not in use.

Escape features a large padded seat for added comfort, pressure reducing anatomical handgrips to maintain your hand strength and the easiest to apply hand brakes available on any walker.

In addition, **Escape** offers a depth adjustable large padded backrest, spring loaded levers for easy handle height adjustment and innovative Gentle Touch braking system with no cables to catch on objects or get in your way. **Escape** turns a world of obstacles into a future of possibilities.


Standard Features



- Available in two beautiful colours
- Large adjustable backrest
- Folding shopping bag
- Large padded seat for added comfort
- Easy to apply cable-free brakes
- Easy to use frame lock for safe storage and transport
- Lightweight aluminum construction
- Curb climber
- 8" wheels

Escape™

Enjoy the Freedom!



Optional Accessories

Escape offers a variety of high-quality accessories that help enhance the walker comfort and performance. These can be purchased separately.



Cane Holder (cane not included)
500-4100



Cup Holder
500-4200



Slow Down Brake
500-4000



Plastic Tray
500-4300

Escape Specifications

Colour	Model Numbers		
<ul style="list-style-type: none"> Charcoal Red 	500-10191 500-10195	500-10211 500-10215	500-10241 500-10245
Description	Super Low	Low	Standard
Seat height	19" (48 cm)	21" (54 cm)	24" (61 cm)
Handle height range	28 - 36" (71 - 91.5 cm)	29.5 - 37.5" (75 - 95.3 cm)	32 - 40" (81.3 - 101.6 cm)
Distance between handles	18" (45.7 cm)	18" (45.7 cm)	18" (45.7 cm)
Distance Between wheels	19" (48.25 cm)	19" (48.25 cm)	19" (48.25 cm)
Seat size	18.5 x 9.5" (47 x 24 cm)	18.5 x 9.5" (47 x 24 cm)	18.5 x 9.5" (47 x 24 cm)
Overall width	25" (63.5 cm)	25" (63.5 cm)	25" (63.5 cm)
Overall length	25.6" (65 cm)	25.6" (65 cm)	25.6" (65 cm)
Depth folded	11" (28 cm)	11" (28 cm)	11" (28 cm)
Wheel diameter	8" (200 mm)	8" (200 mm)	8" (200 mm)
Weight capacity of walker	275 lb (125 kg)	275 lb (125 kg)	275 lb (125 kg)
Weight capacity of shopping bag	11 lb (5 kg)	11 lb (5 kg)	11 lb (5 kg)
Walker Weight	15.4 lb (7 kg)	15.6 lb (7.1 kg)	15.8 lb (7.2 kg)



Triumph Mobility Inc.
www.triumphmobility.com | info@triumphmobility.com | 1.855.546.0711