


Invacare® Replacement Sling (without Commode opening)

Medium	Large	Extra Large	2 Extra Large
R110	R111	R117	-
R112	R113	-	R140
W: 30"	W: 34"	W: 36"	W: 42"
H: 47"	H: 52"	H: 57"	H: 52"
Best Fits: 100-195LB	Best Fits: 150-300LB	Best Fits: 200-400LB	Best Fits: 290-600LB
Capacity: 450LB	Capacity: 450LB	Capacity: 450LB	Capacity: 600LB



Specifications are a guideline only and are subject to change without notice.

Invacare® Replacement Sling (with Commode opening)



Medium	Large	Extra Large	2 Extra Large
R114	R115	R116	R141
W: 32"	W: 35"	W: 37"	W: 42"
H: 48"	H: 50"	H: 56"	H: 50"
Best Fits: 100-195LB	Best Fits: 150-300LB	Best Fits: 210-400LB	Best Fits: 290-600LB
Capacity: 450LB	Capacity: 450LB	Capacity: 450LB	Capacity: 600LB



Specifications are a guideline only and are subject to change without notice.