

Steel Rollator with 6 Wheels, Knockdown

HCPCS: E0143 + E0156



- Strong, durable steel frame supports up to 350 lbs.
- Height adjustable arms (31"-37") allow for a customized fit
- Padded seat allows user to remain comfortable during periods of extended sitting
- Padded backrest adds additional comfort and support when user is in seated position
- Convenient storage pouch under seat for easy access to your personal items
- Ergonomic hand grip with easy to use loop-lock brakes provides comfort and control
- 6" wheels with soft grip tires ideal for indoor and outdoor use
- A compact and foldable design allows you to easily transport your rollator anywhere
- Tool-free quick and easy assembly
- Weighs just 17 lbs.

PRODUCTS**Description**

R800KD-BL	KD Rollator, 6" 4Whl,Steel,Blue,1/ea
R800KD-RD	KD Rollator, 6" 4Whl,Steel,Red,1/ea

SPECIFICATIONS**Description**

Caster Size	6"
Dimensions	22.25" W x 24.8" L x 31.5" H, Basket:12"W x 5" H x 7" L
FoldedSize	6.5" (Depth)
Handle Height Range	31" - 37"
Inside Hand Grip Width	18"
Number of Wheels	4
Seat Depth	12"
Seat To Floor Height	20"
Shipping Weight	20 lbs.
Weight Capacity	350 lbs.
Width inside back legs	18"