

First Class School Chair

HCPCS: E1399-1,E1399



FIGURE A




FIGURE B

- Seat height and depth adjustable solid seat with a 15° anterior or posterior tilt
- Armrests provide lateral limitations as well as a stable resting surface for arms
- Flip-up armrests (Figure A) allow easy transfers and they are height adjustable
- The tray accessory mounts onto the armrest
- Padded Hip Belt 2-point hip belt provides pelvic support and maintains proper hip positioning

PRODUCTS	Description
FC 2000N	FirstClass SchoolChair,Sm,NewStyle,1/ea
FC 4000N	FirstClass SchoolChair,Lg,NewStyle,1/ea

SPECIFICATIONS	Description
Armrest Length	10" 10"
Backrest Height (FC 2000N FC 4000N)	15" 17"
Distance Between Hip Guides (FC 2000N FC 4000N)	8" - 15" 10.5" - 19"
Distance Between Laterals (FC 2000N FC 4000N)	6" - 13" 10" - 17"
Distance Between Armrests (FC 2000N FC 4000N)	13.5" 17.5"
Headrest Added Height (FC 2000N FC 4000N)	4.5" - 8" 4.5" - 8"
Hip Guides (FC 2000N FC 4000N)	3" (H) x 10" (D) 3" (H) x 10" (D)
Lateral (FC 2000N FC 4000N)	4" (H) x 5" (D) 4" (H) x 5" (D)
Product Weight	13 lbs. 18 lbs.
Seat Depth	9.5"-14.5" 13"-17"
Seat To Floor Height	10-14" 15 - 23"
Seat Width	12" 16"
Seat to Footplate (FC 2000N FC4000N)	9.5" - 12.5" 9.5" 17.5"
Seat to bottom of Lateral Pads (FC 2000N FC 4000N)	3.5" - 8.25" 3" - 11.5"
Seat-to-Armrest Height (FC 2000N FC 4000N)	5" - 10" 4" - 13.5"
Shipping Carton Weight (TR 1200)	16.5 lbs. 22.5 lbs.
Weight Capacity	200 lbs
Weight Capacity (13046 13046-HD)	100 lbs. 200 lbs.

ACCESSORIES	Description
FC 2028	Footrest, FC2000, 1pr/bx, BX

ACCESSORIES**Description**

FC 4028	Footrest, FC4000,1/pr/bx, BX
FC 2024	Dining Tray,FC2000,1/bx, BX
FC 4024	Dining Tray, FC 4000,1/bx, BX
FC 8000N	Multi-Axis Headrest, 1/ea, EA
FC 2030N	Anti Tippers, FC2000N, 1pr, PR
FC 2029N	Legs w/Casters,FC2000,1 set/bx, BX
FC 4029N	Legs w/Casters,FC4000,1 set/bx, BX
FC 4030N	Anti Tippers, FC 4000N, 1pr, PR
FC 8027N	Hip Guide, NewStyle, FC, 1pr, PR
FC 8025N	Support Kit, Sm, FC2000N, 1set/bx, BX
FC 8025 LN	Support Kit, Lg, FC4000N, 1set/bx, BX

PARTS**Description**

FC 2002N	Back, FC 2000N, 1/ea, EA
FC 4002N	Back, FC 4000N, 1/ea, EA
FC 8018N	Armrest Compl, FC 2000N / FC 4000N, 1/pr, PR
FC 8060	Screws for Arm Rest, New Style, 1/pkg, PKG
FC ARKNOB	ArmrestLockingKnobs,FC2000N/4000N, 1/pr, PR
FC 8019N	Arm Pad, FC 2000N / FC 4000N, 1/ea, EA
FC 2001N	Back & Accessory Bracket, FC 2000N, 1/ea, EA
FC 4001N	Back & Accessory Bracket, FC 4000N, 1/ea, EA
FC 2011N	Pelvic Belt, FC 2000N, 1/ea, EA
FC 4011N	Pelvic Belt, FC 4000N, 1/ea, EA
FC 2003N	Seat, FC 2000N, 1/ea, EA
FC 4003N	Seat, FC 4000N, 1/ea, EA
FC 0003N	Leg, FC 4000N, 1/ea, EA
FC 0014N	Leg Tip, FC 2000N / FC 4000N, 1/ea, EA
FC KNOBN	Knob Set, FC 2000N / FC 4000N, 4/pkg, PKG