



INVACARE® TOP END®
CROSSFIRE™
ALL TERRAIN
WHEELCHAIR



The **Invacare® Top End® Crossfire™ All Terrain Wheelchair** is the monster truck of wheelchairs and is designed to leave some serious tracks. Do you like to hike, fish, camp, hunt, kayak, garden or go to the beach? You can do it all with the Crossfire All Terrain Wheelchair. Built of 6061 T6 Aluminum, it is equipped with enormous 25" x 2.75" knobby wheels and huge 8" x 2" front casters for extreme traction. The Crossfire All Terrain Wheelchair is engineered with mega ground clearance, plus forward and side-to-side stability so that you can conquer the toughest terrain.

Features

6061 T6 Aluminum frame
Adjustable center-of-gravity
Camber 0-9°
Lock down folding back
Adjustable back height
Adjustable tension seat upholstery
Adjustable tension back upholstery

24" or 25" high performance Sun® wheels
Hard anodized handrims
24" or 25" x 2.75" knobby tires
8" x 2" front casters
Adjustable height and angle footplate
Powder coat paint
250 lb. weight capacity

Model #: CRFAT



1. 24" or 25" x 2.75" knobby wheels and tires
2. 8" x 2" front casters
3. Engineered for front and side-to-side stability
4. Center-of-gravity adjustment
5. Anti-tip option
6. Optional push handles, wheel locks and anti-tip bars



Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Top End

www.topendwheelchair.com
Sports and Recreation Products
4501 63rd Circle North
Pinellas Park, Florida
33781-5914
(727) 522-8677
(800) 532-8677
Fax: (727) 522-1007

©2014

Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols™ and®.
All trademarks are owned
by or licensed to Invacare
Corporation unless
otherwise noted. Natural
Fit is a registered trademark
of Three Rivers Out Front.
Spinergy is a registered
trademark of Spinergy, Inc.
Specifications are subject
to change without
notification.
Form No. 12-504
Rev. 02/14
140182

Options

Spinergy® Spox wheels
Quad release axles
Titanium axles
Coated pushrims
Titanium handrims
Natural Fit® handrims
Fixed back
Push handles
Carbon fiber sideguards

Removeable sideguards
Folding sideguards
Wheel locks
Positioning straps
Alignment gauge
Footplate risers
2" cushion
Anti-tip bars



Yes, you can.®