



Multi-Positioning Seat

For the versatility you need!

The Multi-Positioning Seat offers comfortable and proper positioning support. Made to be used almost anywhere – from home and school to bus and in the community. The modular cushion seating surface allows for individualized adjustment and years of growth. Seamless, latex-free cushions are impermeable to fluids and soft, yet durable. It withstands cracking and cleans easily with common disinfectants. Attachment straps can be used to attach this seat to standard chairs. Customize the MPS by selecting a Standard Headrest or a Headrest with Laterals. Now Available in three color combinations: Aqua Blue & Chocolate, Lilac & Chocolate, and Chocolate & Chocolate. Latex free. 🍷

Choose from 3 different systems:

The Multi-Positioning Seat (MPS), MPS Mobile Floor Kit and MPS Mobile Push Chair Kit. Each includes:

- Adjustable seat and back cushions
- Seat shell
- Adjustable Footrest
- LATCH System for car seat usage
- 5-year manufacturer's warranty

Multi-Positioning Seat (MPS). Headrest not included.

Aqua Blue/Chocolate	Lilac/Chocolate	Chocolate/Chocolate	Size of MPS Seat with Footrest
9288-02	9288-24	08-1451-509	Small
9288-07	9288-29	08-1451-608	Small Extend
9288-13	9288-35	08-1451-574	Large
9288-18	9288-40	08-1451-665	Large Extended

MPS Mobile Floor Kit. Includes MPS and Mobile Floor Base. Headrest not included.

Aqua Blue/Chocolate	Lilac/Chocolate	Chocolate/Chocolate	Size of MPS Mobile Floor Kit
9288-00	9288-22	08-1451-475	Small
9288-05	9288-27	08-1451-483	Small Extended
9288-11	9288-33	08-1451-616	Large
9288-16	9288-38	08-1451-640	Large Extended

MPS Push Chair Kit. Includes MPS and Push Chair Base (9288-44). Headrest not included.

Aqua Blue/Chocolate	Lilac/Chocolate	Chocolate/Chocolate	Size of MPS Push Chair Kit
9288-01	9288-23	08-1451-491	Small
9288-06	9288-28	08-1451-533	Small Extended
9288-12	9288-34	08-1451-624	Large
9288-17	9288-39	08-1451-657	Large Extended

Headrests. (Required) Fits both sizes of the MPS.


Aqua Blue	Lilac	Chocolate	Description
9288-62	9288-63	08-1451-814	Standard Headrest*
9288-64	9288-65	08-1451-822	Headrest with Laterals*

*As of March 2006, a headrest must be purchased separately.

Footrests. Included with the MPS, MPS Mobile Floor Kit and MPS Push Chair Kit.

Chocolate	Description
9288-58	Small Footrest
9288-59	Large Footrest

Shown here with Standard Headrest (sold separately)



9288-65
The standard headrest or the headrest with laterals are required and must be ordered separately. Headrest with laterals shown here.

Use the Special Tomato MPS system everywhere!

- 🏠 **At home**, secures to most chairs by using attachment straps
- 🚌 **In cars**, use the crash-tested seat with LATCH & Tether system
- 🏫 **At school**, secure the Special Tomato in a classroom chair, or use it with the Mobile Base as your all-in-one classroom chair
- 🚗 **On the school bus**, transport the Special Tomato in Push Chair (9288-44, shown on the facing page)
- ✈️ **In an Airplane**, both sizes approved by FAA for air travel

MPS Accessories



MPS Push Chair Base. Umbrella-style folding stroller for easy transport. Folds down to 48" L x 7" W x 7" D. Supports 150 lbs. and weighs only 20 lbs. 30° fixed seat tilt. Designed for use with standard tie-down systems for bus or van transport and meets ANSI/RESNA WC 19 regulations. 🍷
9288-44



MPS Sun Shade for Push Chair Base. Provides protection from the sun and a convenient window for monitoring!
9288-45



Tomato Skins. An optional black cushion cover that custom-fits the Special Tomato MPS. Made of soft, durable fabric to absorb moisture and provide added comfort. Not approved for use in a car.
5507-46 Small
5507-47 Large



Mobile Floor Base. Built-in tray storage on back. Four swivel casters with rear locks. 10°-25° range of tilt-in space. 🍷
Chocolate Size
9288-04 Small
9288-15 Large



Headrest sold separately

five year warranty



MPS Extended Seat Cushion. As your child outgrows the standard cushion, simply purchase the extended cushion.
Aqua Blue Lilac Chocolate Size
9288-08 9288-30 08-1451-541 Small
9288-19 9288-41 08-1451-558 Large



MPS Tray. Depth adjustable. 60° tilt adjustment. Tray easily attaches to back of Mobile Floor Base for storage.
Chocolate Size
9288-03 Small
9288-14 Large

MPS Push Chair Base Accessories
9288-46 8" Pneumatic Tires Upgrade-Set of 4
9288-47 Push Chair Travel Bag

Seating System Sizing Information

Special Tomato Seats	Small		Large	
	Standard Seat	Extended Seat	Standard Seat	Extended Seat
Target User Age	2-5 years	5-8 years	8-11 years	11-14 years
Maximum User Wt.	80 lbs.	80 lbs.	150 lbs.	150 lbs.
Seat Weight	22 lbs.	24 lbs.	26 lbs.	28 lbs.
Seat & Floor Base Wt.	40 lbs.	42 lbs.	44 lbs.	46 lbs.
Seat Width	9"	10"	12"	13"
Seat Depth	8"	10"	13"	15"
Seat to Footrest	6"-16"	6"-16"	8"-18"	8"-18"
Top of Back Cushion	15"-19"	15"-19"	17"-21"	17"-21"
Top of Head Cushion	22"-27"	22"-27"	24"-29"	24"-29"
Hip Flex Angle	-15° to +15°	-15° to +15°	-15° to +15°	-15° to +15°
Floor to Seat Ht.	15"	15"	18"	18"
Tray	24"W x 17"D x 2"H	24"W x 17"D x 2"H	24"W x 17"D x 2"H	24"W x 17"D x 2"H
MPS Seat Only	15"W x 19"D x 42"H	15"W x 19"D x 42"H	19"W x 23"D x 47"H	19"W x 23"D x 47"H
MPS Floor Kit	19"W x 37"D x 37"H	19"W x 37"D x 37"H	22"W x 41"D x 42"H	22"W x 41"D x 42"H
MPS Push Chair Kit	25"W x 38"D x 42"H	25"W x 38"D x 42"H	25"W x 42"D x 47"H	25"W x 42"D x 47"H