

Kanga Tilt-in-Space Wheelchair

HCPCS:



FIGURE A



FIGURE B

- Newly designed reinforced frame
- Available with new optional (sold separately) rehab style seating components with multiple adjustments
- Newly designed height adjustable & removable "T" style armrests
- Frame can grow to the next size without additional components
- Folding frame tilt-in-space wheelchair with 45° of tilt
- Height & angle adjustable back canes
- Width adjustable solid seat pan included
- Removable anti tippers included
- Transit brackets included for securing unoccupied wheelchair in vehicle

PRODUCTS

Description

KGP-8010-2GR	Kanga Tilt in Space, 10" W, Adj Depth
KGP-8012-2GR	Kanga Tilt in Space, 12" W, Adj Depth
KGP-8014-2GR	Kanga Tilt in Space, 14" W, Adj Depth

SPECIFICATIONS

Description

Weight Capacity	200 lbs
-----------------	---------

ACCESSORIES

Description

KGP-FR-AA	Kanga Ped Footrests, Adj Angle, Pr., EA
KGP-ELR-2G	Kanga Pediatric Elevating Legrests, Pr., EA
KG-RW20	Kanga 20" Rear Wheel Assembly, Pr., EA
1300-921112	13009-SBA 13003-SA Seating Pkg, 12in, EA
1300-921114	13009-SBA 13003-SA Seating Pkg, 14in, EA
1300-921110	13009-SBA 13003-SA Seating Pkg, 10in, EA
4008-401010	4008-4010 Headrest, Offset Hardware, 10in, EA
4008-401012	4008-4010 Headrest, Offset Hardware, 12in, EA
4008-401014	4008-4010 Headrest, Offset Hardware, 14in, EA
11060-10110	11060-101 Abductor, Flip Down, 10in, EA
11060-10112	11060-101 Abductor, Flip Down, 12in, EA
11060-10114	11060-101 Abductor, Flip Down, 14in, EA
11072-34410	11072-344 Thigh Supports, Rem., 10in Pr, EA
11072-34412	11072-344 Thigh Supports, Rem., 12in Pr, EA
11072-34414	11072-344 Thigh Supports, Rem., 14in Pr, EA
3040-021010	3040-0210 Trunk Supports, Rem., 10in Pr, EA
3040-021012	3040-0210 Trunk Supports, Rem., 12in Pr, EA

ACCESSORIES**Description**

3040-021014	3040-0210 Trunk Supports, Rem., 14in Pr, EA
12011-10	Upper Extremity Support, Polycarb, 10", EA
12011-12	Upper Extremity Support, Polycarb, 12", EA