

Roscoe Mobility Knee Scooter

The Knee Scooter offers increased comfort and ease of mobility to patients who cannot put weight on their foot or ankle, but want to remain active. With its four large wheels and padded knee platform, the Knee Scooter is perfect for indoor and outdoor use. The offset post on the aluminum knee platform allows for right-to-left adjustment in seconds - no tools required. The Knee Scooter is equipped with an easy folding mechanism that collapses using a thumb release lever and fits beneath the raised knee platform for ease of storage and transportation.

Black Finish



White Finish



NEW!
Burgundy Finish



Product Features

- Scooter can be steered for increased maneuverability
- Ideal for individuals recovering from an ankle or foot injury
- Large 8-inch clear PU wheels with sealed bearings
- Adjustable hand brake ensures safety
- Padded knee platform with easy-adjustment mechanism
- Padded handlebar grips
- Limited lifetime warranty on frame, 90 days on all other parts

Product Specifications

- Knee Platform Measurements: 7" (W) x 14" (L)
- Frame Finish: powder-coated black, white or burgundy
- Adjustable Handle Height: 31" - 41"
- Knee Pad Height Adjustment: 17 1/2" - 22"
- Product Weight: 22.8 lb.
- Weight Capacity: 350 lb.



**Knee Scooter
Replacement
Parts Available!**

New and Improved Maintenance Free Axle and Axle Bearings!



Old Rear Axle

NEW Rear Axle



Adjustable Hand Brake



Handle Bar can be Adjusted

Item #	Description	Patient Height Range	Color	Qty
ROS-KSB	Knee Scooter with 8-Hole Stem	4'11" to 6'6"	Black	each
ROS-KSW	Knee Scooter with 8-Hole Stem	4'11" to 6'6"	White	each
ROS-KSBG	Knee Scooter with 8-Hole Stem	4'11" to 6'6"	Burgundy	each



To Order:
(800) 3-ROSCOE
fax: (440) 572-4261
roscoemedical.com

Roscoe Mobility

Gemini Scooter

The Gemini Scooter is ideal for patients recovering from a foot, ankle, hip or knee injury or illness. The large gel padded seat and convenient front leg peg allows the patient to rest in comfort while elevating the non-weight bearing leg. The Gemini converts from a seated scooter to a knee scooter with the optional knee platform (sold separately) in less than one minute, no tools required. Equipped with an easy-folding mechanism, the Gemini Scooter collapses using a thumb release lever for ease of storage and transportation. It is an ideal rental item as it allows providers to generate more revenue with the versatility of one scooter and carry less inventory.

Indications:

- Knee Scooter: Foot/ankle injuries such as fractures, dislocation, heel spurs
- Seated Scooter: Hip, knee and foot/ankle injuries
- Post-operative foot/ankle surgery
- Diabetic foot/toe ulcers
- Neuromuscular problems or stroke
- Arthritic conditions


Product Features:

- Padded seat allows patient to rest in comfort while elevating the non-weight bearing leg on front foot peg
- Patient has the option of using the gel padded seat or padded knee platform (sold separately) depending on the injury
- Folds compactly for easy storage and transport
- Converts from a seated scooter to a knee scooter in less than one minute, no tools required (knee platform sold separately)
- Quick-release basket included
- Large 8-inch rubber wheels with sealed bearings; great for indoor or outdoor use
- Smooth movement with no hopping or skipping
- Unique steering system with 5 foot turning radius
- Padded handlebar grips for added comfort
- Dual braking system
- Parking break for additional safety
- Triple coated painting process for better durability
- Limited lifetime warranty on frame
90 days on all other parts

Product Specifications:

- Knee Platform Measurements: 7"(W) x 14"(L)
- Knee Pad Height Adjustment: 18" - 23 1/2"
- Adjustable Seat Height: 25" - 34 1/2", with one inch increments
- Adjustable Handle Height: 31 3/4" - 41", with one inch increments
- Product Weight: 23.8 lb.
- Weight Capacity: 250 lb.
- HCPCS Code: E0118* crutch substitute





Seat Extender
Now available with a wider range of adjustment to accommodate patient height!



Item #	Description	Patient Height Range	Color	Qty
30188	Gemini Scooter with Padded Seat	4'11" to 6'6"	Black	each
30190	Optional Knee Platform for Gemini Scooter	4'11" to 6'6"	Black	each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



To Order:
(800) 3-ROSCOE
fax: (440) 572-4261
roscoemedical.com