

Roscoe Mobility

Roscoe Economy Knee Scooter

The Roscoe Economy Knee Scooter offers an affordable choice in mobility for those suffering from an ankle or foot injury. The Roscoe Economy Knee Scooter can easily be folded down for transportation and storage. The handle bars and knee platform are height adjustable to accommodate a wide range of users. The frame and knee pad are manufactured for worry free use.



Product Features

- Scooter can be steered for increased maneuverability.
- Ideal for individuals recovering from an ankle or foot injury.
- Quick Release Steering Assembly folds for easy storage.
- Quick Release Seat includes height adjustments to accommodate different patient heights.
- Padded Knee Platform with easy adjustment mechanism.
- Oversize wheels are ideal for indoor and outdoor use.
- Adjustable braking system.
- Limited Lifetime warranty on frame; 90 days on Non-Durable components.

Product Specifications

- Product Measurement: 31 3/4" (H) x 16 1/8" (W) x 30" (D)
- Knee Platform Measurement: 14" (L) x 7" (W)
- Knee Pad Height Adjustments: 13 3/4" – 19 1/4"
- Frame Finish: Sterling Grey Liquid Paint
- Large 7 1/4" wheels with sealed bearings
- Adjustable Handle Height: 31 3/4" – 39 3/4"
- Product Weight: 19.5 lbs.
- Weight Capacity: 250 lbs.



Lockable Hand Brake



Height Adjustable Knee Platform



Padded Knee Platform

Item #	Description	Patient Height Range	Color	Qty	UPC
ROS-KS2	Roscoe Economy Knee Scooter	4'11" to 6'6"	Sterling Grey	each	092237614624



To Order:
 (800) 3-ROSCOE
 fax: (440) 572-4261
 roscoemedical.com