

Viper Wheelchair

HCPCS: K0004,E1236+K0195+E0990



FIGURE A



FIGURE B

- Viper weighs 32 lbs (Excluding front rigging)
- Easier to propel and transport
- Frame style eliminates seat guides and allows for custom backs and accessories
- Back height adjustable from 18" 20"
- Carbon steel black powder-coated frame double coated for an attractive, chip-proof, maintainable finish
- Durable, flame retardant nylon upholstery resists mildew and bacteria
- Composite Mag-style wheels with aluminum hand rims are lightweight and maintenance free
- Urethane rear tires and casters offer superior performance, smoother ride and are lightweight
- 8" Front casters have three height adjustments
- Nylon, padded armrests provide added patient comfort
- Flip-back arms allow for easy transfer
- Composite steel lined front receiver is angled so the arm will lock into place every time
- The front mounted release lever can be flipped up or down to unlock flip-up arm, depending on the patients needs (Figure A - B)
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Dual axle provides easy transition of seat height to hemi level
- Comes with swing-away footrests or elevating leg rests with tool-free adjustable length riggings
- Composite foot plates with heel loops are durable and lightweight
- Comes standard with flip-up anti-tippers
- Comes with push-to-lock wheel locks
- Available in 12" and 14" seat widths

PRODUCTS

	Description
L412DDA-ELR	Viper,12"FlpBk.DskArm,ELR,1/cs
L412DDA-SF	Viper,12"FlpBk.DskArm,SF,1/cs
L414DDA-ELR	Viper,14"FlpBk.DeskArm,ELR,1/cs
L414DDA-SF	Viper,14"FlpBk.DskArm,SF,1/cs

SPECIFICATIONS

	Description
Back Height	16" - 18"
Carton Shipping Weight	39 lbs or less
Height with 12 and 14 Seats	32"
Height with 16 , 18 and 20 Seat	36"
Length with Riggings with 12 and 14 Seats	36"
Length with Riggings with 16 , 18 and 20 Seat	42"
Product Weight with 12 Seat	30 lbs
Product Weight with 14 Seat	31 lbs

SPECIFICATIONS	Description
Product Weight with 16 Seat	31.5 lbs
Product Weight with 20 Seat	34 lbs
Product Weight with 18 Seat	32 lbs
Seat Depth with 12 Seat	12"
Seat Depth with 14 Seat	14"
Seat Depth with 16 , 18 and 20 Seat	16"
Seat To Floor Height	13.5" - 15.5" 17.5" - 19.5"
Warranty	Limited Lifetime
Weight Capacity	For 12" and 14" seat: 250 lbs For 16", 18" and 20" seat: 300 lbs
Width Closed	12"
Width Open	20" 22.5" 22.5" 24.5" 28"

ACCESSORIES	Description
STDS860	Forwrd Stablzrs/Univ.Fit , 1pr/bx, BX
STDS1818N	Seat Rail Ext.Kit,18"x18", 1/ea, BX
STDS2018N	Seat RailExt Kit,20x18,Blk.Nylon,1/ea, EA
STDS840	RearWheel,TreadedUrethane, Viper GT 1/ea, EA
STDS807	Anti-Tips,Adj.w/Whls,Blk,1pr/bx, BX
STDS1034	Wheelchair Cane/Crutch Holder,1/ea, EA
WASR	Limb Support, Right,1/cs, CS
WASL	Limb Support,Left,1/cs, CS
STDS801	Brake Handle Ext.,8", 1/ea, EA
STDS820	Telescoping IV Pole,W/C, 1/bx, BX
STDS834	OverheadAnti-Theft Bar,Single,1/ea, EA
STDS856	Seat Belt,60"Velcro-type,1/ea, EA
STDS851	Seat Belt, Velcro Close,1/ea, EA
STDS855	Seat Belt,60"" Auto-type,1/ea, EA
STDS850	Seat Belt, Clasp Type,1/ea, EA

PARTS	Description
STDS3J4502	Footrest Adj.Lever,SlvrSprt.1/ea, EA
STDS1011	Wheel Bearing, 8mm x 22mm 1/ea, EA
STDS1001	Lightweight Wheel 24inAlum Handrail,1/ea, EA
STDS1002-1	Axle for 8" caster, 1/ea, EA
STDS1002	Caster,8in Sentra,Viper,Astaire,1/ea, EA
STDS1060	Wheelchair Fork Bearing,1/ea, EA
STDS831	Heel Loops, 1pr, BX
STDS1026	Front Bumper for Ft.Rst/LegRest,1/ea, EA
STDS871	Hand Grip,Sentra,1/ea, EA
STDS1051	Dust Cap, fork stem, 1/ea, EA
STDS4011	BrakeTip, 1/ea, EA
STDS2358BLK	Seat Uphol. 18"x16",Black,1/bx, BX
STDS2356BLK	Seat Uphol. 16"x16",Black,1/bx, BX
STDS1073	Spacers & Screws for Handrim, 8/set, EA
STDS2A3141	Seat Rail Brckt,Rt.Frnt,Alum.1/ea, EA
STDS2363BLK	Back Uphol 16"x18",Black,1/bx, BX
STDS2359BLK	Back Uphol 18"x16",Black,1/bx, BX
STDS2361BLK	Back Uphol 20"x18",Black,1/bx, BX
STDS3003	Caster Fork,Whch,1/ea, EA
STDS004011	Handrim,Alloy24"Wheel/1ea., EA