S



Series PR501

Comforter The Comforter Series from Golden Technologies is all about finding the perfect size for the perfect fit. The Comforter Series offers five sizes to comfortably fit people from 5'0" to 6'8" with weight capacities ranging from 300 to 375 lbs. Featuring a plush, seamed backrest with lumbar support, a soft, comfortable seat and our user-friendly "grip rite" arm design, the Comforter series is your ideal lift chair solution.







= pocketed coil spring seat std.











- /	Ser.	100		
		BA		
ŀ			A	
			7	
		M		۽
0				501 M shown
0				2
1				=
1				윤

Golden Craftsmanship

Each luxury lift & recline chair is hand crafted in our facility using decades of motion furniture experience.



Added Comfort and Safety

Designed to perform year after year using advanced technology backed by the best warranty.



Smart Engineering

Smooth, quiet lifting system offers a secure lift for the best performance.



Smart Tek

In the event of a malfunction, our patented SmartTek™ diagnostic system enables you to get your lift chair up and running quickly.

	Comforter	Comforter	Comforter	Comforter	Comforter
Model Number	PR-501JP	PR-501S	PR-501M	PR-501L	PR-501T
Size	Junior Petite	Small	Medium	Large	Tall
Weight Capacity	300 LB	300 LB	375 LB	375 LB	375 LB
Recline Positions	3	3	3	3	3
Standard Fabrics	5	5	5	5	5
Fabric Options	ALL	ALL	ALL	ALL	ALL
Chaise Seat	YES	YES	YES	YES	YES
Coil Springs	NO	YES	YES	YES	YES
Back Type	Seam	Seam	Seam	Seam	Seam
Head Rest Cover	YES	YES	YES	YES	YES
Arm Covers	YES	YES	YES	YES	YES
Battery Back Up	YES	YES	YES	YES	YES
Overall Width	30.5"	31.5"	32.5"	34.5"	34.5"
Overall Height	39"	40"	42"	44"	47"
Width Between Arms	19"	20"	21"	23"	23"
Floor to Top of Seat	18"	19"	19.5"	20.5"	22.5"
Seat Depth	18.5"	19.5"	21"	22"	24"
Seat to Top of Back	25"	25"	27"	28.5"	30"
Seat to Top of Arm	6.5"	6"	7"	7.5"	6"
Seat Angle Raised	16°	15°	16°	16°	18°
Floor to Top of Foot Rest Recline	d 22"	23"	23"	23"	23.5"
Floor to Top of Back Reclined	32"	32"	32"	33"	36"
Back Angle Closed	110°	110°	110°	110°	110°
Back Angle Reclined	145°	145°	145°	145°	145°
Foot Rest Angle Closed	90°	90°	90°	90°	90°
Foot Rest Angle Opened	180°	180°	180°	180°	180°
Distance Required from Wall Red	lined 17"	17"	17"	18"	19"
Extension Reclined from Back					
of Seat to End of Foot Rest	38"	39°	42"	43"	45"
Life Time Warranty, Chair Frame,		VEC	VEC	VEC	VEC
Lift Frame, Recline Mechanism	YES	YES	YES	YES	YES
Warranty on Electrical Parts	3 Yrs.	3 Yrs.	3 Yrs.	3 Yrs.	3 Yrs.
Pro Rated Years 4 Thru 7	YES	YES	YES	YES	YES
Warranty on Mechanical Labor	3 Yrs.	3 Yrs.	3 Yrs.	3 Yrs.	3 Yrs.
Warranty on Electrical Labor	1 Yr.	1 Yr.	1 Yr.	1 Yr.	1 Yr.
Removable Back Standard	YES	YES	YES	YES	YES
Smart Tek System Standard	YES	YES	YES	YES	YES

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Battery Backup System

Customizable

Left Hand Control Option Available

Heat and Massage Option Available

Gel Memory Foam Option Available

Foot Rest Extension Option Available

Extra Pocket Option Available

Straight Lift Option Available

User Height (suggested)

Battery Backup System

YES

YES

YES

YES

NO

NO

YES

5'0" or under

Batteries not included. Battery Backup System is designed to bring the chair from the reclined to the seated position in the event of a power outage.



YES

YES

YES

YES

YES

YES

YES

5'11"-6'2"

YES

YES

YES

YES

YES

YES

YES

6'3"-6'8"

Custom Sizing...

YES

YES

YES

YES


YFS

YES

YES

5'4"-5'10"

One size doesn't fit all with Golden Technologies seat lift chairs.





YES

YES

YES

YES

NO

YES

YES

5'0"-5'3"

Complies with UL 73 and CSA C22.2 No. 68

