



# PLUS

“Driving Performance” is more than its appearance in most of Luggie. From its elegant appearance to its safety performance, an unmatched breadth of capabilities reveals the depth of character promised by the tag in its back.

*Luggie*<sup>®</sup>

2  
0  
1  
9



## EFFICIENCY

Improving high efficiency substance to support you long mileages and lower your maintain cost.



## RELIABILITY

Luggie Line is FDA Approval, period.



## SAFETY

Intelligent designs to ensure your riding is safety.

Luggie **Elite** is featured 320lbs, extendable, handle bar, and 16" seat. Plus 3 and Plus 4 are available to Elite line. Equipped with 10.5ah battery.



[360 view](#)

SAME: LED battery indicator, seat height adjustable, armrest, SuperCharger



[360 view](#)


Luggie **Super** is featured 360lbs, extended handle bar, and 17" seat. Luggie Plus 3 and Luggie Plus 4 are available to Super line. Equipped with 16.5ah battery and 8.5ah battery.



PLUS 4 is 3-wheel  
turning radius  
but 4-wheel stability

NO NEED to compromise your  
stability demands for true 4-wheel  
scooter.

Plus 4 is the smallest turning radius  
in the market, only 36".



## Rides with refinement

Riding experienced is improved by the patented Omni-suspension which is coated rose gold color to tell the others apart.

## Corners with confidence

### SMART TURN

HANDS-FREE speed control.

Mounted to a rigid and lightweight subframe, assist precise speed controlling response by turning angle. This is effective to assure the riders feel more safety in varying riding scenario.







## DEVELOP TO OUTPERFORM

Luggie is the most efficient mobility scooter. It supports you the longest range without carry a big or heavy battery. Fully charge takes you only 90 min and give you another 10 mile. It's a real travel scooter, no more battery anxiety.

## BORN TO PERFORM

Luggie is the only folding scooter comes with folding variations to help riders **LEVERAGE** the scooter easily. Smallest footprint to help storage in most of constrain space.



Luggie PLUS 4 supports the **HALF FOLDING** to leverage the scooter into trunk without lifter to save your cost and space.





Enlarge 8" solid tire to the rear.

## FASHION THAT FUNCTIONS

A long wheelbase, big and wide wheels, and wide track add up to space, stability and sportiness for *Luggie Super Plus*.

*Richly taste equipped with your Luggie Plus.*

- *Rose Gold Suspension*
- *Black Tires*
- *Metallic Color*



S U P E R  
P L U S 3



5

TOP SPEED

TURNING RADIUS

37"

25<sub>MI</sub>

RANGE PER CHARGE



Model	Elite PLUS 3	SUPER PLUS 3
Weight Cap/ lbs	320	360
Seat Width/ inch	16	1
Length / inch	39	39
Width / inch	21	23
Top Speed	4mph	5mph
Seat Width	16"	17"
Scooter Weight / lbs. (w/o battery)	58	66
Range / (up to miles) @ 16.5AH li-ion battery	25	25
Turning Radius	36"	36"
Ground Clearance @ side deck/ inch	4.5	4.5
Flat Folding	Yes	Yes
W/armrest	Yes	Yes
Suspension	YES	Yes

# Appearance



Lake Green



Caribbean Blue



Titanium Silver



Mineral Gray



Mustard Olive

S U P E R  
P L U S 4



5

TOP SPEED

TURNING RADIUS

36"

25<sub>MI</sub>

RANGE PER CHARGE



Model	Elite PLUS 4	SUPER PLUS 4
Weight Cap/ lbs	320	360
Seat Width/ inch	16	17
Length / inch	39	39
Width / inch	21	23
Top Speed	4mph	5mph
Scooter Weight / lbs (w/o battery)	58	66
Range / (up to miles) @ 16.5AH li-ion battery	25	25
Turning Radius/ inch	36	36
Ground Clearance/ inch @ side deck	4	4
Flat Folding	Yes	Yes
W/armrest	Yes	Yes
Suspension	YES	Yes

# Appearance



Mustard



Caribbean Blue



Titanium Silver



Mineral Gray



Lake Green

A blue and black motorized scooter with a black seat and handlebars. The scooter is shown from a side profile, facing right. It has a blue frame and black wheels. The seat is black with a small logo on the backrest. The handlebars are black and have a silver-colored stem. The scooter is positioned on a dark grey background.

2019

JOIN US

FREERIDER USA  
8696 Utica Ave Rancho Cucamonga, CA 91730  
Email: [sales@freeriderusa.com](mailto:sales@freeriderusa.com)  
TEL: 909-466-7996  
FAX: 909-466-7990