

# TRAV

## Enjoy more mobility with a lightweight, portable wheelchair

Now, you can go anywhere with the Trav. Compact and lightweight, the Trav makes a great supplement to your primary wheelchair. The Trav handles easily in narrow spaces, so it's perfect for use in restaurants and theaters, and while traveling by car, bus, train or plane. Best of all, the Trav is built to the same high standards of quality and safety as all Nissin wheelchairs, giving you the freedom to go more places while also enjoying peace of mind.



## Compact and lightweight



### The Trav...

- ✓ weighs less than 13 lbs
- ✓ folds easily with the touch of one button
- ✓ folds to meet airline requirements for carry-on bags
- ✓ carry bag\* included for convenient storage and carrying

HCPCS Code E1038


## Anytime, anywhere, one-touch operation

### Quality engineered for comfort, mobility and safety

Nissin engineered the Trav with the help from researchers at a major rehabilitation center in Japan with the needs and desires of their clients in mind. The rear wheels are installed inside the chair's frame to reduce its width, increase maneuverability, and yet provide the comfort of a full-size wheelchair. The thigh belt ensures balanced seating. The flip-up armrest facilitates smoother transfers from the chair's seat. The rear wheels have a step break for the convenience of caregivers. The Trav allows you to enjoy more mobility while feeling safe and secure.

### Specifications

|                 | Open     | Folded |
|-----------------|----------|--------|
| Weight Capacity | 250 lbs. | N/A    |
| Height          | 33"      | 21"    |
| Width           | 15.75"   | 15.75" |
| Length          | 25.6"    | 7.5"   |



Nissin Wheelchairs USA  
 3551 Voyager St., Unit 102  
 Torrance, CA 90503  
 Tel: 310-370-8484  
 (877) NISSIN7  
 Fax: 310-370-8454  
 Email: info@nissinusa.net  
 www.nissinusa.net

Please read all the instructions carefully before using. Specifications are subject to change without notice due to constant product innovation and improvement.  
 \*Carry bag pictured is not the free bag that is included, but is available for purchase.