

Specifications

Colour	Model Numbers	
■ Sky Blue ■ Midnight Blue	600-210006 600-210007	600-190006 600-190007
Overall width	24.75" / 63 cm	24.75" / 63 cm
Folded width	10.5" / 26.7 cm	10.5" / 26.7 cm
Length	28.25" / 71.75 cm	28.25" / 71.75 cm
Handle height range	32.4" – 40.2" 82.3 cm – 102 cm	29.4" – 37.2" 74.7 cm – 94.5 cm
Distance between handles	17.75" / 45 cm	17.75" / 45 cm
Seat Height	21" / 54 cm	19" / 48 cm
Seat width	16" / 40.6 cm	16" / 40.6 cm
Seat depth (rollator)	12.5" / 31.75 cm	12.5" / 31.75 cm
Seat depth (transport chair)	16" / 40.6 cm	16" / 40.6 cm
Weight (Transport)	30 lbs (13.6 kg)	30 lbs (13.6 kg)
Weight (Total)	32 lbs (14.5 kg)	32 lbs (14.5 kg)
Weight Capacity	275 lbs (125 kg)	275 lbs (125 kg)




SKY BLUE



MIDNIGHT BLUE



Limited lifetime warranty
 Please visit:
triumphmobility.com
 for additional information



SCAN ME!

Contact us

Triumph Mobility
 151 Bentley Ave, Unit 4
 Nepean, Ontario, K2E 6T7
 Tel: +1 (855) 546 0711
 Fax: +1 (866) 225 3059
 Email: info@triumphmobility.com
www.triumphmobility.com



TriumphMobility

©2020Triumph Mobility Inc. All rights reserved.
 Triumph™ and Prestige™ are trademarks of
 Triumph Mobility Inc. The Triumph Prestige
 rollator and transport chair is covered by pending
 patent and industrial design protection.

Printed in Canada
 File name: R_202204_TriumphPrestige_brochure_EN

Proudly designed in Canada




TRIUMPH
 prestige.

The Prestige Lifestyle

- ✓ No tools needed
- ✓ No fussy instructions

Finally, a support to keep pace with your busy life.



ROLLATOR MODE

TRANSPORT CHAIR MODE



TRIUMPH
prestige



Perfect Flexibility

10.5 in



The Prestige locks closed at **10.5 in** at its widest point – perfect for tucking it away in a restaurant or at home. Its economic size also makes it easy to fit in a car trunk, weighing only **13.6 kg**



The Prestige transforms from rollator to transport chair in seconds: **swivel the flexible padded back rest, flip the handgrips around**, you're ready to move!



Don't worry about wearing it out: the frame has a **lifetime warranty** and all other parts are guaranteed for 2 years.



≤ 11 lbs (5 kg)

The basket underneath makes the Prestige perfect for shopping, securely attached at six different points and strong enough to **hold up to 11 lbs.**

What else?

Comfortable



- The padded backrest adjusts to the natural angle in both modes
- A padded seat cushion with a roomy seat depth of 12.5" in rollator mode and 16" in transport chair mode
- Ergonomic, angled and adjustable handgrips
- Padded armrest at a natural position of 9.5" from seat to arm

Stylish



- Two fashionable colours to choose from
- Wires all internal for a sleek look with nothing to snag

Safe



- Soft tires provide a natural suspension and balance
- Drum brakes provide reliable all-weather braking
- Sturdy aluminum frame design