

ROBOTEER

R E B O O T Y O U R L I F E


EXCLUSIVELY DISTRIBUTED BY: [meducare](https://www.meducare.com)




Robooter E60/E60-A

ELECTRIC WHEELCHAIR

The next generation of performance has arrived



[meducare](https://www.meducare.com)



www.meducare.com



Adjustable suspension
A smoother ride has arrived
Comfort all day



Omni-wheel designed for low noise and vibration
360 turning means no k-turns
Works great in all road conditions
including gravel, mud, grass, etc



Fast charging, long range



Ergonomic Seat and Back



○ Pearl White



● Stone Black



● Lake Green

PARAMETER

PRODUCT MODEL	E60	E60-A
Folding mode	Manual folding	Manual folding
Maximum speed	5 mph	5 mph
Maximum barrier height	2.3"	2.3"
Maximum crossing width	4"	4"
Slope climbing performance	9°	9°
Total weight	71.4 lbs.	75.4 lbs.
Backrest angle	Not adjustable	Adjustable
Battery weight	7.5 lbs.	7.5 lbs.
Specifications of front and rear wheels	10-inch front wheel, 10-inch rear wheel	10-inch front wheel, 10-inch rear wheel
Maximum load	< 330 lbs.	< 330 lbs.
Battery capacity	20 Ah/25 Ah	20 Ah/25 Ah
Overall dimension (length * width * height)	38.9"x25.2"x36.7"	41.4"x25.2"x37.6"
Folding shape dimension (length * width * height)	30.3"x25.2"x19.7"	33.7"x25.2"x20.7"
Seat wide	16.5"	16.5"
Theoretical driving distance	14.3 miles(20 Ah)、18 miles(25 Ah)	14.3 miles(20 Ah)、18 miles(25 Ah)