

# Durable medical equipment **product guide**



Innovative products and support



**Elevate care  
with medical  
equipment  
that goes  
the distance**

High-quality engineering, long-lasting performance, and a comprehensive and unique product offering makes Medline your best choice for durable medical products that enhance a user's lifestyle and quality of care.












We're committed to delivering outstanding customer value and providing high-quality health care products. With strategically located distribution centers across the country and a large fleet we can get the products you need, when you need them.

Your Medline Representative is eager to assist you. From product selection to smart inventory suggestions, you can count on us.



**Find many new, industry-leading products throughout this catalog.** Look for the items with this symbol: **New**

**What's inside:**

|   |                                     |    |
|---|-------------------------------------|----|
|  | Crutches .....                      | 6  |
|  | Knee Walkers .....                  | 8  |
|  | Canes .....                         | 9  |
|  | Rollators .....                     | 12 |
|  | Walkers .....                       | 22 |
|  | Wheelchairs .....                   | 28 |
|  | Transport Chairs .....              | 40 |
|  | Bath Safety .....                   | 44 |
|  | Homecare Beds and Accessories ..... | 56 |
|  | Patient Lifts .....                 | 60 |
|  | Mobility Rack .....                 | 63 |

# The Best of Medline's Featured Products



## Faux Teak Shower Benches

A great combination of convenience, functionality, and elegance.

Page 55



## Offset Folding Canes

No more hanging your walking cane on doorknobs, or bending over to pick it up.

Page 11



## All-Terrain Knee Walkers

Get mobile easily and comfortably with this great alternative to crutches.

Page 8



## Upright Rollators

Give those with limited mobility the freedom to walk comfortably upright without leaning over.

Page 18



## Upright Walkers

Take the step toward mobility and maintain upright posture.

Page 24



## Aluminum Commode/Shower Chair

Easy for caregiver to transport patient from bedside to shower.

Page 47



**byACRE Overland Rollators**  
A companion for exploring and adventuring.

Page 19



**Euro Style Shower Chair with Back**  
Shower in style with a new-look shower chair featuring Microban® protection.

Page 49



**Seated Scooters**  
Mobility, convenience and portability you need during leg or foot rehabilitation.

Page 8



**Trigger Walker**  
Lightweight, easy to fold without moving hands from supports, rear glide breaks and two storage pouches.

Page 24



**Transfer Bench Shower Curtain**  
A functional and elegant solution for privacy and safety for tub transfer bench users.

Page 51



**Transfer Bench**  
When issues with pain, strength and flexibility make it difficult to bathe safely, this is the perfect solution. Pairs well with our shower curtain.


Page 51

\*These antimicrobial properties are built-in to protect the product. This product does not protect users or others against bacteria, viruses, germs or other disease organisms.

# Crutches

Crutches are the traditional mobility aid for when one leg is almost entirely non-weightbearing. All axillary (under the arm) models can be adjusted easily to both key measures: overall height and hand height.

1. Select a pair of crutches that's appropriate to the weight and height of the subject.
2. When sizing for axillary crutches, first measure from the ground to near the bottom of the axilla, or armpit. The crutch cushion should be 2" to 3" below the axilla. When standing straight, the cushion should not be touching the axilla. This measurement will be the overall crutch height. (See Dimension A)
3. Keeping arms at your side with a natural bend at the elbow, measure from the ground to the center of the wrist. This measurement will be the hand grip height. (See Dimension B)



## Premium Red-Dot® Axillary Crutches

- Lightweight and durable aluminum, with Red-Dot handgrip for quick and easy sizing

| Item No.  | Wt. Cap. | Size  | Patient Ht. | HCPCS | Pkg.    |
|-----------|----------|-------|-------------|-------|---------|
| G90-214-8 | 300 lb   | Tall  | 5'10"-6'6"  | E0114 | 8 pr/cs |
| G91-214-8 | 300 lb   | Adult | 5'2"-5'10"  | E0114 | 8 pr/cs |
| G92-214-8 | 300 lb   | Youth | 4'6"-5'2"   | E0114 | 8 pr/cs |

## Basic Axillary Crutches

- Economical design made of lightweight and durable aluminum

| Item No.    | Wt. Cap. | Size  | Patient Ht. | HCPCS | Pkg.     |
|-------------|----------|-------|-------------|-------|----------|
| MDS50514-10 | 250 lb   | Tall  | 5'10"-6'6"  | E0114 | 10 pr/cs |
| MDS51514-10 | 250 lb   | Adult | 5'2"-5'10"  | E0114 | 10 pr/cs |
| MDS52514-10 | 250 lb   | Youth | 4'6"-5'2"   | E0114 | 10 pr/cs |

## Standard Axillary Crutches

- Lightweight and durable, with an I-beam to help reduce wobbliness

### Steel

| Item No.  | Wt. Cap. | Size  | Patient Ht. | HCPCS | Pkg.    |
|-----------|----------|-------|-------------|-------|---------|
| MDS80534S | 350 lb   | Tall  | 5'10"-6'6"  | E0114 | 6 pr/cs |
| MDS80535S | 350 lb   | Adult | 5'2"-5'10"  | E0114 | 6 pr/cs |
| MDS80536S | 350 lb   | Youth | 4'6"-5'2"   | E0114 | 6 pr/cs |

### Aluminum

| Item No.  | Wt. Cap. | Size  | Patient Ht. | HCPCS | Pkg.    |
|-----------|----------|-------|-------------|-------|---------|
| MDSV80534 | 300 lb   | Tall  | 5'10"-6'6"  | E0114 | 8 pr/cs |
| MDSV80535 | 300 lb   | Adult | 5'2"-5'10"  | E0114 | 8 pr/cs |
| MDSV80536 | 300 lb   | Youth | 4'6"-5'2"   | E0114 | 8 pr/cs |



G91-214-8

MDSV80535

MDS51514-10

All Crutches in retail packaging.

## Child-Size Axillary Crutches

- Sized for youngsters and made of lightweight and durable aluminum

| Item No.  | Wt. Cap. | Size  | Patient Ht. | HCPCS | Pkg.    |
|-----------|----------|-------|-------------|-------|---------|
| MDS80337Z | 250 lb   | Child | 4'0"-4'6"   | E0114 | 2 pr/cs |
| G53314-8  | 275 lb   | Child | 3'9"-4'7"   | E0114 | 8 pr/cs |

## 3-in-1 Quick-Fit Axillary Crutches

- Adjusts to fit all three sizes: youth, adult and tall adult
- Red-Dot easy hand-grip adjustments
- I-beam height adjustment
- Made of lightweight and durable aluminum

| Item No.   | Wt. Cap. | Size      | Patient Ht. | HCPCS | Pkg.    |
|------------|----------|-----------|-------------|-------|---------|
| MDS80540-8 | 300 lb   | Universal | 4'7"-6'7"   | E0114 | 8 pr/cs |

## Premium Bariatric Axillary Crutches

- Heavy-duty steel for demanding healthcare use

| Item No.   | Wt. Cap. | Size  | Patient Ht. | HCPCS | Pkg. |
|------------|----------|-------|-------------|-------|------|
| MDS80334XW | 650 lb   | Tall  | 5'10"-6'6"  | E0114 | 1 pr |
| MDS80335XW | 650 lb   | Adult | 5'2"-5'10"  | E0114 | 1 pr |

## Standard Bariatric Axillary Crutches

- Heavy-duty steel for demanding healthcare use

| Item No. | Wt. Cap. | Size  | Patient Ht. | HCPCS | Pkg.    |
|----------|----------|-------|-------------|-------|---------|
| G60314B  | 550 lb   | Tall  | 5'10"-6'6"  | E0114 | 2 pr/cs |
| G61314B  | 550 lb   | Adult | 5'2"-5'10"  | E0114 | 2 pr/cs |

## Single-Tube Crutches

- Heavy-duty steel for demanding healthcare use

| Item No. | Wt. Cap. | Size      | Patient Ht. | HCPCS | Pkg.    |
|----------|----------|-----------|-------------|-------|---------|
| MDS80542 | 400 lbs. | Universal | 4'6"-6'6"   | -     | 8 pr/cs |

## Premium Forearm Crutches

- Lightweight and durable aluminum construction

| Item No. | Wt. Cap. | Size  | Patient Ht. | HCPCS | Pkg. |
|----------|----------|-------|-------------|-------|------|
| G05160   | 300 lb   | Tall  | 5'10"-6'6"  | E0110 | 1 pr |
| G05161   | 300 lb   | Adult | 5'0"-6'2"   | E0110 | 1 pr |
| G05162Y  | 300 lb   | Youth | 4'2"-5'2"   | -     | 1 pr |
| G05163C  | 200 lb   | Child | 3'2"-4'6"   | -     | 1 pr |

## Standard Forearm Crutches

- Lightweight and durable aluminum construction

| Item No.  | Wt. Cap. | Size  | Patient Ht. | HCPCS | Pkg. |
|-----------|----------|-------|-------------|-------|------|
| MDS805160 | 250 lb   | Tall  | 5'10"-6'6"  | E0110 | 1 pr |
| MDS805161 | 250 lb   | Adult | 5'0"-6'2"   | E0110 | 1 pr |
| MDS805162 | 250 lb   | Youth | 4'2"-5'2"   | E0110 | 1 pr |

## Crutch Replacement Accessories

| Item No.   | Description                              | HCPCS | Pkg.    |
|------------|--|-------|---------|
| MDS80265W  | Super Crutch Tip, XL, Gray               | A4637 | 6 pr/cs |
| MDS80266RW | Crutch Tip, Gray                         | A4637 | 6 pr/cs |
| G00841     | Super Crutch Tip, 7/8", Black, Guardian® | A4637 | 8 pr/cs |
| G00842     | Super Crutch Tip, 7/8", Gray, Guardian   | A4637 | 8 pr/cs |



# Knee Walkers

Knee walkers are a great alternative to crutches for below-the-knee injuries. We offer a range of mobility aids to fit all heights, shapes and sizes.



## Standard Knee Walkers

- Easy to maneuver with dual hand brakes for superb stopping control
- Can be used with right or left leg

|            | Item No.    | Description      | Wt. Capacity      | Knee pad Height | HCPCS | Retail Pkg. | Pkg. |
|------------|-------------|------------------|-------------------|-----------------|-------|-------------|------|
| <b>New</b> | MDS86000EB  | 8" Wheels, Blue  | 350 lb (158.8 kg) | 14.5"-21"       | E0118 | N           | 1 ea |
| <b>New</b> | MDS86000EJR | 6" Wheels, Red   | 250 lb (113.4 kg) | 14"-17.5"       | E0118 | N           | 1 ea |
|            | MDS86000G4  | 8" Wheels, Black | 300 lb (136.1 kg) | 17"-23"         | E0118 | N           | 1 ea |
| <b>New</b> | MDS86000SS  | Seated Scooter   | 300 lb (136.1 kg) | 21"-28.5"       | -     | N           | 1 ea |

## All-Terrain Knee Walkers

|            | Item No. | Description       | Wt. Capacity      | Knee pad Height | HCPCS | Retail Pkg. | Pkg. |
|------------|----------|-------------------|-------------------|-----------------|-------|-------------|------|
|            | ATKNEE   | 12" Wheels, Black | 350 lb (158.8 kg) | 20"-26"         | E0118 | N           | 1 ea |
| <b>New</b> | ATKNEEJR | 10" Wheels, Black | 350 lb (158.8 kg) | 17.25"-24"      | E0118 | N           | 1 ea |



ATKNEE




MDS86000SS



# Canes

Standard canes with a single tip should be used for those looking for help with balance. Quad canes are useful for those looking to bear weight as they have a wider base. When sizing a cane, follow these steps:

1. Select a cane that's appropriate to the weight and height of the user.
2. With arms at side and a natural bend at the elbow, measure from the ground to the center of the wrist. This measurement will determine the ideal cane height.
3. Loosen the collar on the cane and press the push-button lock, slide the lower tube into or out of the body of the cane until you get to the wrist-height measurement, and let the push button "click" into the closest hole to the measurement. Tighten the collar (check and retighten the collar regularly).



## Premium Offset Handle Canes

- Offset handle design centers the user's weight over the strongest part of the cane
- Ergonomically designed to provide comfort and support for partial weight-bearing users
- Push-button, height-adjustable to find the perfect height
- Secured with a locking ring to provide quiet, rattle-free use and added safety

| Item No. | Description                              | Wt. Capacity      | Handle Material | Height Range | HCPCS | Retail Pkg. | Pkg. |
|----------|--|-------------------|-----------------|--------------|-------|-------------|------|
| G05355   | Offset Handle Cane, Vinyl Handle, Chrome | 300 lb (136.1 kg) | Vinyl           | 31"-40"      | E0100 | Y           | 6/cs |
| G05356   | Offset Handle Cane, Vinyl Handle, Black  | 300 lb (136.1 kg) | Vinyl           | 31"-40"      | E0100 | Y           | 6/cs |

## Offset Canes

- Single-point canes
- With foam handles and wrist straps

| Item No.     | Description | Wt. Capacity      | Handle Material | Height Range | HCPCS | Retail Pkg. | Pkg. |
|--------------|-------------|-------------------|-----------------|--------------|-------|-------------|------|
| MDS86420     | Black       | 300 lb (136.1 kg) | Foam            | 29"-38"      | E0100 | Y           | 6/cs |
| MDS86420BRZW | Bronze      | 300 lb (136.1 kg) | Foam            | 29"-38"      | E0100 | Y           | 6/cs |
| MDS86420CHRW | Chrome      | 300 lb (136.1 kg) | Foam            | 29"-38"      | E0100 | Y           | 6/cs |



G05355



MDS86420



MDS86420BRZW

## New Gel Grip Offset Canes

- Our cost-effective fashion canes are mobility aids and fashion accessories
- Colorful options provide users with the right cane for different occasions
- Comfortable gel handles and convenient wrist straps

| Item No.     | Description      | Wt. Capacity      | Handle Material | Height Range | HCPCS | Retail Pkg. | Pkg. |
|--------------|------------------|-------------------|-----------------|--------------|-------|-------------|------|
| MDS86440BLH  | Gel Grip, Blue   | 300 lb (136.1 kg) | Gel             | 29"-38"      | E0100 | Y           | 1 ea |
| MDS86440PRH  | Gel Grip, Purple | 300 lb (136.1 kg) | Gel             | 29"-38"      | E0100 | Y           | 1 ea |
| MDS86440CHRH | Gel Grip, Chrome | 300 lb (136.1 kg) | Gel             | 29"-38"      | E0100 | Y           | 1 ea |
| MDS86440BKH  | Gel Grip, Black  | 300 lb (136.1 kg) | Gel             | 29"-38"      | E0100 | Y           | 1 ea |

## Quad Canes

- Provide a secure grip with four feet
- Allow greater weight distribution for stability than single-point canes

| Item No.    | Description        | Wt. Capacity      | Handle Material | Height Range | HCPCS | Retail Pkg. | Pkg. |
|-------------|--------------------|-------------------|-----------------|--------------|-------|-------------|------|
| MDS86222W   | Small Base, Black  | 300 lb (136.1 kg) | Foam            | 29"-38"      | E0105 | Y           | 2/cs |
| MDS86222CHR | Small Base, Chrome | 300 lb (136.1 kg) | Foam            | 29"-38"      | E0105 | Y           | 2/cs |
| G05345S     | Small Base, Chrome | 300 lb (136.1 kg) | Vinyl           | 29"-38"      | E0105 | Y           | 2/cs |
| MDS86228    | Large Base, Black  | 300 lb (136.1 kg) | Foam            | 29"-38"      | E0105 | Y           | 2/cs |
| MDS86228CHR | Large Base, Chrome | 300 lb (136.1 kg) | Foam            | 29"-38"      | E0105 | Y           | 2/cs |

## Folding T-Handle Cane

- Folding capabilities make this cane perfect for travel
- Fits in a suitcase or carry-on luggage

| Item No.  | Description | Wt. Capacity      | Handle Material | Height Range | HCPCS | Retail Pkg. | Pkg. |
|-----------|-------------|-------------------|-----------------|--------------|-------|-------------|------|
| MDS86430W | Black       | 250 lb (113.4 kg) | Vinyl           | 32"-36"      | E0100 | Y           | 6/cs |



MDS86440BLH



MDS86440CHRH



MDS86222W



MDS86430W

New  
Gel  
Grips

Folds down  
for easy  
transport  
and storage



## Standard Aluminum Canes

- Standard handle lets users hook the cane over their arm while opening doors or performing other tasks
- Foam grip provides a comfortable secure grip

| Item No.    | Description | Wt. Capacity      | Handle Material | Height Range | HCPCS | Retail Pkg. | Pkg. |
|-------------|-------------|-------------------|-----------------|--------------|-------|-------------|------|
| MDS86415    | Black       | 300 lb (136.1 kg) | Foam            | 30"-38"      | E0100 | Y           | 6/cs |
| MDS86415CHR | Chrome      | 300 lb (136.1 kg) | Foam            | 30"-38"      | E0100 | Y           | 6/cs |
| G05370      | Chrome      | 300 lb (136.1 kg) | Vinyl           | 30"-39"      | E0100 | Y           | 6/cs |

## Bariatric Canes

- All-steel construction with dual-tubing insert for added strength and security

| Item No.                 | Description                     | Wt. Capacity      | Handle Material | Height Range | HCPCS | Retail Pkg. | Pkg. |
|--------------------------|---------------------------------|-------------------|-----------------|--------------|-------|-------------|------|
| MDS86420XW               | Offset Handle, Black            | 500 lb (226.8 kg) | Vinyl           | 29"-38"      | E0100 | Y           | 1 ea |
| <b>New</b> MDS86420XWCHR | Offset Single Point, Chrome     | 500 lb (226.8 kg) | Vinyl           | 29"-38"      | E0100 | Y           | 1 ea |
| <b>New</b> MDS86420XWTL  | Offset Handle Cane, Tall, Black | 500 lb (226.8 kg) | Vinyl           | 37"-48"      | E0100 | Y           | 1 ea |
| MDS86222XW               | Small Base Quad, Black          | 500 lb (226.8 kg) | Vinyl           | 29"-38"      | E0105 | Y           | 1 ea |
| MDS86228XW               | Large Base Quad, Black          | 500 lb (226.8 kg) | Vinyl           | 29"-38"      | E0105 | Y           | 1 ea |

## **New** Hybrid Canes

- Gel grip for user comfort
- Folding capabilities make this cane perfect for travel
- Four-point tip allows cane to be self-standing

| Item No. | Description     | Wt. Capacity      | Handle Material | Height Range | HCPCS | Retail Pkg. | Pkg. |
|----------|-----------------|-------------------|-----------------|--------------|-------|-------------|------|
| 4UCANE   | Aluminum, Black | 350 lb (158.8 kg) | Gel Grip        | 32"-36"      | E0100 | Y           | 1 ea |
| 4UCANEB  | Aluminum, Blue  | 350 lb (158.8 kg) | Gel Grip        | 32"-36"      | E0100 | Y           | 1 ea |

## **New** TriGlide™ Canes

| Item No. | Description | Wt. Capacity      | Handle Material | Height Range | HCPCS | Retail Pkg. | Pkg. |
|----------|-------------|-------------------|-----------------|--------------|-------|-------------|------|
| TGLCANE  | Black       | 350 lb (158.8 kg) | Gel             | 32"-37"      | -     | Y           | 1 ea |
| TGLCANEB | Blue        | 350 lb (158.8 kg) | Gel             | 32"-37"      | -     | Y           | 1 ea |



TGLCANE



MDS86420XW



MDS86228XW



4UCANE

## Canes Accessories

| Item No.   | Description        | Retail Pkg. | Pkg. |
|------------|--------------------|-------------|------|
| MDS86425W  | Quad Cane Tips     | Y           | 6/cs |
| MDS86426W  | Standard Cane Tips | Y           | 6/cs |
| 4PTCANETIP | 4-Point Cane Tip   | Y           | 4/cs |



4UCANEB


# Rollators

Versatility meets personality with our superior group of rollators for the customer looking for more. From essential basics to premier collections, Medline takes your customer's needs and lifestyle into consideration. Engineered with design features that solve specific problems for the user in unique colors and additional accessories. Discover both functional and luxury in a product category that historically only focused on utility.

Rollators provide users with both a mobility aid and a seat to rest on when tired. Most rollators allow handle adjustment for a proper fit and comfort.

To accurately size the rollator handle height keep arms at your side with a natural bend at the elbow. Align your wrist to the handle as shown in the image.

Follow your rollator's manual on how to adjust the handle as rollator handle height adjustment mechanisms differ from one rollator to another.



## Basic Steel Rollators

- Durable powder-coated steel frame
- Storage bag under the seat is included
- Folds for easy storage/transport

| Item No.      | Description                      | Wt. Capacity      | Seat Height | Wheel Size | Height Range | Hand Grip Width | HCPCS       | Retail Pkg. | Pkg. |
|---------------|----------------------------------|-------------------|-------------|------------|--------------|-----------------|-------------|-------------|------|
| MDS86850ESKD  | Burgundy                         | 350 lb (158.8 kg) | 20.5"       | 6"         | 31"-35"      | 17"             | E0143+E0156 | Y           | 1 ea |
| MDS86850ES    | Burgundy                         | 350 lb (158.8 kg) | 20.5"       | 6"         | 31"-35"      | 17.5"           | E0143+E0156 | N           | 1 ea |
| MDS86850EBS   | Blue                             | 350 lb (158.8 kg) | 20.5"       | 6"         | 31"-35"      | 17.5"           | E0143+E0156 | N           | 1 ea |
| MDS86830EBLS  | Black                            | 300 lb (136.1 kg) | 20.5"       | 6"         | 31"-35"      | 17.5"           | E0143+E0156 | N           | 1 ea |
| MDS86850EGS   | Green                            | 300 lb (136.1 kg) | 20.5"       | 6"         | 31"-35"      | 17.5"           | E0143+E0156 | N           | 1 ea |
| MDS86840EBS   | Light Blue                       | 300 lb (136.1 kg) | 20.5"       | 6"         | 31"-35"      | 19.75"          | E0143+E0156 | N           | 1 ea |
| MDS86850ES8   | Burgundy                         | 350 lb (158.8 kg) | 20.5"       | 8"         | 31"-35"      | 17.5"           | E0143+E0156 | N           | 1 ea |
| MDS86850EBS8  | Blue                             | 350 lb (158.8 kg) | 20.5"       | 8"         | 31"-35"      | 17.5"           | E0143+E0156 | N           | 1 ea |
| MDS86830EBLS8 | Black                            | 300 lb (136.1 kg) | 20.5"       | 8"         | 31"-35"      | 17.5"           | E0143+E0156 | N           | 1 ea |
| MDS86840EBS8  | Light Blue                       | 300 lb (136.1 kg) | 20.5"       | 8"         | 31"-35"      | 19.75"          | E0143+E0156 | N           | 1 ea |
| MDS86850EGS8  | Green                            | 300 lb (136.1 kg) | 20.5"       | 8"         | 31"-35"      | 17.5"           | E0143+E0156 | N           | 1 ea |
| MDS86850BAG   | Black Bag for MDS86850 Rollators |                   |             |            |              |                 |             | N           | 1 ea |



MDS86840EBS



MDS86850ES8

## Basic Steel Rollators with 8" Wheels and Basket

- Easily removable wire basket
- Flip up backrest for easy storage

| Item No.     | Description | Wt. Capacity      | Seat Height | Wheel Size | Height Range | Hand Grip Width | HCPCS | Retail Pkg. | Pkg. |
|--------------|-------------|-------------------|-------------|------------|--------------|-----------------|-------|-------------|------|
| MDS86860ERS8 | Red         | 300 lb (136.1 kg) | 21.25"      | 8"         | 31"-35"      | 18.5"           | -     | N           | 1 ea |
| MDS86860EBS8 | Blue        | 300 lb (136.1 kg) | 21.25"      | 8"         | 31"-35"      | 18.5"           | -     | N           | 1 ea |

## Basic Steel Rollators with 6" Wheels and Microban\*

- Antimicrobial\* protection on handles and seats to protect the product

| Item No.     | Description | Wt. Capacity      | Seat Height | Wheel Size | Height Range | Hand Grip Width | HCPCS | Retail Pkg. | Pkg. |
|--------------|-------------|-------------------|-------------|------------|--------------|-----------------|-------|-------------|------|
| MDS86850ESAB | Teal        | 300 lb (136.1 kg) | 20"         | 6"         | 30.25"-34.5" | 17.5"           | -     | N           | 1 ea |
| MDS86850ESAR | Burgundy    | 300 lb (136.1 kg) | 20"         | 6"         | 30.25"-34.5" | 17.5"           | -     | N           | 1 ea |

\*These antimicrobial properties are built-in to protect the product. This product does not protect users or others against bacteria, viruses, germs or other disease organisms.

## Basic Aluminum Rollators

- Lightweight durable aluminum frame
- Storage bag under the seat is included
- Folds for easy storage/transport

| Item No.     | Description | Wt. Capacity      | Seat Height | Wheel Size | Height Range | Hand Grip Width | HCPCS       | Retail Pkg. | Pkg. |
|--------------|-------------|-------------------|-------------|------------|--------------|-----------------|-------------|-------------|------|
| MDS86850EB   | Blue        | 250 lb (113.4 kg) | 20.5"       | 6"         | 31"-35"      | 17.5"           | E0143+E0156 | N           | 1 ea |
| MDS86850E    | Burgundy    | 250 lb (113.4 kg) | 20.5"       | 6"         | 31"-35"      | 17.5"           | E0143+E0156 | Y           | 1 ea |
| MDS86850EP   | Purple      | 250 lb (113.4 kg) | 20.5"       | 6"         | 31"-35"      | 17.5"           | E0143+E0156 | N           | 1 ea |
| MDS86850EG   | Green       | 250 lb (113.4 kg) | 20.5"       | 6"         | 31"-35"      | 17.5"           | E0143+E0156 | N           | 1 ea |
| MDS86830EBL  | Black       | 300 lb (136.1 kg) | 20.5"       | 6"         | 31"-35"      | 17.5"           | E0143+E0156 | N           | 1 ea |
| MDS86830EBL8 | Black       | 300 lb (136.1 kg) | 20.5"       | 8"         | 31.5"-37.5"  | 17.5"           | E0143+E0156 | N           | 1 ea |

## Basic Jr. Rollators with 6" Wheels

| Item No.     | Description | Wt. Capacity      | Seat Height | Wheel Size | Height Range | Hand Grip Width | HCPCS       | Retail Pkg. | Pkg. |
|--------------|-------------|-------------------|-------------|------------|--------------|-----------------|-------------|-------------|------|
| MDS86850EJR  | Burgundy    | 250 lb (113.4 kg) | 20.5"       | 6"         | 29"-33"      | 17.5"           | E0143+E0156 | N           | 1 ea |
| MDS86850EBJR | Blue        | 250 lb (113.4 kg) | 20.5"       | 6"         | 29"-33"      | 17.5"           | E0143+E0156 | N           | 1 ea |



### Superlight Rollators with 6" Wheels

- Lightweight aluminum frame
- Weighs 12.6 lbs. making it easy to transport

| Item No.      | Description | Wt. Capacity      | Seat Height | Height Range | Hand Grip Width | HCPCS       | Retail Pkg. | Pkg. |
|---------------|-------------|-------------------|-------------|--------------|-----------------|-------------|-------------|------|
| MDS86825SLR   | Red         | 250 lb (113.4 kg) | 18"-23"     | 29.5"-36"    | 17"             | E0143+E0156 | Y           | 1 ea |
| MDS86825SLBLK | Black       | 250 lb (113.4 kg) | 18"-23"     | 29.5"-36"    | 17"             | E0143+E0156 | Y           | 1 ea |
| MDS86825SLB   | Smoky Blue  | 250 lb (113.4 kg) | 18"-23"     | 29.5"-36"    | 17"             | E0143+E0156 | Y           | 1 ea |

### Universal Steel Rollator with 6" Wheels

- Durable steel frame
- Foldable for easy transport or store when not in use

| Item No.     | Description | Wt. Capacity      | Seat Height | Height Range | Hand Grip Width | HCPCS | Retail Pkg. | Pkg. |
|--------------|-------------|-------------------|-------------|--------------|-----------------|-------|-------------|------|
| MDS86825SLRS | Red         | 350 lb (158.8 kg) | 18"-23"     | 29.5"-36"    | 17"             | -     | Y           | 1 ea |

### Luxe Wide Rollator with 7" Wheels

- Lightweight aluminum frame
- Wider seat
- Height adjustable handle and seat

| Item No.  | Description | Wt. Capacity      | Seat Height | Height Range | Hand Grip Width | HCPCS         | Retail Pkg. | Pkg. |
|-----------|-------------|-------------------|-------------|--------------|-----------------|---------------|-------------|------|
| MDS86835W | White       | 300 lb (136.1 kg) | 19.5"-24"   | 29"-33"      | 19"             | E0143 + E0156 | N           | 1 ea |



Superlight options weigh only 12.6 lbs!

MDS86825SLB



Wider seat

MDS86835W

## Bariatric Rollators with 8" Wheels

- Heavy-duty bariatric rollators
- Wider distance between handles for users needing extra room

| Item No.      | Description                    | Wt. Capacity      | Seat Height | Height Range | Hand Grip Width | HCPCS       | Retail Pkg. | Pkg. |
|---------------|--------------------------------|-------------------|-------------|--------------|-----------------|-------------|-------------|------|
| MDS86800XWS   | Steel, Burgundy                | 500 lb (226.8 kg) | 22"         | 31.5"–37.5"  | 23.5"           | E0149+E0156 | N           | 1 ea |
| MDS86800XW    | Aluminum, Burgundy             | 400 lb (181.4 kg) | 22"         | 31.5"–37.5"  | 23.5"           | E0149+E0156 | N           | 1 ea |
| MDS86800XWBSK | Basket for MDS86800XW Rollator |                   |             |              |                 | –           | N           | 1 ea |

## New Premium Bariatric Rollator with 8" Wheels

- Premium wheels for added maneuverability

| Item No.  | Description     | Wt. Capacity      | Seat Height | Height Range | Hand Grip Width | HCPCS | Retail Pkg. | Pkg. |
|-----------|-----------------|-------------------|-------------|--------------|-----------------|-------|-------------|------|
| MDS868XWR | Steel, Burgundy | 500 lb (226.8 kg) | 22"         | 31.5"–37.25" | 23.5"           | –     | N           | 1 ea |

## Deluxe Comfort Rollator with 8" Wheels

- Lightweight frame
- Memory foam seat
- Under storage bag and front pouch
- Cup and cane holder

| Item No.   | Description          | Wt. Capacity      | Seat Height | Height Range | Hand Grip Width | HCPCS | Retail Pkg. | Pkg. |
|------------|----------------------|-------------------|-------------|--------------|-----------------|-------|-------------|------|
| MDS86826SC | Aluminum, Smoky Teal | 300 lb (136.1 kg) | 23"         | 32.5"–37.5"  | 18.5"           | –     | N           | 1 ea |



Durable premium wheels for various terrain

MDS868XWR



Cane and cup holder included

MDS86826SC

### Guardian Rollators with 8" Wheels

- Aluminum frame
- Padded seat and curved padded backrest provide enhanced support and comfort

| Item No. | Description | Wt. Capacity      | Seat Height | Height Range | Hand Grip Width | HCPCS       | Retail Pkg. | Pkg. |
|----------|-------------|-------------------|-------------|--------------|-----------------|-------------|-------------|------|
| G07887B  | Blue        | 300 lb (136.1 kg) | 20"         | 31"-36.5"    | 17.5"           | E0143+E0156 | N           | 1 ea |
| G07887R  | Rose        | 300 lb (136.1 kg) | 20"         | 31"-36.5"    | 17.5"           | E0143+E0156 | N           | 1 ea |

### Momentum Rollators with 6" Wheels

- Sporty aesthetic with durable steel frame with large padded seat
- Seat flips up to enable user to walk inside
- All-terrain wheels make it easier to get around
- Product weighs 18.8 lbs.

| Item No.  | Description | Wt. Capacity      | Seat Height | Height Range | Hand Grip Width | HCPCS | Retail Pkg. | Pkg. |
|-----------|-------------|-------------------|-------------|--------------|-----------------|-------|-------------|------|
| MDS86870R | Red         | 250 lb (113.4 kg) | 20.5"       | 34.25"-41.5" | 18"             | -     | N           | 1 ea |
| MDS86870W | White       | 250 lb (113.4 kg) | 20.5"       | 34.25"-41.5" | 18"             | -     | N           | 1 ea |
| MDS86870G | Grey        | 250 lb (113.4 kg) | 20.5"       | 34.25"-41.5" | 18"             | -     | N           | 1 ea |

Seat flips up to allow user to walk inside the rollator



G07887B



MDS86870G



## Empower Rollators with 8" Wheels

- Aluminum frame with memory foam seat
- Padded backrest flips forward and backward enabling the user to sit comfortably on either side
- Seat flips up to enable user to walk inside

| Item No.     | Description          | Wt. Capacity       | Seat Height | Height Range | Hand Grip Width | HCPCS | Retail Pkg. | Pkg. |
|--------------|----------------------|--------------------|-------------|--------------|-----------------|-------|-------------|------|
| MDS86845BM   | Blue                 | 300 lb (136.1 kg)  | 20.5"-24.5" | 34.25"-41.5" | 18"             | -     | Y           | 1 ea |
| MDS86845BLKM | Black                | 300 lb (136.1 kg)  | 20.5"-24.5" | 34.25"-41.5" | 18"             | -     | Y           | 1 ea |
| MDS86845SM   | Silver               | 300 lb. (136.1 kg) | 20.5"-24.5" | 34.25"-41.5" | 18"             | -     | Y           | 1 ea |
| MDS86845BK   | Black with Microban* | 300 lb (136.1 kg)  | 20.5"-24.5" | 34.25"-41.5" | 18"             | -     | Y           | 1 ea |
| MDS86845RD   | Red with Microban*   | 300 lb (136.1 kg)  | 20.5"-24.5" | 34.25"-41.5" | 18"             | -     | Y           | 1 ea |

\*These antimicrobial properties are built-in to protect the product. This product does not protect users or others against bacteria, viruses, germs or other disease organisms.

**Weighs only 17 lbs**

Convenient loop-style, easy-squeeze brakes

Padded elbow rests

Spacious memory foam seat is 35% larger than most rollator seats

Adjustable seat height

Roomy neoprene bag for storage

Large, stable 8" wheels for support and maneuverability

Wheels nest to fold flatter

MDS86845BM

Handle height adjusts for the perfect fit

Handles fold in

Backrest flips for comfortable seating in either direction

Seat flips up for better stability and accessibility

### Three-Wheel Rollators

- Easy to maneuver with a tray, basket and storage bag

| Item No.     | Description     | Wt. Capacity      | Seat Height | Height Range | Hand Grip Width | HCPCS   | Retail Pkg. | Pkg. |
|--------------|-----------------|-------------------|-------------|--------------|-----------------|---------|-------------|------|
| MDS86840TRIB | Blue, 8" Wheels | 275 lb (124.7 kg) | 20.5"       | 31.75"–35.5" | 14.5"           | Pending | N           | 1 ea |
| MDS86840TRIR | Red, 8" wheels  | 275 lb (124.7 kg) | 20.5"       | 31.75"–35.5" | 14.5"           | Pending | N           | 1 ea |



MDS86840TRIB

### New Upright Rollator

- Easy maneuverability and helps to promote good posture

| Item No.   | Description | Wt. Capacity      | Seat Height | Height Range | Hand Grip Width | HCPCS | Retail Pkg. | Pkg. |
|------------|-------------|-------------------|-------------|--------------|-----------------|-------|-------------|------|
| MDSUPROLLW | White       | 300 lb (136.1 kg) | 21.5"       | 41"–46"      | 14.5"           | –     | N           | 1 ea |



MDSUPROLLW

### European Style Rollators

- European-style rollators feature a lightweight aluminum frame
- Fold horizontally for easy storage
- Storage bag included

| Item No.               | Description            | Wt. Capacity       | Seat Height | Height Range | Color     | HCPCS | Retail Pkg. | Pkg. |
|------------------------|------------------------|--------------------|-------------|--------------|-----------|-------|-------------|------|
| MDS86890EUB            | Simplicity 2 Rollator  | 300 lb. (136.1 kg) | 19"         | 32"–36"      | Blue      | –     | N           | 1 ea |
| MDS86890EUS            | Simplicity 2 Rollator  | 300 lb (136.1 kg)  | 19"         | 32"–36"      | Silver    | –     | N           | 1 ea |
| <b>New</b> MDS86870ANB | Lightweight, Knockdown | 300 lb (136.1 kg)  | 21"         | 34"–38"      | Navy Blue | –     | N           | 1 ea |
| <b>New</b> MDS86870SHW | With Suspension Shocks | 300 lb (136.1 kg)  | 22"         | 35"–39"      | White     | –     | N           | 1 ea |
| <b>New</b> MDS86880R   | Euro-Style             | 300 lb (136.1 kg)  | 21"         | 30"–35"      | Red       | –     | N           | 1 ea |



MDS86890EUB



MDS86870ANB



MDS86870SHW



MDS86880R

**New byACRE Overland All-Terrain Carbon Fiber Rollator**

- Durable pneumatic tires and superb suspension
- Weighs 14.78 lbs.
- Carbon fiber frame is 5x stronger than steel
- Ergonomic handles have integrated brake lines and are height adjustable in 1' increments
- Includes organizer bag

| Item No. | Description    | Wt. Capacity      | Seat Height | Height Range | Seat Width | HCPCS | Retail Pkg. | Pkg. |
|----------|----------------|-------------------|-------------|--------------|------------|-------|-------------|------|
| BYA1500G | Defender Green | 330 lb (149.7 kg) | 24"         | 32.3"-37.4   | 18.5"      | -     | N           | 1 ea |



**New byACRE Overland All-Terrain Carbon Fiber Rollator Accessories**

| Item No. | Description        | Retail Pkg. | Pkg. |
|----------|--------------------|-------------|------|
| BYAOPUMP | Mini Pump          | N           | 1 ea |
| BYAOBR   | Backrest           | N           | 1 ea |
| BYAOTVL  | Travel Bag         | N           | 1 ea |
| BYAODP   | Daypack            | N           | 1 ea |
| BYAOPB   | Pocket Bag         | N           | 1 ea |
| BYAOCNHR | Cane Holder, Right | N           | 1 ea |
| BYAOCNHL | Cane Holder, Left  | N           | 1 ea |



## byACRE Ultralight Carbon Fiber Rollators

- World's lightest rollators at 10.6 lbs.
- Carbon fiber frame is 5x stronger than steel
- Ergonomic handles have integrated brake lines and are height adjustable in 1" increments
- All byACRE rollators include organizer bag

World's  
lightest  
rollator

### Compact

| Item No.   | Description | Wt. Capacity      | Seat Height | Height Range | Seat Width | HCPCS | Retail Pkg. | Pkg. |
|------------|-------------|-------------------|-------------|--------------|------------|-------|-------------|------|
| BYA100SMBK | Black       | 285 lb (129.3 kg) | 22"         | 31"-36"      | 16"        | -     | N           | 1 ea |
| BYA100SMR  | Red         | 285 lb (129.3 kg) | 22"         | 31"-36"      | 16"        | -     | N           | 1 ea |
| BYA100SMW  | White       | 285 lb (129.3 kg) | 22"         | 31"-36"      | 16"        | -     | N           | 1 ea |

### Regular

| Item No.   | Description | Wt. Capacity      | Seat Height | Height Range | Seat Width | HCPCS | Retail Pkg. | Pkg. |
|------------|-------------|-------------------|-------------|--------------|------------|-------|-------------|------|
| BYA100MDBK | Black       | 285 lb (129.3 kg) | 24"         | 31.5"-37.4"  | 16.5"      | -     | N           | 1 ea |
| BYA100MDR  | Red         | 285 lb (129.3 kg) | 24"         | 31.5"-37.4"  | 16.5"      | -     | N           | 1 ea |
| BYA100MDW  | White       | 285 lb (129.3 kg) | 24"         | 31.5"-37.4"  | 16.5"      | -     | N           | 1 ea |

### Comfort

| Item No.   | Description | Wt. Capacity      | Seat Height | Height Range | Seat Width | HCPCS | Retail Pkg. | Pkg. |
|------------|-------------|-------------------|-------------|--------------|------------|-------|-------------|------|
| BYA100LGBK | Black       | 285 lb (129.3 kg) | 24"         | 31.5"-37.4"  | 18.5"      | -     | N           | 1 ea |
| BYA100LGR  | Red         | 285 lb (129.3 kg) | 24"         | 31.5"-37.4"  | 18.5"      | -     | N           | 1 ea |
| BYA100LGW  | White       | 285 lb (129.3 kg) | 24"         | 31.5"-37.4"  | 18.5"      | -     | N           | 1 ea |



BYA100SMBK



BYA100MDW



BYA100LGR

Folds for easy  
portability



## byACRE Ultralight Carbon Fiber Rollator Accessories

| Item No. | Description                       | Retail Pkg. | Pkg. |
|----------|-----------------------------------|-------------|------|
| BYAGRBG  | Grocery Bag, Grey                 | N           | 1 ea |
| BYAWKNBG | Weekend Bag, Grey                 | N           | 1 ea |
| BYATVLBG | Travel Bag, Grey                  | N           | 1 ea |
| BYACNHR  | Cane Holder, Right                | N           | 1 ea |
| BYACNHL  | Cane Holder, Left                 | N           | 1 ea |
| BYAOBGR  | Organizer Bag, Small/Medium, Grey | N           | 1 ea |
| BYAOBGL  | Organizer Bag, Large, Grey        | N           | 1 ea |
| BYASMBR  | Backrest, Small/Medium            | N           | 1 ea |
| BYALBR   | Backrest, Large                   | N           | 1 ea |



BYAGRBG



BYAWKNBG



BYATVLBG



BYAOBGR




BYASMBR

# Walkers

Performance and style are found in our superior group of walkers, engineered with design features that solve specific mobility issues.

Proper walker sizing can help prevent poor posture and injury from misuse or falling. Medline's walkers use push-button mechanisms, so they can be adjusted to various heights.

1. Select a walker that's appropriate to the weight of the patient.
2. Keeping arms at your side with a natural bend at the elbow, measure from the ground to the center of the wrist.
3. Adjust all four legs to the same height so that the handgrips match the same height as the wrist.



## Hemi Walkers

- Lightweight design allows the hemi sidestepper to be maneuvered easily—even with one arm
- Bi-level hand grip offers extra assistance

| Item No. | Description                 | Wt. Capacity      | Handle Material | Height Range | Hand Grip Width | HCPCS | Retail Pkg. | Pkg. |
|----------|-----------------------------|-------------------|-----------------|--------------|-----------------|-------|-------------|------|
| G07770X1 | Adult Hemi Walker           | 250 lb (113.4 kg) | PVC             | 32"–36"      | N/A             | E0135 | N           | 1 ea |
| G07770XW | Bariatric Adult Hemi Walker | 500 lb (226.8 kg) | PVC             | 32"–36"      | N/A             | E0148 | N           | 1 ea |
| G07771   | Youth Hemi Walker           | 250 lb (113.4 kg) | PVC             | 28.5"–32"    | N/A             | E0135 | N           | 2/cs |

## Two-button Folding Walkers\*

- Durable walker with two-button folding system
- Push-button, height-adjustable legs

| Item No.     | Description         | Wt. Capacity      | Handle Material | Height Range | Hand Grip Width | HCPCS          | Retail Pkg. | Pkg. |
|--------------|---------------------|-------------------|-----------------|--------------|-----------------|----------------|-------------|------|
| MDS864104    | Aluminum, No Wheels | 300 lb (136.1 kg) | PVC             | 32"–39"      | 18"             | E0135          | N           | 4/cs |
| MDS864104B   | Steel, No Wheels    | 400 lb (181.4 kg) | PVC             | 32"–39"      | 18"             | E0135          | N           | 4/cs |
| MDS86410W4   | Aluminum, 3" Wheels | 300 lb (136.1 kg) | PVC             | 32"–38"      | 17.5"           | E0143          | N           | 4/cs |
| MDS86410W54  | Aluminum, 5" Wheels | 300 lb (136.1 kg) | PVC             | 32"–38"      | 17.5"           | E0143          | N           | 4/cs |
| MDS86410W54B | Steel, 5" Wheels    | 400 lb (181.4 kg) | PVC             | 32.5"–39.5"  | 17.5"           | E0143 or E0149 | N           | 4/cs |



G07770X1



MDS864104



MDS86410W54B

\*All walkers come unassembled

### Youth Two-Button Folding Walkers\*

- Youth size durable walker with two-button folding system
- Comfortable grip with nonmarking tips



MDS86410JW54

| Item No.      | Description         | Wt. Capacity      | Handle Material | Height Range | Hand Grip Width | HCPCS | Retail Pkg. | Pkg. |
|---------------|---------------------|-------------------|-----------------|--------------|-----------------|-------|-------------|------|
| MDS86410J4    | Aluminum, No Wheels | 300 lb (136.1 kg) | PVC             | 25"-32"      | 17.5"           | E0135 | N           | 4/cs |
| MDS86410J4B   | Steel, No Wheels    | 400 lb (181.4 kg) | PVC             | 25"-32"      | 17.5"           | -     | N           | 4/cs |
| MDS86410JW54  | Aluminum, 5" Wheels | 300 lb (136.1 kg) | PVC             | 25"-32"      | 17.5"           | E0143 | N           | 4/cs |
| MDS86410JW54B | Steel, 5" Wheels    | 400 lb (181.4 kg) | PVC             | 25"-32"      | 17.5"           | E0149 | N           | 4/cs |

### Bariatric Folding Walkers\*

- Extra-wide frames and higher weight capacities provide a comfortable fit for larger patients
- Includes two-button folding system
- All bariatric folding walkers come assembled



MDS86410XWW

| Item No.     | Description         | Wt. Capacity      | Handle Material | Height Range | Hand Grip Width | HCPCS | Retail Pkg. | Pkg. |
|--------------|---------------------|-------------------|-----------------|--------------|-----------------|-------|-------------|------|
| MDS86410XW   | Aluminum, No Wheels | 500 lb (226.8 kg) | PVC             | 31"-41"      | 20.5"           | E0148 | N           | 1 ea |
| MDS86410XWB  | Steel, No Wheels    | 650 lb (294.8 kg) | PVC             | 31"-41"      | 20"             | -     | N           | 1 ea |
| MDS86410XWW  | Aluminum, 5" Wheels | 500 lb (226.8 kg) | PVC             | 33.5"-43.25" | 20.5"           | E0149 | N           | 1 ea |
| MDS86410XWWB | Steel, 5" Wheels    | 600 lb (272.2 kg) | PVC             | 31"-41"      | 20"             | E0149 | N           | 1 ea |

### Folding Paddle Walkers\*

- Ideal design for patients with manual dexterity limitations
- The paddle mechanism can be flipped-up or pushed-down to fold the walker



MDS86410P

| Item No.     | Description               | Wt. Capacity      | Handle Material | Height Range | Hand Grip Width | HCPCS | Retail Pkg. | Pkg. |
|--------------|---------------------------|-------------------|-----------------|--------------|-----------------|-------|-------------|------|
| MDS86410P    | Aluminum, No Wheels       | 300 lb (136.1 kg) | PVC             | 32"-39"      | 16.5"           | E0135 | N           | 4/cs |
| MDS86410P5   | Aluminum, 5" Wheels       | 300 lb (136.1 kg) | PVC             | 32"-39"      | 16.5"           | E0143 | N           | 4/cs |
| MDS86410P5J  | Junior, 5" Wheels         | 300 lb (136.1 kg) | PVC             | 25.5"-32.5"  | 16.5"           | E0143 | N           | 4/cs |
| MDS86410KDBW | Blue, Aluminum, 5" Wheels | 350 lb (158.8 kg) | PVC             | 32.5"-39.5"  | 16.75"          | E0149 | Y           | 1 ea |

\*All walkers come unassembled

## Upright Folding Walkers\*

- Upright folding walkers help users with posture compared to standard walkers
- Comfortable grip with mark-resistant tips

| Item No.    | Description               | Wt. Capacity      | Handle Material | Height Range | Hand Grip Width | HCPCS | Retail Pkg. | Pkg. |
|-------------|---------------------------|-------------------|-----------------|--------------|-----------------|-------|-------------|------|
| MDSUPWALK   | Aluminum, 5" Wheels       | 300 lb (136.1 kg) | PVC             | 46"-53"      | 20"             | -     | N           | 1 ea |
| MDSUPWALKR  | Aluminum, 5" Wheels, Red  | 300 lb (136.1 kg) | PVC             | 46"-53"      | 20"             | -     | N           | 1 ea |
| MDSUPWALKJR | Aluminum, Youth 5" Wheels | 300 lb (136.1 kg) | Cushion         | 35.5"-42.5"  | 18"             | -     | N           | 1 ea |

## Folding Trigger Walker\*

- Trigger mechanism can be pulled up to fold the walker
- Wheeled glidebrakes to help user brake easily
- Includes a bag for easy storage when travelling
- Two storage pouches included

| Item No.     | Description         | Wt. Capacity      | Handle Material | Height Range | Hand Grip Width | HCPCS | Retail Pkg. | Pkg. |
|--------------|---------------------|-------------------|-----------------|--------------|-----------------|-------|-------------|------|
| MDS86410TRIG | Aluminum, 5" Wheels | 300 lb (136.1 kg) | PVC             | 29.5"-37"    | 17"             | -     | N           | 1 ea |



MDS86410KDBW



MDSUPWALK



MDS86410TRIG

\*All walkers come unassembled



## Guardian Walkers

- All Guardian walkers come assembled
- Easy-fold sides for easy movement through narrow spaces
- Comfortable grip and mark-resistant tips

| Item No. | Description                             | Wt. Capacity      | Handle Material | Height Range  | Hand Grip Width | HCPCS | Retail Pkg. | Pkg. |
|----------|---|-------------------|-----------------|---------------|-----------------|-------|-------------|------|
| G30757W  | Aluminum, 5" Wheels, Adult              | 300 lb (136.1 kg) | PVC             | 32.5"-39.5"   | 16.5"           | E0143 | N           | 4/cs |
| G30758W  | Aluminum, 5" Wheels, Youth              | 300 lb (136.1 kg) | PVC             | 27.25"-34.25" | 16.5"           | E0143 | N           | 4/cs |
| G30754B  | Bariatric, Steel, No Wheels, Adult      | 650 lb (294.8 kg) | PVC             | 32"-39"       | 22"             | E0148 | N           | 2/cs |
| G30755P  | Aluminum, No Wheels, Adult              | 300 lb (136.1 kg) | PVC             | 32.5"-39.5"   | 16.5"           | E0135 | N           | 4/cs |
| G30756P  | Aluminum, No Wheels, Youth              | 300 lb (136.1 kg) | PVC             | 27"-34"       | 16.5"           | E0135 | N           | 4/cs |
| G07783M  | Aluminum, 5" Wheels, Pediatric, Tweener | 175 lb (79.4 kg)  | PVC             | 24"-36"       | 9.5"            | -     | N           | 1 ea |
| G07749   | Aluminum, No Wheels, Pediatric, Child   | 200 lb. (90.7 kg) | PVC             | 24.5"-28.5"   | 12.75"          | E0135 | N           | 1 ea |
| G07750   | Aluminum, No Wheels, Pediatric, Toddler | 175 lb (79.4 kg)  | PVC             | 20.5"-24.5"   | 10.75"          | E0135 | N           | 1 ea |
| G07767   | Aluminum, No Wheels, Heavy Duty, Adult  | 400 lb (181.4 kg) | PVC             | 32.5"-39.5"   | 16.5"           | E0148 | N           | 1 ea |



G30754B



G30755P



G07749

\*All MDS walkers come unassembled

## Walker Accessories

| Item No.                | Description                                     | Retail Pkg. | Pkg.    |
|-------------------------|---|-------------|---------|
| G07729M                 | 7/8" Footpiece, 3" Wheels and Glide Caps        | Y           | 1 pr    |
| G07731-8                | 1" Footpiece, 3" Wheels and Glide Caps          | Y           | 1 pr    |
| G07731G                 | Footpiece, 3" Wheels                            | Y           | 3 pr/cs |
| G07722-8                | Guardian, 5" Front Walker Wheels                | Y           | 1 pr    |
| MDS86615W               | 3" Casters for Walkers                          | Y           | 2 pr/cs |
| MDS86615W5              | Footpiece, 5" Wheels and Glide Caps             | Y           | 1 pr    |
| MDS86615W5W             | Footpiece, 5" Wheels and Glide Caps             | Y           | 6 pr/cs |
| G07724                  | Pediatric Extension Legs, 5" Wheels             | Y           | 2 pr/cs |
| G07725                  | Walker Footpiece Extension Kit, 5" Wheels       | Y           | 1 kit   |
| G07825-8                | Swivel Footpiece, 3" Wheels                     | Y           | 1 pr    |
| G07722-8B               | Fixed 5" Bariatric Dual Wheels                  | Y           | 1 pr    |
| G07722G                 | Fixed 5" Bariatric Dual Wheels                  | Y           | 3 pr/cs |
| G07723M                 | Footpiece, 5" Wheels                            | Y           | 1 pr    |
| MDS86615W               | Footpiece, 3" Wheels                            | Y           | 1 pr    |
| G07820-8M               | For All Guardian Walkers                        | Y           | 1 pr    |
| MDS86615W5              | Footpiece, 5" Wheels                            | Y           | 2 pr/cs |
| G07732G                 | Glidebrakes                                     | Y           | 3 pr/cs |
| MDS86615T               | Rubber Tips                                     | Y           | 1 pr    |
| G07900CS                | Glide Caps                                      | Y           | 24/cs   |
| G07904                  | Glide Caps                                      | Y           | 6/cs    |
| MDS86615C               | 1" Glide Caps for Walkers                       | Y           | 1 pr    |
| MDS86615CW              | 1" Glide Caps for Walkers                       | Y           | 6/cs    |
| MDS86615SKI             | Ski Gliders for Walkers                         | Y           | 6/cs    |
| G07712                  | 4-Piece Footpiece Extension                     | Y           | 2 pr/pk |
| MDS86615L               | 10" Leg Replacement                             | Y           | 1 pr    |
| G07850M                 | Flip Tray for Guardian Walkers                  | Y           | 3/cs    |
| G07850MB                | Bariatric Walker Flip Tray                      | Y           | 3/cs    |
| G07702                  | Walker Platform Attachment for Guardian Walkers | Y           | 1 ea    |
| G07703M                 | Walker Platform Attachment Youth                | Y           | 1 ea    |
| MDS86615P               | Platform Attachment for MDS Aluminum Walkers    | Y           | 1 ea    |
| G07741                  | Walker Front Carry Pouch                        | Y           | 3/cs    |
| G07715                  | Carrying Basket                                 | Y           | 2/cs    |
| G07872                  | Walker Seat                                     | Y           | 2/cs    |
| MDS86615K               | Walker Basket                                   | Y           | 2/cs    |
| MDS86615PB              | Platform Attachment for Basic Walkers           | Y           | 1 ea    |
| MDSWBALL                | Walker Ball Glides                              | Y           | 6 pr/cs |
| MDS86615BAG             | 2-Button Walker Bag                             | Y           | 4/cs    |
| MDSRUBSKI               | Walker Glides                                   | Y           | 6/cs    |
| MDS86615W5T             | 5" Long Extension Wheels                        | Y           | 1 pr    |
| <b>New</b> MDSUCPCANEHW | Combo Cup and Cane Holder                       | N           | 6/cs    |



# Walker Accessories Comparison Chart

● Fits    ■ Does not fit

|   | Single   | Paddle   | Non-Standard       | Basic  | Easy-Care                                | 2-Button  | Single   | Bariatric  |
|---|----------|--|--------------------|--|--|---|----------|--|
|   | MDS86617 | MDS86617<br>MDS86410P<br>MDS86410P5<br>MDS86410P5J | G07767<br>MDS86204 | MDS86410JW54B<br>MDS86410W54B<br>MDS864104B<br>MDS86410J4B | G30755P<br>G30757W<br>G30756P<br>G30758W | MDS864104<br>MDS86410W54<br>MDS86410J4<br>MDS86410JW54<br>MDS86410W4<br>MDS86410JW4 | MDS86615 | G30754B<br>MDS86410XWB<br>MDS86410XWWB<br>MDS86410XWWW<br>MDS86410XW |
| <b>G07731-8</b><br>3" Walker Wheels, 8-Hole, 1 pr         | ●        | ■ ● ● ● ●  | ■ ■                | ● ● ● ●  | ● ● ● ●                                  | ● ● ● ● ● ●   | ●        | ■ ● ● ● ● ■  |
| <b>G07731G</b><br>3" Walker Wheels, 3/cs                  | ●        | ■ ● ● ● ●  | ■ ■                | ● ● ● ●  | ● ● ● ●                                  | ● ● ● ● ● ●   | ●        | ■ ● ● ● ● ■  |
| <b>G07729M</b><br>3" Walker Wheels, 1 pr                  | ■        | ■ ■ ■ ■ ■  | ■ ■                | ■ ■ ■ ■ ■  | ■ ■ ■ ■ ■                                | ■ ■ ■ ■ ■ ■ ■   | ■        | ■ ■ ■ ■ ■ ■ ■  |
| <b>MDS86615W</b><br>3" Walker Wheels with Glidecaps, 2/st | ●        | ■ ● ● ● ●  | ■ ■                | ● ● ● ●  | ● ● ● ●                                  | ● ● ● ● ● ●   | ●        | ■ ● ● ● ● ■  |
| <b>G07722-8</b><br>5" Walker Wheels, 8-Holes, 1 pr        | ●        | ■ ● ● ● ●  | ■ ■                | ● ● ● ●  | ● ● ● ●                                  | ● ● ● ● ● ●   | ●        | ■ ● ● ● ● ■  |
| <b>G07722G</b><br>5" Wheels for Walker, 3/cs              | ●        | ■ ● ● ● ●  | ■ ■                | ● ● ● ●  | ● ● ● ●                                  | ● ● ● ● ● ●   | ●        | ■ ● ● ● ● ■  |
| <b>MDS86615W5</b><br>5" Walker Wheels, 2/st               | ●        | ■ ● ● ● ●  | ■ ■                | ● ● ● ●  | ● ● ● ●                                  | ● ● ● ● ● ●   | ●        | ■ ● ● ● ● ■  |
| <b>MDS86615W5T</b><br>5" Walker Wheels, 2/st              | ●        | ■ ● ● ● ●  | ■ ■                | ● ● ● ●  | ■ ■ ■ ■ ■                                | ■ ■ ■ ■ ■ ■ ■   | ■        | ■ ● ● ● ● ■  |
| <b>G07724</b><br>5" Walker Wheels (7/8" Frame), 1 kt      | ■        | ■ ■ ■ ■ ■  | ■ ■                | ■ ■ ■ ■ ■  | ■ ■ ■ ■ ■                                | ■ ■ ■ ■ ■ ■ ■   | ■        | ■ ■ ■ ■ ■ ■ ■  |
| <b>G07725</b><br>5" Walker Wheels, 1 kt                   | ●        | ■ ● ● ● ●  | ■ ■                | ● ● ● ●  | ■ ● ● ●                                  | ● ● ● ● ● ●   | ●        | ■ ● ● ● ● ■  |
| <b>G07722-8B</b><br>5" Duel Wheels, Bariatric, 1 pr       | ●        | ■ ● ● ● ●  | ■ ■                | ● ● ● ●  | ■ ■ ■ ■ ■                                | ■ ■ ■ ■ ■ ■ ■   | ■        | ● ● ● ● ● ●  |
| <b>G07825-8</b><br>3" Walker Wheels, Swivel, 1 pr         | ●        | ■ ● ● ● ●  | ■ ■                | ● ● ● ●  | ● ● ● ●                                  | ● ● ● ● ● ●   | ●        | ■ ● ● ● ● ■  |
| <b>G07811</b><br>5" Walker Wheels, 1 pr                   | ●        | ■ ● ● ● ●  | ■ ■                | ● ● ● ●  | ■ ■ ■ ■ ■                                | ■ ■ ■ ■ ■ ■ ■   | ■        | ■ ● ■ ● ■ ● ■  |
| <b>MDS86615W5</b><br>5" Walker Swivel Wheels, 2/st        | ●        | ■ ● ● ● ●  | ■ ■                | ● ● ● ●  | ■ ■ ■ ■ ■                                | ■ ■ ■ ■ ■ ■ ■   | ■        | ■ ● ● ● ● ■  |
| <b>G07732G</b><br>Brake Footpiece Set, 3/cs               | ●        | ■ ● ● ● ●  | ■ ■                | ● ● ● ●  | ● ● ● ●                                  | ● ● ● ● ● ●   | ●        | ■ ● ● ● ● ■  |
| <b>G07900CS</b><br>Glidecaps, 24/cs                       | ●        | ● ● ● ● ●  | ● ●                | ● ● ● ●  | ● ● ● ●                                  | ● ● ● ● ● ●   | ●        | ● ● ● ■ ● ● ●  |
| <b>G07904</b><br>Glidecaps, 6/cs                          | ●        | ● ● ● ● ●  | ● ●                | ● ● ● ●  | ● ● ● ●                                  | ● ● ● ● ● ●   | ●        | ● ● ● ■ ● ● ●  |
| <b>G07712</b><br>Extension Footpieces, 1/pk               | ●        | ■ ● ● ● ●  | ■ ■                | ● ● ● ●  | ● ● ● ●                                  | ● ● ● ● ● ●   | ●        | ■ ● ● ■ ● ● ■  |
| <b>MDS86615L</b><br>Walker Leg Replacement, 2/st          | ●        | ■ ● ● ● ●  | ● ●                | ● ● ● ●  | ● ● ● ●                                  | ● ● ● ● ● ●   | ●        | ■ ● ● ■ ● ● ■  |
| <b>G07741</b><br>Walker Pouch, 3/cs                       | ●        | ● ● ● ● ●  | ● ●                | ● ● ● ●  | ● ● ● ●                                  | ● ● ● ● ● ●   | ●        | ■ ● ● ● ● ●  |
| <b>G07715</b><br>Carrying Basket, 2/cs                    | ■        | ■ ■ ■ ■ ■  | ■ ■                | ■ ■ ■ ■ ■  | ● ● ● ●                                  | ■ ■ ■ ■ ■ ■ ■   | ■        | ■ ■ ■ ■ ■ ■ ■  |



# Wheelchairs

When choosing a wheelchair, there are many factors to consider including Arm Style, Leg Style, Weight Capacity and Seat Dimension. We've issued the charts in a way that makes it easier to select the right one.

**Arm Style:** Many users require the armrests in order to safely perform transfers in and out of a wheelchair. Individuals who can stand also use the armrests for support when getting in and out of the wheelchair. Properly fitted armrests when used help in maintaining a correct and functional seated posture. Arm Style choices include:

- Desk-length swing back and removable
- Desk-length swing-back and height adjustable
- Desk-length height-adjustable swing-back and removable
- Full-length removable
- Full-length removable and height adjustable
- Full-length swing-back and removable
- Full-length height-adjustable swing-back and removable

**Leg Style:** The main purpose of these leg rests is to redistribute the pressure on the user and increase comfort. There are two options for leg style:

- Swing-Away
- Elevating

**Weight Capacity:** Wheelchairs come with varying weight limits. Basic wheelchairs have weight limits of 250–350 lbs, lightweight or pediatric wheelchairs have weight limits of 200–250 lbs, and extra wide and bariatric wheelchairs have weight limits of 500–850 lbs.

**Seat Dimensions:** Seat widths typically range from 16 to 20 inches, with an average width of about 18 inches. Take seat depth into consideration as it will make a difference in overall comfort and usage. Seat depth is the distance from the back of the seat to the front edge and should be measured when the user is sitting straight up.




# Wheelchair Sizing Guide

Determine the best fit for your wheelchair in six easy steps.

The most important step to consider when choosing a wheelchair is a proper fit. It can make the difference between recovering well and causing an injury.

## Improper fit can cause:

- Posture issues
- Breathing problems
- Sores
- Discomfort
- Pelvic issues
- Hip problems



## 1 Determine Arm Type and Height

- **Full-length arm:** Gives support to pivot transfers when the patient pushes off to stand.
- **Desk-length arm:** Makes it easy to move around tables and desks.
- **Height:** Measure from the elbow to the seat of the chair while holding arms up with elbows bent at a 90-degree angle. Use adjustable arms when possible.

## 2 Determine Footrest Style

- **Elevated leg rest:** Raises legs
- **Articulation leg rest:** Extends longer for taller patients.
- **Length:** Measure from back of the knee to the heel of the foot.

## 3 Determine Wheelchair Weight and Weight Limit

- **Weight capacity:** Choose wheelchair weight based on the patient's weight.
- **Upper body strength:** Choose lighter wheelchairs for weaker patients.

## 4 Determine Seat Size

- **Seat width:** Measure from hip to hip in a straight line and then add two inches.
- **Seat depth:** Measure from the back of the hip to the back of the knee while seated and then subtract one inch.

## 5 Determine Back Height

- **Standard back height:** Measure from collarbone to the seat while seated.
- **Taller back height:** Require additional upper trunk support to be installed.

## 6 Determine Floor to Seat Height

- **Movement:** Check if the patient uses feet to propel or move.
- **Height:** Measure the distance between the back of the knee and the heel.

### K1 Guardian Wheelchairs

- Durable hammertone frame; available in nylon or vinyl upholstery
- Hemi-height adjustable
- Accommodates 2" seat extender

#### Nylon

| Arm Style              | Leg Style  | Item No.  | Frame Color | Wt. Capacity      | Seat Dimensions | HCPCS                     | Retail Pkg. | Pkg. |
|------------------------|------------|-----------|-------------|-------------------|-----------------|---------------------------|-------------|------|
| Desk-Length Swing Back | Swing-Away | K1166N22S | Black       | 300 lb (136.1 kg) | 16" x 16"       | K0001 OR<br>K0002 + E0951 | N           | 1 ea |
|                        |            | K1186N22S | Black       | 300 lb (136.1 kg) | 18" x 16"       | K0001 OR<br>K0002 + E0951 | N           | 1 ea |
|                        |            | K1206N22S | Black       | 300 lb (136.1 kg) | 20" x 16"       | K0001 OR<br>K0002 + E0951 | N           | 1 ea |
|                        | Elevating  | K1166N22E | Black       | 300 lb (136.1 kg) | 16" x 16"       | K0001 OR<br>K0002 + E0990 | N           | 1 ea |
|                        |            | K1186N22E | Black       | 300 lb (136.1 kg) | 18" x 16"       | K0001 OR<br>K0002 + E0990 | N           | 1 ea |
|                        |            | K1206N22E | Black       | 300 lb (136.1 kg) | 20" x 16"       | K0001 OR<br>K0002 + E0990 | N           | 1 ea |
| Full-Length Removable  | Swing-Away | K1186N13S | Black       | 300 lb (136.1 kg) | 18" x 16"       | K0001 OR<br>K0002 + E0951 | N           | 1 ea |
|                        | Elevating  | K1186N13E | Black       | 300 lb (136.1 kg) | 18" x 16"       | K0001 OR<br>K0002 + E0990 | N           | 1 ea |

#### Vinyl

| Arm Style              | Leg Style  | Item No.  | Frame Color | Wt. Capacity      | Seat Dimensions | HCPCS                    | Retail Pkg. | Pkg. |
|------------------------|------------|-----------|-------------|-------------------|-----------------|--------------------------|-------------|------|
| Desk-Length Swing Back | Swing-Away | K1166V22S | Black       | 300 lb (136.1 kg) | 16" x 16"       | K0006                    | N           | 1 ea |
|                        |            | K1186V22S | Black       | 300 lb (136.1 kg) | 18" x 16"       | K0006                    | N           | 1 ea |
|                        |            | K1206V22S | Black       | 300 lb (136.1 kg) | 20" x 16"       | K0006 + E2201            | N           | 1 ea |
|                        | Elevating  | K1166V22E | Black       | 300 lb (136.1 kg) | 16" x 16"       | K0006 + K0195            | N           | 1 ea |
|                        |            | K1186V22E | Black       | 300 lb (136.1 kg) | 18" x 16"       | K0006 + K0195            | N           | 1 ea |
|                        |            | K1206V22E | Black       | 300 lb (136.1 kg) | 20" x 16"       | K0006 + E2201 +<br>K0195 | N           | 1 ea |
| Full-Length Removable  | Swing-Away | K1186V11S | Black       | 300 lb (136.1 kg) | 18" x 16"       | -                        | N           | 1 ea |
|                        | Elevating  | K1186V11E | Black       | 300 lb (136.1 kg) | 18" x 16"       | -                        | N           | 1 ea |



## K2 Guardian Wheelchairs

- Durable Hammertone frame
- Ships in hemi height; includes seat belt, anti-tippers and brake extenders

### Nylon

| Arm Style                                   | Leg Style | Item No.   | Frame Color | Wt. Capacity      | Seat Dimensions | HCPCS                 | Retail Pkg. | Pkg. |
|---|-----------|------------|-------------|-------------------|-----------------|-----------------------|-------------|------|
| Full-Length Removable and Height Adjustable | Elevating | K2166N42E1 | Black       | 300 lb (136.1 kg) | 16" x 16"       | K0002 + E0973 + K0195 | N           | 1 ea |
|   |           | K2186N42E1 | Black       | 300 lb (136.1 kg) | 18" x 16"       | K0002 + E0973 + K0195 | N           | 1 ea |
|   |           | K2206N42E1 | Black       | 300 lb (136.1 kg) | 20" x 16"       | K0002 + E0973 + K0195 | N           | 1 ea |

### New Guardian Extra-Wide Chairs

- Durable Hammertone frame and wipeable vinyl upholstery
- Heavy-duty frame with 500lb weight capacity
- Hemi-height adjustable

### Vinyl

| Arm Style                  | Leg Style  | Item No.  | Frame Color | Wt. Capacity      | Seat Dimensions | HCPCS | Retail Pkg. | Pkg. |
|----------------------------|------------|-----------|-------------|-------------------|-----------------|-------|-------------|------|
| Removable Desk-Length Arms | Swing-Away | KX208V21S | Black       | 500 lb (226.8 kg) | 20" x 18"       | -     | N           | 1 ea |
|                            |            | KX228V21S | Black       | 500 lb (226.8 kg) | 22" x 18"       | -     | N           | 1 ea |
|                            |            | KX248V21S | Black       | 500 lb (226.8 kg) | 24" x 18"       | -     | N           | 1 ea |
|                            |            | KX268V21S | Black       | 500 lb (226.8 kg) | 26" x 18"       | -     | N           | 1 ea |
|                            | Elevating  | KX208V21E | Black       | 500 lb (226.8 kg) | 20" x 18"       | -     | N           | 1 ea |
|                            |            | KX228V21E | Black       | 500 lb (226.8 kg) | 22" x 18"       | -     | N           | 1 ea |
|                            |            | KX248V21E | Black       | 500 lb (226.8 kg) | 24" x 18"       | -     | N           | 1 ea |
|                            |            | KX268V21E | Black       | 500 lb (226.8 kg) | 26" x 18"       | -     | N           | 1 ea |

Accommodates 2" seat extenders



K2166N42E1

Includes anti-tippers



MDS806550PLUS

# The K3 wheelchair comes standard with the accessories your customers are looking for.

## New K3 Guardian Wheelchairs

- Durable Hammertone frame and easy-to-clean nylon upholstery
- Swing-back arms are also removable for easy transfers and storage
- Hemi-height adjustable

Height-adjustable swing-back arms for comfort and easy side transfers

Seat belt prevents sliding or falling out of the wheelchair

Choose from elevating or swing away foot rests

2" (5.1cm) seat depth extender for customer comfort


Break extenders allow the user to independently lock the breaks



K3186N44E1



Anti-tippers for added safety



Heel loops for support and to prevent sliding (swing-away footrest models only)



## New K3 Guardian Wheelchairs

- Durable Hammertone frame and easy-to-clean nylon upholstery
- Swing-back arms are also removable for easy transfers and storage
- Hemi-height adjustable

| Arm Style   | Leg Style  | Item No.    | Frame Color | Wt. Capacity      | Seat Dimensions | HCPCS                         | Retail Pkg. | Pkg. |
|---|------------|-------------|-------------|-------------------|-----------------|-------------------------------|-------------|------|
| Desk-Length<br>Swing-Back<br>and Removable                            | Swing-Away | K3166N24S   | Black       | 300 lb (136.1 kg) | 16" x 16"       | K0006                         | N           | 1 ea |
|   |            | K3186N24S   | Black       | 300 lb (136.1 kg) | 18" x 16"       | K0006                         | N           | 1 ea |
|   |            | K3206N24S   | Black       | 300 lb (136.1 kg) | 20" x 16"       | K0006 + E2201                 | N           | 1 ea |
|   | Elevating  | K3166N24E   | Black       | 300 lb (136.1 kg) | 16" x 16"       | K0006 + K0195                 | N           | 1 ea |
|   |            | K3186N24E   | Black       | 300 lb (136.1 kg) | 18" x 16"       | K0006 + K0195                 | N           | 1 ea |
|   |            | K3206N24E   | Black       | 300 lb (136.1 kg) | 20" x 16"       | K0006 + K0195 + E2201         | N           | 1 ea |
| Desk-Length,<br>Height-<br>Adjustable,<br>Swing-Back<br>and Removable | Swing-Away | K3166N34S   | Black       | 300 lb (136.1 kg) | 16" x 16"       | K0006 + E0973                 | N           | 1 ea |
|   |            | K3186N34S   | Black       | 300 lb (136.1 kg) | 18" x 16"       | K0006 + E0973                 | N           | 1 ea |
|   |            | K3206N34S   | Black       | 300 lb (136.1 kg) | 20" x 16"       | K0006 + E2201 + E0973         | N           | 1 ea |
|   | Elevating  | K3166N34E   | Black       | 300 lb (136.1 kg) | 16" x 16"       | K0006 + K0195 + E0973         | N           | 1 ea |
|   |            | K3186N34E   | Black       | 300 lb (136.1 kg) | 18" x 16"       | K0006 + K0195 + E0973         | N           | 1 ea |
|   |            | K3206N34E   | Black       | 300 lb (136.1 kg) | 20" x 16"       | K0006 + K0195 + E2201 + E0973 | N           | 1 ea |
| Full-Length,<br>Swing Back<br>and Removable                           | Swing-Away | K3166N14S   | Black       | 300 lb (136.1 kg) | 16" x 16"       | K0006                         | N           | 1 ea |
|   |            | K3186N14S   | Black       | 300 lb (136.1 kg) | 18" x 16"       | K0006                         | N           | 1 ea |
|   |            | K3206N14S   | Black       | 300 lb (136.1 kg) | 20" x 16"       | K0006 + E2201                 | N           | 1 ea |
|   | Elevating  | K3166N14E   | Black       | 300 lb (136.1 kg) | 16" x 16"       | K0006 + K0195                 | N           | 1 ea |
|   |            | K3186N14E   | Black       | 300 lb (136.1 kg) | 18" x 16"       | K0006 + K0195                 | N           | 1 ea |
|   |            | K3206N14E   | Black       | 300 lb (136.1 kg) | 20" x 16"       | K0006 + K0195 + E2201         | N           | 1 ea |
| Full-Length,<br>Height-<br>Adjustable,<br>Swing Back<br>and Removable | Swing-Away | K3166N44S   | Black       | 300 lb (136.1 kg) | 16" x 16"       | K0006 + E0973                 | N           | 1 ea |
|   |            | K3186N44S   | Black       | 300 lb (136.1 kg) | 18" x 16"       | K0006 + E0973                 | N           | 1 ea |
|   |            | K3206N44S   | Black       | 300 lb (136.1 kg) | 20" x 16"       | K0006 + E2201 + E0973         | N           | 1 ea |
|   | Elevating  | K3166N44E   | Black       | 300 lb (136.1 kg) | 16" x 16"       | K0006 + K0195 + E0973         | N           | 1 ea |
|   |            | K3186N44E   | Black       | 300 lb (136.1 kg) | 18" x 16"       | K0006 + K0195 + E0973         | N           | 1 ea |
|   |            | K3206N44E   | Black       | 300 lb (136.1 kg) | 20" x 16"       | K0006 + K0195 + E2201 + E0973 | N           | 1 ea |
|   |            | K3166N44E1* | Black       | 300 lb (136.1 kg) | 16" x 16"       | K0006 + K0195                 | N           | 1 ea |
|   |            | K3186N44E1* | Black       | 300 lb (136.1 kg) | 18" x 16"       | K0006 + K0195                 | N           | 1 ea |
|   |            | K3206N44E1* | Black       | 300 lb (136.1 kg) | 20" x 16"       | K0006 + K0195 + E2201         | N           | 1 ea |

\*Includes anti-tippers, seat extension, brake extension and seat belt.



## K4 Lightweight Wheelchairs

- Durable, lightweight TIG-welded frame folds down for easy transport

### Nylon

| Arm Style               | Leg Style  | Item No.     | Frame Color       | Wt. Capacity      | Seat Dimensions       | HCPCS                  | Retail Pkg. | Pkg. |
|-------------------------|------------|--------------|-------------------|-------------------|-----------------------|------------------------|-------------|------|
| Desk-Length, Swing Back | Swing-Away | MDS806500N   | Black             | 300 lb (136.1 kg) | 16" x 16"             | K0004                  | N           | 1 ea |
|                         |            | MDS806500    | Black             | 300 lb (136.1 kg) | 18" x 16"             | K0004                  | N           | 1 ea |
|                         |            | MDS806560    | Black             | 350 lb (158.8 kg) | 20" x 16"             | K0004+E2201            | N           | 1 ea |
|                         | Elevating  | MDS806570    | Black             | 350 lb (158.8 kg) | 22" x 16"             | K0004 + E2201          | N           | 1 ea |
|                         |            | MDS806550N   | Black             | 300 lb (136.1 kg) | 16" x 16"             | K0004 + K0195 or E0990 | N           | 1 ea |
|                         |            | MDS806550    | Black             | 300 lb (136.1 kg) | 18" x 16"             | K0004                  | N           | 1 ea |
|                         |            | MDS806565    | Black             | 350 lb (158.8 kg) | 20" x 16"             | K0004+E2201+K0195      | N           | 1 ea |
| Full-Length Removable   | MDS806575  | Black        | 350 lb (158.8 kg) | 22" x 16"         | K0004 + E2201 + K0195 | N                      | 1 ea        |      |
|                         | Swing-Away | MDS806500FLA | Black             | 300 lb (136.1 kg) | 18" x 16"             | -                      | N           | 1 ea |
|                         | Elevating  | MDS806550FLA | Black             | 300 lb (136.1 kg) | 18" x 16"             | -                      | N           | 1 ea |

## K4 Basic Wheelchairs

- Durable Grey frame
- Hemi height adjustable

| Arm Style | Leg Style | Item No.      | Frame Color | Wt. Capacity       | Seat Dimensions | HCPCS                                   | Retail Pkg. | Pkg. |
|-----------|-----------|---------------|-------------|--------------------|-----------------|---|-------------|------|
|           |           | MDS806500E    | Silver      | 300 lb (136.1 kg)  | 18" x 16"       | K0004                                   | N           | 1 ea |
|           |           | MDS806500NE18 | Silver      | 300 lb (136.1 kg)  | 16" x 18"       | K0004                                   | N           | 1 ea |
|           |           | MDS806560E    | Silver      | 300 lb (136.1 kg)  | 20" x 18"       | K0006 + E2201                           | N           | 1 ea |
|           | Elevating | MDS806550NE   | Silver      | 300 lb (136.1 kg)  | 16" x 16"       | K0004 + K0195                           | N           | 1 ea |
|           |           | MDS806550E    | Silver      | 300 lb (136.1 kg)  | 18" x 16"       | K0004+K0195                             | N           | 1 ea |
|           |           | MDS806500E18  | Silver      | 300 lb (136.1 kg)  | 18" x 18"       | -                                       |             | 1 ea |
|           |           | MDS806565E    | Silver      | 300 lb. (136.1 kg) | 20" x 18"       | K0003 or K0006 + E2201 + K0195 or E0990 | N           | 1 ea |
|           |           |               |             |                    |                 |   |             |      |



## K4 Plus Wheelchairs

- Durable Hammertone frame
- Includes anti-tippers

| Arm Style                                     | Leg Style  | Item No.      | Frame Color | Wt. Capacity      | Seat Dimensions | HCPCS                 | Retail Pkg. | Pkg. |
|---|------------|---------------|-------------|-------------------|-----------------|-----------------------|-------------|------|
| Desk-Length, Swing-Back and Height Adjustable | Swing-Away | MDS806500PLUS | Black       | 300 lb (136.1 kg) | 18" x 16"       | -                     | N           | 1 ea |
|   |            | MDS806560PLUS | Black       | 300 lb (136.1 kg) | 20" x 18"       | K0004 + E2201 + E0971 | N           | 1 ea |
|   | Elevating  | MDS806550PLUS | Black       | 300 lb (136.1 kg) | 18" x 16"       | -                     | N           | 1 ea |

## New Guardian K4 Wheelchairs

- Durable Hammertone frame and easy to clean nylon upholstery
- Hemi-height adjustable
- Quick-release rear wheels

### Nylon

| Arm Style  | Leg Style  | Item No.  | Frame Color | Wt. Capacity      | Seat Dimensions | HCPCS | Retail Pkg. | Pkg. |
|--|------------|-----------|-------------|-------------------|-----------------|-------|-------------|------|
| Desk-Length, Swing Back and Removable                    | Swing-Away | K4166N24S | Black       | 300 lb (136.1 kg) | 16" x 16"       | -     | N           | 1 ea |
|  |            | K4186N24S | Black       | 300 lb (136.1 kg) | 18" x 16"       | -     | N           | 1 ea |
|  |            | K4206N24S | Black       | 300 lb (136.1 kg) | 20" x 16"       | -     | N           | 1 ea |
|  | Elevating  | K4166N24E | Black       | 300 lb (136.1 kg) | 16" x 16"       | -     | N           | 1 ea |
|  |            | K4206N24E | Black       | 300 lb (136.1 kg) | 18" x 16"       | -     | N           | 1 ea |
|  |            | K4186N24E | Black       | 300 lb (136.1 kg) | 20" x 16"       | -     | N           | 1 ea |
| Desk-Length, Height-Adjustable, Swing Back and Removable | Swing-Away | K4186N34S | Black       | 300 lb (136.1 kg) | 18" x 16"       | -     | N           | 1 ea |
|  | Elevating  | K4186N34E | Black       | 300 lb (136.1 kg) | 18" x 16"       | -     | N           | 1 ea |



K4186N24S

## Guardian Recliners

- Infinite position reclining up to 180°
- Durable hammertone frame

### Nylon

| Arm Style              | Leg Style | Item No.  | Frame Color | Wt. Capacity      | Seat Dimensions | HCPCS                 | Retail Pkg. | Pkg. |
|------------------------|-----------|-----------|-------------|-------------------|-----------------|-----------------------|-------------|------|
| Desk-Length, Removable | Elevating | KR168N21E | Black       | 300 lb (136.1 kg) | 16" x 18"       | K0006 + E0990 + E1226 | N           | 1 ea |
|                        |           | KR188N21E | Black       | 300 lb (136.1 kg) | 18" x 18"       | K0006 + E0990 + E1226 | N           | 1 ea |
|                        |           | KR208N21E | Black       | 300 lb (136.1 kg) | 20" x 18"       | K0006 + E0990 + E1226 | N           | 1 ea |
| Full-Length Removable  | Elevating | KR188N11E | Black       | 300 lb (136.1 kg) | 18" x 18"       | K0006 + E0990 + E1226 | N           | 1 ea |

## Reclining Wheelchairs

- Infinite position reclining from 90°–140°
- Durable carbon steel frame
- Hemi-height adjustable

### Vinyl

| Arm Style              | Leg Style | Item No.  | Frame Color | Wt. Capacity      | Seat Dimensions | HCPCS  | Retail Pkg. | Pkg. |
|------------------------|-----------|-----------|-------------|-------------------|-----------------|--|-------------|------|
| Desk-Length, Removable | Elevating | MDS808350 | Chrome      | 300 lb (136.1 kg) | 16" x 17"       | K0003 + E0966 + E1226 + K0195 + E0971                  | N           | 1 ea |
|                        |           | MDS808450 | Chrome      | 300 lb (136.1 kg) | 18" x 17"       | K0003 + E0966 + E1226 + K0195 + E0971                  | N           | 1 ea |
|                        |           | MDS808550 | Chrome      | 300 lb (136.1 kg) | 20" x 17"       | K0003 + E0966 + E1226 + K0195 + E0971 + E2201          | N           | 1 ea |
|                        |           | MDS808650 | Chrome      | 300 lb (136.1 kg) | 22" x 17"       | K0007 + E1226 + E2201 + E0971 + E0966 + K0195 or E0990 | N           | 1 ea |

## Hybrid Wheelchairs

- Combination wheelchair and transport chair
- Quick release large wheel

### Nylon

| Arm Style              | Leg Style  | Item No.     | Frame Color | Wt. Capacity      | Seat Dimensions | HCPCS                           | Retail Pkg. | Pkg. |
|------------------------|------------|--------------|-------------|-------------------|-----------------|---------------------------------|-------------|------|
| Desk-Length, Removable | Swing-Away | MDS806250NH2 | Blue        | 300 lb (136.1 kg) | 16" x 16"       | K0004 or E1038                  | N           | 1 ea |
|                        |            | MDS806250H2  | Blue        | 300 lb (136.1 kg) | 18" x 16"       | K0001 or K0002                  | N           | 1 ea |
|                        | Elevating  | MDS806300H2  | Blue        | 300 lb (136.1 kg) | 18" x 16"       | K0004 or E1038 + K0195 or E0990 | N           | 1 ea |



### Excel Wheelchairs

- Durable carbon steel frame
- Hemi-height adjustable

#### Vinyl

| Arm Style              | Leg Style  | Item No.      | Frame Color | Wt. Capacity      | Seat Size | HCPCS                           | Retail Pkg. | Pkg. |
|------------------------|------------|---------------|-------------|-------------------|-----------|---------------------------------|-------------|------|
| Desk-Length, Removable | Swing-Away | MDS806250N    | Chrome      | 300 lb (136.1 kg) | 16" x 16" | K0001 or K0002                  | N           | 1 ea |
|                        |            | MDS806250D    | Chrome      | 300 lb (136.1 kg) | 18" x 16" | K0001 or K0002 + E0990 or K0195 | N           | 1 ea |
|                        |            | MDS806400     | Chrome      | 300 lb (136.1 kg) | 20" x 16" | K0001 + E2201                   | N           | 1 ea |
|                        | Elevating  | MDS806300N    | Chrome      | 300 lb (136.1 kg) | 16" x 16" | K0001 or K0002 + K0195          | N           | 1 ea |
|                        |            | MDS806300D    | Chrome      | 300 lb (136.1 kg) | 18" x 16" | K0001 or K0002                  | N           | 1 ea |
| Full-Length, Permanent | Swing-Away | MDS806450     | Chrome      | 300 lb (136.1 kg) | 20" x 16" | K0001 + E2201 + K0195           | N           | 1 ea |
|                        | Elevating  | MDS806150D    | Chrome      | 300 lb (136.1 kg) | 18" x 16" | K0001 or K0002                  | N           | 1 ea |
|                        |            | MDS806200N    | Chrome      | 300 lb (136.1 kg) | 16" x 16" | K0002 + K0195 or E0990          | N           | 1 ea |
| Full-Length, Removable | Elevating  | MDS806100D    | Chrome      | 300 lb (136.1 kg) | 18" x 16" | K0001 or K0002                  | N           | 1 ea |
|                        |            | MDS806400FLA  | Chrome      | 300 lb (136.1 kg) | 20" x 16" | -                               | N           | 1 ea |
|                        | Elevating  | MDS806300DFLA | Chrome      | 300 lb (136.1 kg) | 18" x 16" | K0001 or K0002                  | N           | 1 ea |

### Excel Extra-Wide Wheelchairs

- Durable carbon steel frame
- Hemi-height adjustable

#### Vinyl

| Arm Style              | Leg Style  | Item No.     | Frame Color | Wt. Capacity      | Seat Dimensions | HCPCS                  | Retail Pkg. | Pkg. |
|------------------------|------------|--------------|-------------|-------------------|-----------------|------------------------|-------------|------|
| Desk-Length, Removable | Swing-Away | MDS806700    | Chrome      | 500 lb (226.8 kg) | 20" x 18"       | K0007                  | N           | 1 ea |
|                        |            | MDS806800    | Chrome      | 500 lb (226.8 kg) | 22" x 18"       | K0007                  | N           | 1 ea |
|                        |            | MDS806900    | Chrome      | 500 lb (226.8 kg) | 24" x 18"       | K0007+E2202            | N           | 1 ea |
|                        | Elevating  | MDS806750    | Chrome      | 500 lb (226.8 kg) | 20" x 18"       | K0007 + K0195 or E0990 | N           | 1 ea |
|                        |            | MDS806850    | Chrome      | 500 lb (226.8 kg) | 22" x 18"       | K0007 + K0195 or E0990 | N           | 1 ea |
| Full-Length, Removable | Swing-Away | MDS806950    | Chrome      | 500 lb (226.8 kg) | 24" x 18"       | K0007 + E2202 + K0195  | N           | 1 ea |
|                        |            | MDS806700FLA | Chrome      | 500 lb (226.8 kg) | 20" x 18"       | K0007                  | N           | 1 ea |
|                        |            | MDS806800FLA | Chrome      | 500 lb (226.8 kg) | 22" x 18"       | K0007                  | N           | 1 ea |
|                        | Elevating  | MDS806900FLA | Chrome      | 500 lb (226.8 kg) | 24" x 18"       | K0007 + E2202          | N           | 1 ea |
|                        |            | MDS806750FLA | Chrome      | 500 lb (226.8 kg) | 20" x 18"       | K0007 + K0195 or E0990 | N           | 1 ea |
|                        |            | MDS806850FLA | Chrome      | 500 lb (226.8 kg) | 22" x 18"       | K0007 + K0195 or E0990 | N           | 1 ea |
|                        | Elevating  | MDS806950FLA | Chrome      | 500 lb (226.8 kg) | 24" x 18"       | K0007 + E2202 + K0195  | N           | 1 ea |



MDS806300D



MDS806900

### Shuttle Bariatric Wheelchairs

- Durable carbon steel frame
- Hemi-height adjustable

| Arm Style              | Leg Style | Item No.  | Frame Color | Wt. Capacity      | Seat Dimensions | HCPCS                                  | Retail Pkg. | Pkg. |
|------------------------|-----------|-----------|-------------|-------------------|-----------------|--|-------------|------|
| Desk-length, removable | Elevating | MDS809650 | Black       | 850 lb (385.5 kg) | 24" x 20"       | K0007 + K0195 or E0990 + E2202 + E2203 | N           | 1 ea |
|                        |           | MDS809750 | Black       | 850 lb (385.5 kg) | 26" x 20"       | K0007 + E2202 + E2203 + K0195 or E0990 | N           | 1 ea |
|                        |           | MDS809850 | Black       | 850 lb (385.5 kg) | 28" x 20"       | K0007 + K0108 + K0195 or E0990 + E2203 | N           | 1 ea |
|                        |           | MDS809950 | Black       | 850 lb (385.5 kg) | 30" x 20"       | K0007 + K0108 + K0195 or E0990 + E2203 | N           | 1 ea |

### Pediatric Wheelchairs

- Adjustable telescoping handles
- Hemi-height adjustable

#### Vinyl

| Arm Style               | Leg Style | Item No.    | Frame Color | Wt. Capacity      | Seat Dimensions | HCPCS                          | Retail Pkg. | Pkg. |
|-------------------------|-----------|-------------|-------------|-------------------|-----------------|--------------------------------|-------------|------|
| Desk-Length, Swing-Back | Elevating | MDS806140PD | Yellow      | 250 lb (113.4 kg) | 14" x 12"       | E1236 + E0971 + E0990 or K0195 | N           | 1 ea |

#### Nylon

| Arm Style               | Leg Style  | Item No.  | Frame Color | Wt. Capacity      | Seat Dimensions | HCPCS | Retail Pkg. | Pkg. |
|-------------------------|------------|-----------|-------------|-------------------|-----------------|-------|-------------|------|
| Desk-Length, Swing-Back | Swing-Away | KPD2N22S  | Yellow      | 250 lb (113.4 kg) | 12" x 12"       | E1236 | N           | 1 ea |
|                         |            | KPD4N22S  | Yellow      | 250 lb (113.4 kg) | 14" x 12"       | E1236 | N           | 1 ea |
|                         |            | KPD2N22S1 | Blue        | 250 lb (113.4 kg) | 12" x 12"       | E1236 | N           | 1 ea |
|                         |            | KPD4N22S1 | Blue        | 250 lb (113.4 kg) | 14" x 12"       | E1236 | N           | 1 ea |



MDS809750



MDS806140PD



KPD2N22S



KPD4N22S1

## Wheelchair Accessories

| Item No.      | Description   | Pkg.  |
|---------------|---|-------|
| MDS85190      | Oxygen Tank Holder and IV Pole for Wheelchair   | 1 ea  |
| WCA6100       | Oxygen Tank Holder for Adult Wheelchair, Holds D and E Tanks                                      | 1 ea  |
| MDS85189      | Rear Anti-Tip Device for Excel 2000 Wheelchair  | 1 pr  |
| STMSM22       | Anti-Rollback Device  | 1 ea  |
| MDS85189H     | Rear Anti-Tip Device for K1 Basic Wheelchair  | 1 pr  |
| MDS85181U     | Oxygen Tank Holder for K2 Basic, Excel 2000 and Extra-Wide Wheelchairs, Holds Standard Size Tanks | 1 ea  |
| MDS85183      | Telescoping IV Pole for Excel Extra-Wide Wheelchair   | 1 ea  |
| MDS85189WK4   | Anti-Tip Rear Wheel for K4 and Recliner   | 1 pr  |
| WCABELT       | Wheelchair Safety Belt  | 1 ea  |
| MDS85157PR    | Black Handle Extension Wheel Lock 5/8" x 8"   | 1 pr  |
| MDS85197A     | Adjustable Anti-Theft Device for Wheelchair   | 1 ea  |
| MDSK41RMDRIVE | Single-Arm Drive attachment for K4 Basic  | 1 ea  |
| WCA806991     | Tinnerman Nut for Wheelchair  | 1 pr  |
| WCA6200       | Adjustable Cup Holder   | 1 ea  |
| MDS85196      | Anti-Fold Device  | 1 ea  |
| MDS85195      | Super Hemi Height Kit for Wheelchairs with Standard Wheels  | 1 set |
| MDS851965SH   | Anti-Fold Device for Shuttle Extra-Wide Wheelchair  | 1 ea  |
| MDS85197SH    | Overhead Antitheft Device for Shuttle Extra-Wide Wheelchair                                       | 1 ea  |
| MDS85180      | 5-in-1 IV/O2/Anti-Theft Device  | 1 ea  |
| STMSM22W      | Anti-Rollback Device  | 1 ea  |
| MDS85189SH    | Rear Anti-Tipper for Shuttle Extra-Wide Wheelchair  | 1 pr  |
| MDS85189K     | Rear Anti-Tipper for Kidz Pediatric Wheelchair  | 1 pr  |
| WCAA0060      | K1/K2 Nylon 2" Seat Depth Extender  | 1 ea  |
| WCAA0061      | K1/K2 Nylon 2" Seat Depth Extender  | 1 ea  |
| WCAA0070      | K1 & K2 Anti-Tipper   | 1 pr  |
| WCAA0069      | K3 Anti-Tipper  | 1 pr  |
| WCAA0066      | K3 16" Seat Extender  | 1 ea  |
| WCAA0067      | K3 18" Seat Extender  | 1 ea  |
| WCAA0068      | K3 20" Seat Extender  | 1 ea  |



MDS85183

# Transport Chairs

Transport chairs make it easier for caregivers to help individuals get around comfortably and easily. Compared to standard wheelchairs, a transport chair requires a caregiver to push the chair.

A transport chair is much lighter than a wheelchair and come in a variety of colors and materials such as steel and aluminum. They are more compact, commonly come in 17" or 19" seat width to fit everyone's needs. The wheel sizes vary from 8" to 12" for durability and are great for light use and getting from one place to another which makes getting around town easier.



MDS808200B  
p. 41



### Basic Steel Transport Chairs with 8" Wheels

- Durable steel frame with black upholstery
- Includes detachable footrests and a seat belt for safety
- Back folds down for convenient transporting

| Item No.                | Description | Wt. Capacity      | Leg Style  | Arm Style              | Seat Width | HCPCS | Retail Pkg. | Pkg. |
|-------------------------|-------------|-------------------|------------|------------------------|------------|-------|-------------|------|
| MDS808200BN             | Black       | 250 lb (113.4 kg) | Swing-away | Full-Length, Permanent | 17"        | -     | N           | 1 ea |
| MDS808200B              | Black       | 250 lb (113.4 kg) | Swing-away | Full-Length, Permanent | 19"        | E0138 | N           | 1 ea |
| MDS808200WAT            | Red         | 300 lb (136.1 kg) | Swing-away | Desk-Length, Permanent | 19"        | -     | N           | 1 ea |
| MDS808200               | Chrome      | 300 lb (136.1 kg) | Swing-away | Full-Length, Permanent | 19"        | E0138 | N           | 1 ea |
| MDS808210KDR            | Burgundy    | 300 lb (136.1 kg) | Swing-away | Desk-Length, Permanent | 19"        | -     | Y           | 1 ea |
| <b>New</b> MDS808200KD  | Red         | 300 lb (136.1 kg) | Swing-away | Desk-Length, Permanent | 19"        | -     | Y           | 1 ea |
| <b>New</b> MDS808200KDT | Teal        | 300 lb (136.1 kg) | Swing-away | Full-Length, Permanent | 19"        | -     | Y           | 1 ea |

### Aluminum Transport Chairs with 8" Wheels

- Powder-coated aluminum frame with black upholstery
- Includes seat belt for safety and padded arms

| Item No.     | Description | Wt. Capacity      | Leg Style  | Arm Style              | Seat Width | HCPCS | Retail Pkg. | Pkg. |
|--------------|-------------|-------------------|------------|------------------------|------------|-------|-------------|------|
| MDS808200ABE | Blue        | 300 lb (136.1 kg) | Swing-away | Full-Length, Permanent | 19"        | E1038 | N           | 1 ea |
| MDS808200ARE | Red         | 300 lb (136.1 kg) | Swing-away | Full-Length, Permanent | 19"        | E1038 | N           | 1 ea |



MDS808200KD



MDS808200KDT



MDS808200ABE

## Aluminum Transport Chairs with 12" Wheels

- Powder coated aluminum frame with Black upholstery
- Includes seat belt for safety and padded arms
- Larger back tire for maneuvering through obstacles
- Hand brakes on handles

| Item No.                 | Description | Wt. Capacity      | Leg Style  | Arm Style              | Seat Width | HCPCS | Retail Pkg. | Pkg. |
|--------------------------|-------------|-------------------|------------|------------------------|------------|-------|-------------|------|
| MDS808210ABE             | Blue        | 300 lb (136.1 kg) | Swing-Away | Full-Length, Permanent | 19"        | E1038 | N           | 1 ea |
| MDS808210ARE             | Red         | 300 lb (136.1 kg) | Swing-Away | Full-Length, Permanent | 19"        | E1038 | N           | 1 ea |
| <b>New</b> MDS808210ABEX | Blue        | 300 lb (136.1 kg) | Swing-Away | Full-Length, Permanent | 19"        | -     | N           | 1 ea |

## Excel® Translator Combination Rollator/Transport Chairs

- Durable steel frame with Black upholstery
- Includes easy-to-adjust, push-button, swing-away footrests
- Footrests lock to sides of rollator when not in use
- Height-adjustable push grips and comfortable hand brake

| Item No.     | Description | Wt. Capacity      | Leg Style  | Arm Style              | Seat Width | HCPCS           | Retail Pkg. | Pkg. |
|--------------|-------------|-------------------|------------|------------------------|------------|-----------------|-------------|------|
| MDS808200TRR | Red         | 250 lb (113.4 kg) | Swing-Away | Desk Length, Permanent | 19"        | -               | N           | 1 ea |
| MDS808200TR  | Blue        | 250 lb (113.4 kg) | Swing-Away | Desk Length, Permanent | 19"        | E0143+<br>E0156 | N           | 1 ea |



MDS808210ARE



MDS808200TR

### Ultralight Transport Chairs

- Lightweight aluminum frame with Black upholstery
- Carrying handles and compact design allow for easy and convenient storage and transport

| Item No.      | Description | Wt. Capacity      | Leg Style  | Arm Style              | Seat Width | HCPCS | Retail Pkg. | Pkg. |
|---------------|-------------|-------------------|------------|------------------------|------------|-------|-------------|------|
| MDS808200F3R  | Red         | 300 lb (136.1 kg) | Swing-Away | Desk-Length, Permanent | 19"        | -     | Y           | 1 ea |
| MDS808200F3B  | Blue        | 300 lb (136.1 kg) | Swing-Away | Desk-Length, Permanent | 19"        | -     | Y           | 1 ea |
| MDS808200F3BK | Black       | 300 lb (136.1 kg) | Swing-Away | Desk-Length, Permanent | 19"        | -     | N           | 1 ea |
| MDS808200F3S  | Silver      | 300 lb (136.1 kg) | Swing-Away | Desk-Length, Permanent | 19"        | -     | N           | 1 ea |

### Bariatric Transport Chairs

- Durable frame with Black upholstery
- Seat belt, handbrakes and wheel lock offer better control

| Item No.     | Description                           | Wt. Capacity      | Leg Style  | Arm Style               | Seat Width | HCPCS | Retail Pkg. | Pkg. |
|--------------|---------------------------------------|-------------------|------------|-------------------------|------------|-------|-------------|------|
| TCB228N21S   | Bariatric Steel, Blue                 | 500 lb (226.8 kg) | Swing-Away | Desk Length, Swing Back | 22"        | -     | N           | 1 ea |
| MDS80822B    | Bariatric Steel, Red                  | 500 lb (226.8 kg) | Swing-Away | Desk-Length             | 22"        | E1039 | N           | 1 ea |
| MDS808200BAR | Bariatric, Lightweight, Aluminum, Red | 400 lb (181.4 kg) | Elevating  | Elevating               | 22"        | E1039 | N           | 1 ea |



# Bath Safety

Injury and health conditions can impair mobility and contribute to safety challenges in the bathroom. To help provide security and stability, Medline offers bath safety products that can be used throughout the bathroom.

## Commodes

A bedside commode can be placed next to the bed, making it easy and accessible anytime. Commode liners inside the commode bucket make clean-up easier and sanitary.

### 3-in-1 Commodes

- Use as a bedside commode or over toilet
- Folding commode reduces storage space by 35%

| Item No.   | Description  | Wt. Capacity      | Seat Height | Seat Depth | Inside Arm Width | Overall Width | HCPCS | Retail Pkg. | Pkg. |
|------------|--|-------------------|-------------|------------|------------------|---------------|-------|-------------|------|
| MDS89664   | Grey, Fixed, Standard Seating                      | 350 lb (158.8 kg) | 15"–21"     | 15"        | 19"              | 24"           | E0163 | N           | 4/cs |
| G1-200X1   | White, Folding, Standard Seating                   | 350 lb (158.8 kg) | 16"–22"     | 14.5"      | 19"              | 21"           | E0163 | N           | 1 ea |
| G1-200X4   | White, Folding, Standard Seating                   | 350 lb (158.8 kg) | 16"–22"     | 14.5"      | 19"              | 21"           | E0163 | N           | 4/cs |
| G1-200RX1  | White, Folding, Standard Seating                   | 350 lb (158.8 kg) | 16"–22"     | 14.5"      | 19"              | 21"           | E0163 | Y           | 1 ea |
| G1-200LX1  | White, Folding, Elongated Seating                  | 350 lb (158.8 kg) | 16"–22"     | 16.5"      | 19"              | 21"           | E0163 | N           | 1 ea |
| G1-200LX4  | White, Folding, Elongated Seating                  | 350 lb (158.8 kg) | 16"–22"     | 16.5"      | 19"              | 21"           | E0163 | N           | 4/cs |
| G1-200MX1  | Dark Grey Folding Standard Seating with Microban*  | 350 lb (158.8 kg) | 16"–22"     | 14.5"      | 19"              | 21"           | –     | N           | 1 ea |
| G1-200LMX1 | Dark Grey Folding Elongated Seating with Microban* | 350 lb (158.8 kg) | 16"–22"     | 16.5"      | 19"              | 21"           | –     | N           | 1 ea |

\*These antimicrobial properties are built-in to protect the product. This product does not protect users or others against bacteria, viruses, germs or other disease organisms.



G1-200X1



G1-200LX1

Drop arm holds up to 350 lbs

### Drop Arm Commodes

- Armrests conveniently swing away for easy patient transfer
- Padded options available for increased comfort

| Item No.   | Description                                    | Wt. Capacity      | Seat Height | Seat Depth | Inside Arm Width | Overall Width | HCPCS | Retail Pkg. | Pkg. |
|------------|--|-------------------|-------------|------------|------------------|---------------|-------|-------------|------|
| G1-301DX1  | White, Drop Arm, Elongated                     | 350 lb (158.8 kg) | 19"-23"     | 16.5"      | 20"              | 23"           | E0165 | N           | 1 ea |
| G1-301DX2  | White, Drop Arm, Elongated                     | 350 lb (158.8 kg) | 19"-23"     | 16.5"      | 20"              | 23"           | E0165 | N           | 2/cs |
| G1-301DPX1 | White, Padded, Drop Arm, Elongated             | 350 lb (158.8 kg) | 19"-23"     | 16.5"      | 20"              | 23"           | E0165 | N           | 1 ea |
| G1-303DWX1 | White Drop Arm Elongated, Four Locking Casters | 350 lb (158.8 kg) | 17"-22.5"   | 17.5"      | 16.5"            | 23"           | -     | N           | 1 ea |
| G98204     | Chrome Padded Seat and Back with Drop Arm      | 350 lb (158.8 kg) | 20"-25"     | 18"-20"    | 18.5"            | 22"           | E0165 | N           | 1 ea |

Easy to fold and budget friendly

### Economy Commodes

- Basic commode without seat lid or plastic armrests
- Available with standard or elongated seat

| Item No.   | Description                      | Wt. Capacity      | Seat Height | Seat Depth | Inside Arm Width | Overall Width | HCPCS | Retail Pkg. | Pkg. |
|------------|----------------------------------|-------------------|-------------|------------|------------------|---------------|-------|-------------|------|
| G1-202EX4  | White, Folding Standard Seating  | 350 lb (158.8 kg) | 16"-22"     | 14.5"      | 19"              | 21"           | E0163 | N           | 4/cs |
| G1-202ELX4 | White, Folding Elongated Seating | 350 lb (158.8 kg) | 16"-22"     | 16.5"      | 19"              | 21"           | E0163 | N           | 4/cs |



G1-301DX1



G1-303DWX1

## Bariatric Commodes

- Fixed and folding commodes have sturdy steel frames
- Perfect for high weight capacity needs

| Item No.   | Description                             | Wt. Capacity        | Seat Height | Seat Depth | Inside Arm Width | Overall Width | HCPCS | Retail Pkg. | Pkg. |
|------------|---|---------------------|-------------|------------|------------------|---------------|-------|-------------|------|
| G1-201LBX1 | White, Bariatric Steel Folding Commode  | 650 lb (294.8 kg)   | 16"-22"     | 16.5"      | 24"              | 30.5"         | E0168 | N           | 1 ea |
| G1-400BX2  | White, Bariatric Steel Fixed Bench Seat | 650 lb (294.8 kg)   | 17"-20.5"   | 18.5"      | 24"              | 29.5"         | E0168 | N           | 2 ea |
| G1-400BDX1 | White, Bariatric Steel Drop Arm Commode | 1,000 lb (453.6 kg) | 17.5"-21.5" | 18.5"      | 25"              | 31"           | E0168 | N           | 1 ea |



G1-202ELX4



G1-201LBX1



G1-400BX2



G1-400BDX1

## Aluminum Commodes

- Rust-resistant aluminum
- Multi-functional bathroom usage

| Item No.    | Description   | Wt. Capacity        | Seat Height | Seat Depth | Inside Arm Width | Overall Width | HCPCS | Retail Pkg. | Pkg. |
|-------------|---|---------------------|-------------|------------|------------------|---------------|-------|-------------|------|
| G1-501LWX1  | Aluminum Elongated Commode, with Four Locking Casters           | 300 lb (136.1 kg)   | 21.5"       | 16.5"      | 18"              | 22"           | -     | N           | 1 ea |
| G1-502PWX1  | Aluminum Padded Shower Chair Commode, with Four Locking Casters | 300 lb (136.1 kg)   | 20.5"       | 16"        | 18.5"            | 21"           | E0163 | N           | 1 ea |
| G1-503WCPX1 | Aluminum Shower Commode Wheelchair                              | 275 lb (124.7 kg)   | 20"-22"     | 15.5"      | 19"              | 25"           | -     | N           | 1 ea |
| G1-504WCPX1 | Aluminum Shower Commode Transport Chair                         | 275 lb (124.7 kg)   | 20"-22"     | 15.5"      | 19"              | 23"           | -     | N           | 1 ea |
| G1-505LX1   | Aluminum 3-In-1 Folding Commode                                 | 350 lb (158.8 kg)   | 15.5"-21.5" | 16.5"      | 17"              | 22"           | -     | N           | 1 ea |
| G1-506DWX1  | Aluminum Drop Arm Commode, with Four Locking Casters            | 300 lb (136.1 kg)   | 18.5"-23"   | 16.5"      | 17.5"            | 22"           | -     | N           | 1 ea |
| G1-507BWX1  | Aluminum Bariatric Commode, with Four Locking Casters           | 1,000 lb (453.6 kg) | 21"         | 21"        | 25.5"            | 28.5"         | -     | N           | 1 ea |



G1-502PWX1



G1-503WCPX1



G1-506DWX1



G1-507BWX1

## Buckets

- Replacement commode buckets
- Available with or without lid and handle

| Item No.                | Description  | HCPCS | Retail Pkg. | Pkg.  |
|-------------------------|--|-------|-------------|-------|
| MDS80306B               | Grey Bucket with Lid and Metal Handle                            | E0167 | N           | 6/cs  |
| MDS80306EB              | Grey Bucket Only   | E0167 | N           | 18/cs |
| G1-004X18               | Grey Bucket Lid Only   | -     | N           | 18/cs |
| G1-005X6                | Grey Bucket with Lid and Plastic Handle                          | -     | N           | 6/cs  |
| <b>New</b> G1-400DBBUCK | Grey Bucket with Lid and Metal Handle Replacement for G1-400BDX1 | -     | N           | 1 ea  |

## **New** Commode Liners

- Disposable commode liners make clean up easy

| Item No.      | Description  | HCPCS | Retail Pkg. | Pkg.       |
|---------------|--|-------|-------------|------------|
| MDS89664LINER | Commode Liner with Absorbent Pad, 12/roll                            | -     | Y           | 6 rolls/cs |
| MDS89664LINEH | Commode Liner with Absorbent Pad, 12/roll                            | -     | Y           | 1 roll     |
| G1-003X6      | Commode Liner with Absorbent Pad, 12/roll, with Box Dispenser Cutout | -     | Y           | 6 rolls/cs |
| G1-001X1      | Commode Liner Wall Mount for G1-003X6                                | -     | N           | 1 ea       |
| G1-002X1      | Commode Liners, No Absorbent Pad, 50/roll                            | -     | N           | 1 roll     |



MDS80306B



G1-001X1



MDS89664LINER



# Shower Chairs

Shower chairs and transfer benches can make getting in and out of the tub or shower easier and more convenient. Shower chairs are placed inside the tub or shower, while transfer benches straddle the side of a tub, allowing one to sit outside the tub then move over the surface of the bench into the tub.

## Shower Chairs with Back

- Rust-resistant aluminum chairs are perfect for patients who have difficulty standing in the shower
- Backrest included for added stability, security and comfort

| Item No.             | Description  | Wt. Capacity      | Seat Height   | Seat Depth | HCPCS | Retail Pkg. | Pkg. |
|----------------------|--|-------------------|---------------|------------|-------|-------------|------|
| G2-100KHX1           | Grey Shower Chair with Arms                                | 350 lb (158.8 kg) | 16"-21"       | 13.5"      | -     | N           | 1 ea |
| MDS89745RAH          | White Shower Chair with Arms                               | 350 lb (158.8 kg) | 16"-21"       | 13.5"      | E2045 | Y           | 1 ea |
| MDS89745ELMBH        | Tan Shower Chair with Arms, with Microban*                 | 350 lb (158.8 kg) | 15"-21"       | 15"        | -     | Y           | 1 ea |
| MDS89745KDMBH        | Light Blue Shower Chair, With Microban*                    | 300 lb (136.1 kg) | 13.75"-20.75" | 12"        | -     | Y           | 1 ea |
| G2-101KX1            | Grey Tool-Free Shower Chair                                | 400 lb (181.4 kg) | 13.75"-20.5"  | 12"        | E0240 | N           | 1 ea |
| G2-101KRX1           | Grey Retail Shower Chair, Tool-Free                        | 400 lb (181.4 kg) | 13.75"-20.5"  | 12"        | E0240 | Y           | 1 ea |
| G2-102KX1            | Grey Knockdown Shower Chair                                | 400 lb (181.4 kg) | 15"-19"       | 12"        | E0240 | N           | 1 ea |
| G2-102KX4            | Grey Knockdown Shower Chair                                | 400 lb (181.4 kg) | 15"-19"       | 12"        | E0240 | N           | 4/cs |
| G2-103RX1            | Grey Folding Tool-Free Shower Chair with Back              | 300 lb (136.1 kg) | 14.5"-20.5"   | 15"        | -     | Y           | 1 ea |
| G2-104AX2            | Grey Deluxe Aluminum Shower Chair with Back                | 400 lb (181.4 kg) | 13.5"-20.5"   | 12"        | -     | N           | 2/cs |
| G2-105RX1            | White Premium Plastic Shower Chair with Back               | 350 lb (158.8 kg) | 16"-20.5"     | 15"        | -     | Y           | 1 ea |
| <b>New</b> G2-106KX1 | Grey Knockdown Shower Chair with Cut Out for Easy Cleaning | 300 lb (136.1 kg) | 16.5"-21.5"   | 14.5"      | -     | N           | 1 ea |
| MDSMOMCHAIRG         | Dark Grey Momentum Shower Chair                            | 250 lb (113.4 kg) | 16"-20.5"     | 15"        | -     | Y           | 2/cs |
| MDSMOMCHAIRGH        | Dark Grey Momentum Shower Chair                            | 250 lb (113.4 kg) | 16"-20.5"     | 15"        | -     | Y           | 1 ea |
| G2-111KMX1           | White Shower Chair with Microban*                          | 400 lb (181.4 kg) | 16.5"-20.5"   | 14"        | -     | N           | 1 ea |
| G2-112KMX1           | Blue Shower Chair with Microban*                           | 400 lb (181.4 kg) | 16.5"-20.5"   | 14"        | -     | N           | 1 ea |

\*These antimicrobial properties are built-in to protect the product. This product does not protect users or others against bacteria, viruses, germs or other disease organisms.



G2-100KHX1



G2-102KX1



G2-105RX1

### Shower Chairs without Back

- Rust-resistant aluminum shower chairs for patients who have difficulty standing in the shower
- No backrest for easy maneuverability in tight spaces

| Item No.     | Description                         | Wt. Capacity       | Seat Height  | Seat Depth | HPCPS | Retail Pkg. | Pkg. |
|--------------|-------------------------------------|--------------------|--------------|------------|-------|-------------|------|
| G2-200KH1    | Grey Shower Chair with Arms         | 350 lb (158.8 kg)  | 16"-21"      | 13.5"      | -     | N           | 1 ea |
| MDS89740RWAH | White Shower Chair with Arms        | 350 lb (158.8 kg)  | 16"-21"      | 13.5"      | -     | Y           | 1 ea |
| G2-201KR1    | Grey Retail Shower Chair, Tool-Free | 400 lb (181.4 kg)  | 13.75"-20.5" | 12"        | E0245 | Y           | 1 ea |
| G2-202KX1    | Grey Shower Chair, Knockdown        | 400 lb (181.4 kg)  | 15"-19"      | 12"        | E0245 | N           | 1 ea |
| G2-203FR1    | Grey Folding Shower Chair           | 250 lb. (113.4 kg) | 14.5"-16"    | 12"        | E0245 | Y           | 1 ea |
| G2-204AX2    | Grey Assembled Shower Chair         | 400 lb (181.4 kg)  | 13.5"-20.5"  | 12"        | -     | N           | 2/cs |
| G2-300KR1    | Grey Bath Stool                     | 300 lb (136.1 kg)  | 13.75"-20.5" | 12.5"      | E0245 | Y           | 1 ea |
| G2-300KRP1   | Grey Padded Bath Stool with Swivel  | 300 lb (136.1 kg)  | 16"-22"      | 13.75"     | E0245 | Y           | 1 ea |
| G2-400RX1    | Grey Tub Portable Bath Bench        | 250 lb (113.4 kg)  | 22"          | 12.5"      | E0245 | Y           | 1 ea |



G2-400RX1



G2-300KR1

Great for small spaces



G2-200KH1



G2-201KR1



G3-300KRP1

## Bariatric Shower Chairs

- Heavy-duty reinforced frames
- For high weight capacity needs

| Item No.   | Description                           | Wt. Capacity      | Seat Height | Seat Depth | HCPSC | Retail Pkg. | Pkg. |
|------------|---------------------------------------|-------------------|-------------|------------|-------|-------------|------|
| G2-100BAX1 | Grey Assembled Shower Chair with Back | 650 lb (294.8 kg) | 16.5"–20.5" | 17"        | E0245 | N           | 1 ea |
| G2-102BX1  | Grey Shower Chair with Back           | 550 lb (249.5 kg) | 15"–20"     | 12"        | E0240 | N           | 1 ea |
| G2-202BX1  | Grey Shower Chair without Back        | 550 lb (249.5 kg) | 15"–20"     | 12"        | E0245 | N           | 1 ea |

## Transfer Benches

- Allows patients to enter the bath tub from a seated position
- Padded and bariatric options available to meet various needs

| Item No.              | Description                         | Wt. Capacity      | Seat Height  | Seat Depth | HCPSC | Retail Pkg. | Pkg. |
|-----------------------|-------------------------------------|-------------------|--------------|------------|-------|-------------|------|
| G3-100KBX1            | Grey Transfer Bench                 | 400 lb (181.4 kg) | 18"–22.5"    | 16.25"     | E0247 | N           | 1 ea |
| G3-100KBX2            | Grey Transfer Bench                 | 400 lb (181.4 kg) | 18"–22.5"    | 16.25"     | E0247 | N           | 2/cs |
| G3-100KBRX1           | Grey Transfer Bench                 | 400 lb (181.4 kg) | 18"–22.5"    | 16.25"     | E0247 | Y           | 1 ea |
| G3-101KBX1            | Grey Transfer Bench Commode         | 400 lb (181.4 kg) | 18"–22"      | 16"        | E0247 | N           | 1 ea |
| G3-101KBPX1           | Grey Padded Transfer Bench Commode  | 400 lb (181.4 kg) | 18"–22"      | 16"        | E0247 | N           | 1 ea |
| G3-200KBX1            | Grey Bariatric Transfer Bench       | 500 lb (226.8 kg) | 19.25"–23.5" | 16.5"      | E0248 | N           | 1 ea |
| <b>New</b> G3-300BAX1 | Grey Bariatric Transfer Bench       | 600 lb (272.2 kg) | 20"–24"      | 16.5"      | –     | N           | 1 ea |
| <b>New</b> G3-500KRX1 | White Composite Transfer Bench      | 350 lb (158.8 kg) | 17"–21"      | 17.5"      | –     | Y           |      |
| G3-001RX1             | White Transfer Bench Shower Curtain | N/A               | N/A          | N/A        | –     | Y           | 1 ea |



G2-102BX1



G3-100KBX1

**400-lb weight cap at a fraction of the price**



G3-001RX1

**Great add-on item for your Transfer Bench**

## Raised Toilet Seats

A raised toilet seat installs on top of or in place of the toilet seat, decreasing the distance to the seat when sitting down or standing up.

### New Economy Raised Toilet Seats

- Makes using the toilet easier for those with trouble bending
- Simple design made of polypropylene resin is easy to clean
- Slip-resistant pads and lip inside toilet offer additional stability

| Item No.    | Description   | Wt. Capacity      | HCPCS | Retail Pkg. | Pkg. |
|-------------|---|-------------------|-------|-------------|------|
| G4-111MX1   | Elevated Locking Toilet Seat, with Microban*, 5", White           | 400 lb (181.4 kg) | -     | N           | 1 ea |
| G4-111MRX1  | Elevated Locking Toilet Seat, with Microban*, 5", White           | 400 lb (181.4 kg) | -     | Y           | 1 ea |
| G4-112MX1   | Elevated Locking Toilet Seat with Arms, with Microban*, 5", Blue  | 400 lb (181.4 kg) | -     | N           | 1 ea |
| G4-115MRX1  | Elevated Locking Toilet Seat with Arms, with Microban*, 5", White | 400 lb (181.4 kg) | -     | Y           | 1 ea |
| G4-201PX1   | Padded Elevated Toilet Seat, White                                | -                 | -     | N           | 1 ea |
| G4-301LRX1  | 3.5" Seat Riser with Arms That Can Fit Elongated Toilets          | 400 lb (181.4 kg) | -     | Y           | 1 ea |
| G4-301RX1   | 3.5" Seat Riser with Arms That Can Fit Standard Toilets           | 400 lb (181.4 kg) | -     | Y           | 1 ea |
| G4-300LRX1  | 3.5" Hinged Toilet Seat Riser, Elongated                          | 400 lb (181.4 kg) | -     | Y           | 1 ea |
| G4-300RX1   | 3.5" Hinged Toilet Seat Riser                                     | 400 lb (181.4 kg) | -     | Y           | 1 ea |
| G4-300BRX1  | Bariatric Seat Riser without Arms                                 | 600 lb (272.2 kg) | -     | Y           | 1 ea |
| G4-300BHRX1 | Bariatric Seat Riser with Arms                                    | 600 lb (272.2 kg) | -     | Y           | 1 ea |
| G4-400RX1   | Raised Toilet Seat with Side Locks, 4", White                     | 400 lb (181.4 kg) | -     | Y           | 1 ea |
| G4-500RX1   | Toilet Seat Riser without Arms, without Lock, 5", White           | 400 lb (181.4 kg) | -     | Y           | 1 ea |
| G4-501RX1   | Toilet Seat Riser without Arms, with Lock, 5", White              | 400 lb (181.4 kg) | -     | Y           | 3/cs |
| G4-502RX1   | Elevated Locking Toilet Seat with Arms, 5", White                 | 400 lb (181.4 kg) | -     | Y           | 3/cs |
| G4-503LRX1  | Elongated Elevated Locking Toilet Seat with Arms, 5", White       | 400 lb (181.4 kg) | -     | Y           | 3/cs |
| G7-200RX1   | Toilet Safety Light   | -                 | -     | Y           | 1 ea |

\*These antimicrobial properties are built-in to protect the product. This product does not protect users or others against bacteria, viruses, germs or other disease organisms.



G4-301LRX1



G4-300RX1



G4-300BHRX1



G4-400RX1



G7-200RX1



G4-503LRX1

## Grab Bars

Grab bars allow safe mobility to perform everyday tasks that might otherwise require assistance in or around showers and tubs, making bathing safer.

### New Toilet Safety Rails

| Item No.   | Description  | Wt. Capacity      | Seat Height | HCPCS | Retail Pkg. | Pkg. |
|------------|--|-------------------|-------------|-------|-------------|------|
| G5-100RX1  | Aluminum Foldable Toilet Safety Rails                      | 250 lb (113.4 kg) | 27.5"-32"   | -     | Y           | 1 ea |
| G5-200X1   | Aluminum Foldable Toilet Safety Rails with Padded Armrests | 300 lb (136.1 kg) | 26"-31"     | -     | N           | 1 ea |
| MDS86100RF | Aluminum Foldable Toilet Safety Rails                      | 250 lb (113.4 kg) | 27.5"-32"   | E0243 | Y           | 1 ea |
| G3-0300H   | Aluminum Toilet Safety Rails with Padded Armrests          | 300 lb (136.1 kg) | 26"-31"     | -     | N           | 1 ea |

### New White Enamel Grab Bars

- Enamel-coated steel grab bars mount to the wall
- Can be used horizontally or vertically and include mounting hardware

| Item No.   | Description | Wt. Capacity      | HCPCS | Retail Pkg. | Pkg. |
|------------|-------------|-------------------|-------|-------------|------|
| G5-012WRX1 | 12" Wide    | 300 lb (136.1 kg) | E0241 | Y           | 1 ea |
| G5-016WRX1 | 12" Wide    | 300 lb (136.1 kg) | E0241 | Y           | 1 ea |
| G5-018WRX1 | 12" Wide    | 300 lb (136.1 kg) | E0241 | Y           | 1 ea |
| G5-024WRX1 | 16" Wide    | 300 lb (136.1 kg) | E0241 | Y           | 1 ea |
| G5-032WRX1 | 16" Wide    | 300 lb (136.1 kg) | E0241 | Y           | 1 ea |



G5-012WRX1

### New Knurled Chrome Grab Bars

- Chrome plated steel grab bars mount to the wall
- Can be used horizontally or vertically and include mounting hardware

| Item No.   | Description | Wt. Capacity      | HCPCS | Retail Pkg. | Pkg. |
|------------|-------------|-------------------|-------|-------------|------|
| G5-012CRX1 | 12" Wide    | 300 lb (136.1 kg) | -     | Y           | 1 ea |
| G5-016CRX1 | 12" Wide    | 300 lb (136.1 kg) | -     | Y           | 1 ea |
| G5-018CRX1 | 12" Wide    | 300 lb (136.1 kg) | -     | Y           | 1 ea |
| G5-024CRX1 | 16" Wide    | 300 lb (136.1 kg) | -     | Y           | 1 ea |
| G5-032CRX1 | 16" Wide    | 300 lb (136.1 kg) | -     | Y           | 1 ea |



G5-012CRX1

### New Suction Grab Bars

- Attach to the wall without screws
- Grab bars stay in place with suction cups
- White



| Item No.   | Description                   | Wt. Capacity      | HCPCS | Retail Pkg. | Pkg. |
|------------|-------------------------------|-------------------|-------|-------------|------|
| G5-012SRX1 | With Locking Indicator, 12" L | 250 lb (113.4 kg) | -     | Y           | 1 ea |
| G5-016SRX1 | With Locking Indicator, 16" L | 250 lb (113.4 kg) | -     | Y           | 1 ea |
| G5-019SRX1 | Pivot Grab Bar, 19" L         | 250 lb (113.4 kg) | -     | Y           | 1 ea |

### Tub Grab Bars

- Step-through clamp design for better hand positioning
- Tool-free installation with lined pads that won't mar tub surface
- White

| Item No.              | Description                                    | Wt. Capacity      | HCPCS | Retail Pkg. | Pkg. |
|-----------------------|--|-------------------|-------|-------------|------|
| <b>New</b> G5-300RX1  | Steel, Locks to Side of Tub                    | 250 lb (113.4 kg) | -     | Y           | 1 ea |
| <b>New</b> G5-301RX1  | Composite, Locks to Side of Tub                | 250 lb (113.4 kg) | -     | Y           | 1 ea |
| <b>New</b> G5-301MRX1 | Composite with Microban*, Locks to Side of Tub | 250 lb (113.4 kg) | -     | Y           | 1 ea |
| <b>New</b> G5-302HRX1 | Steel, Height Adjustable, Locks to Side of Tub | 300 lb (136.1 kg) | -     | Y           | 1 ea |
| MDS86321W             | Steel, Locks to Side of Tub                    | 250 lb (113.4 kg) | E0246 | Y           | 2/cs |
| MDS86321WH            | Steel, Locks to Side of Tub                    | 250 lb (113.4 kg) | E0246 | Y           | 1 ea |



\*These antimicrobial properties are built-in to protect the product. This product does not protect users or others against bacteria, viruses, germs or other disease organisms.

## Accessories

### Cast Protectors

- Help keep casts dry when bathing
- Available in adult and child sizes
- Clear

| Item No.   | Description | HCPCS | Retail Pkg. | Pkg. |
|------------|-------------|-------|-------------|------|
| CUR200ALL  | Adult Leg   | -     | Y           | 6/cs |
| CUR200ALLH | Adult Leg   | -     | Y           | 1 ea |
| CUR200AAA  | Adult Arm   | -     | Y           | 6/cs |
| CUR200AAAH | Adult Arm   | -     | Y           | 1 ea |
| CUR100KLL  | Child Leg   | -     | Y           | 6/cs |
| CUR100KLLH | Child Leg   | -     | Y           | 1 ea |
| CUR100KAA  | Child Arm   | -     | Y           | 6/cs |
| CUR100KAAH | Child Arm   | -     | Y           | 1 ea |



CUR200ALL



CUR200AAA

## Shower Heads

- Easy to install White handheld shower heads with diverter valve
- Two styles

| Item No.    | Description                                  | HCPCS | Retail Pkg. | Pkg. |
|-------------|--|-------|-------------|------|
| G7-100RX1   | Deluxe Hand Held Shower, 3 Massage Functions | -     | Y           | 1 ea |
| MDSHSHOWER  | White Hand Held with 6.5 ft Hose             | -     | Y           | 6/cs |
| MDSHSHOWERH | White Hand Held with 6.5 ft Hose             | -     | Y           | 1 ea |



## Momentum Shower Caddies

- Basket helps organize shampoo, conditioner and soap
- Supports up to 5 lbs

| Item No.      | Description               | HCPCS | Retail Pkg. | Pkg. |
|---------------|---------------------------|-------|-------------|------|
| MDSMOMCADDYG  | Wall Mounted Shower Caddy | -     | Y           | 2/cs |
| MDSMOMCADDYGH | Wall Mounted Shower Caddy | -     | Y           | 1 ea |



## Bath Mats

- Slip-resistant heavyweight mat with suction cups on bottom for secure grip
- Antimicrobial\* treatment helps protect bath mat from growth of stains and odors caused by bacteria, mold and mildew

| Item No.    | Description   | HCPCS | Retail Pkg. | Pkg. |
|-------------|---|-------|-------------|------|
| MDSBATHMAT  | Bath Mat with Microban*, White                        | -     | Y           | 4/cs |
| MDSBATHMATH | Bath Mat with Microban*, White                        | -     | Y           | 1 ea |
| MDSMOMMATG  | Momentum Bath Mat with Microban*, 29.5" x 16.8", Grey | -     | Y           | 2/cs |
| MDSMOMMATGH | Momentum Bath Mat with Microban*, 29.5" x 16.8", Grey | -     | Y           | 1 ea |

\*These antimicrobial properties are built-in to protect the product. This product does not protect users or others against bacteria, viruses, germs or other disease organisms.



## Bath Furniture

- Teak look without the constant maintenance
- Frames are aluminum seat and shelf are composite wood

| Item No. | Description                          | Wt. Capacity     | HCPCS | Retail Pkg. | Pkg. |
|----------|--------------------------------------|------------------|-------|-------------|------|
| G2-500X1 | Bathroom Bench with Faux Wood        | 200 lb (90.7 kg) | -     | N           | 1 ea |
| G2-501X1 | Bathroom Corner Bench with Faux Wood | 200 lb (90.7 kg) | -     | N           | 1 ea |



# Home Care Beds and Accessories

Medline offers high-quality, cost-effective homecare beds and accessories for the convenience and satisfaction of patients and caregivers alike.

Meeting or exceeding very stringent dynamic and static loading requirements, each bed is inspected, fully system tested, and approved at maximum weight capacities. Recessed headspring supports and foot assembly reduce the chance of pinching and injuries. With some of the industries lightest beds; transport, setup and storage is easier.

Our beds can be cleaned without having to remove the motor and end panels are easy to wipe down and keep clean. To complete your Medline homecare bed, we offer high-performance mattresses and easy-to-install, durable side rails. Patient mobility, transfers, and movement are all possible with our easy-to-assemble trapezes, bed assist bars, and reachers.

## Basic Beds

- Headspring supports recessed to resist pinching and injuries
- Foot assembly recessed to reduce chance of pinching
- 2-year warranty on electronics, limited lifetime warranty on the frame and welds

| Item No.      | Description                     | Wt. Capacity      | Height Range | HCPCS | Retail Pkg. | Pkg. |
|---------------|---------------------------------|-------------------|--------------|-------|-------------|------|
| MDR107002E    | Semi-Electric Basic Bed         | 450 lb (204.1 kg) | 15"-20"      | E0295 | N           | 1 ea |
| MDR107002E-4* | Semi-Electric Basic Bed, 4-Pack | 450 lb (204.1 kg) | 15"-20"      | E0295 | N           | 1 ea |
| MDR107003E    | Full-Electric Basic Bed         | 450 lb (204.1 kg) | 15"-20"      | E0297 | N           | 1 ea |
| MDR107003E-4* | Full-Electric Basic Bed, 4-Pack | 450 lb (204.1 kg) | 15"-20"      | E0297 | N           | 1 ea |
| MDR107003ELO  | Full-Electric Basic Low Bed     | 450 lb (204.1 kg) | 9.5"-22.5"   | E0297 | N           | 1 ea |

\*Order in multiples of 4.



MDR107003E



## Medlite Beds

- Premium lightweight home care beds
- 5-year warranty on electronics, limited lifetime warranty on the frame and welds (all Medlite beds)

| Item No.    | Description                   | Wt. Capacity      | Height Range | HCPCS | Retail Pkg. | Pkg. |
|-------------|-------------------------------|-------------------|--------------|-------|-------------|------|
| MDR107002L  | Medlite Semi-Electric Bed     | 450 lb (204.1 kg) | 14.5"-22.5"  | E0295 | N           | 1 ea |
| MDR107003L  | Medlite Full-Electric Bed     | 450 lb (204.1 kg) | 14.5"-22.5"  | E0297 | N           | 1 ea |
| MDR107003LO | Medlite Full-Electric Low Bed | 450 lb (204.1 kg) | 9.5"-22.5"   | E0297 | N           | 1 ea |

## Bariatric Bed

- Powerful motors and sturdy steel frame for heavy-duty homecare bed needs
- 2-year warranty on electronics, limited lifetime warranty on the frame and welds

| Item No.  | Description                           | Wt. Capacity      | Height Range | HCPCS | Retail Pkg. | Pkg. |
|-----------|---------------------------------------|-------------------|--------------|-------|-------------|------|
| MDR107004 | Full-Electric Bariatric Home Care Bed | 600 lb (272.2 kg) | 18"-26.5"    | E0301 | N           | 1 ea |



## Side Rails

- Half- and full-length Medline home care bed side rails
- 1-year warranty


| Item No.    | Description                | HCPCS | Retail Pkg. | Pkg. |
|-------------|----------------------------|-------|-------------|------|
| MDS89697    | Clamp-on Half Rail         | E0305 | N           | 1 pr |
| MDS89698N   | Spring Loaded Half Rail    | E0305 | N           | 1 pr |
| MDS89694N   | Spring Loaded Full Rail    | E0310 | N           | 1 pr |
| MDS89696N   | Home-Style Adjustable Rail | -     | N           | 1 pr |
| MDS89695N   | Economy Full Rail          | E0310 | N           | 1 pr |
| MDS89697BAR | Bariatric Half Rail        | E0305 | N           | 1 pr |



### Bed Accessories

- Low bed conversion kits, extension kits and transport aids for Medline home care beds

| Item No.   | Description                           | HCPCS | Retail Pkg. | Pkg. |
|------------|---------------------------------------|-------|-------------|------|
| MDSBEDCART | Bed Cart Kit for Home Care Beds       | -     | N           | 1 ea |
| MDR630EXTN | 4" Extension Kit for Home Care Beds   | -     | N           | 1 ea |
| MDSLWBASIC | Low Bed Conversion Kit for MDR107003E | -     | N           | 1 ea |
| MDSLWBED   | Low Bed Conversion Kit for MDR107003L | -     | N           | 1 ea |



MDSBEDCART



MDSLWBED

### Reacher

- Allows user to reach items both high and low

| Item No.    | Description       | HCPCS | Retail Pkg. | Pkg. |
|-------------|-------------------|-------|-------------|------|
| MDS86031RCE | 31" Reacher, Blue | -     | N           | 6/cs |



MDS86031RCE

### Bed Assist Bars

- Provides assistance getting in and out of bed
- Multiple hand positioning options to meet user's preference

| Item No.      | Description            | Handle material | HCPCS | Retail Pkg. | Pkg. |
|---------------|------------------------|-----------------|-------|-------------|------|
| MDS6800BA     | Bed Assist Bar, Black  | Foam            | E0315 | Y           | 3/cs |
| MDS6800BAM    | M-Shape Handle, White  | Foam            | E0315 | N           | 1 ea |
| MDS6800SWIV   | SwivAssist Bar, Black  | Foam            | E0315 | N           | 3/cs |
| MDS6800SWIVRH | SwivAssist Bar, Retail | Foam            | E0315 | Y           | 1 ea |



MDS6800SWIV



MDS6800BAM



MDS6800BA

## Mattresses

- Cost-effective homecare bed mattresses available in three different materials
- Durable vinyl covers

| Item No.   | Description                 | Wt. Capacity      | HCPCS | Retail Pkg. | Pkg. |
|------------|-----------------------------|-------------------|-------|-------------|------|
| MDR237827  | Innerspring, 36" x 80" x 6" | 250 lb (113.4 kg) | E0271 | N           | 1 ea |
| MDR237827E | Fiber, 36" x 80" x 6"       | 250 lb (113.4 kg) | E0272 | N           | 1 ea |
| MDR230981R | Foam, 36" x 80" x 6"        | 350 lb (158.8 kg) | E0272 | N           | 1 ea |
| MDR23784   | Innerspring, 36" x 84" x 6" | 250 lb (113.4 kg) | E0271 | N           | 1 ea |
| MDR23784E  | Fiber, 36" x 84" x 6"       | 250 lb (113.4 kg) | -     | N           | 1 ea |
| MDR211892R | Foam, 42" x 84" x 6"        | 600 lb (272.2 kg) | E0272 | N           | 1 ea |



MDR237827



MDR237827E



MDR230981R

## Overbed Tables

- Standard overbed tables with various features and material options
- Overbed tables for medical use

| Item No.  | Description                 | Wt. Capacity    | Height Range | HCPCS | Retail Pkg. | Pkg. |
|-----------|-----------------------------|-----------------|--------------|-------|-------------|------|
| MDS104015 | Economy Overbed Table       | 50 lb (22.7 kg) | 30"-42"      | E0274 | N           | 1 ea |
| MDS107015 | Composite Top Overbed Table | 50 lb (22.7 kg) | 30"-42"      | E0274 | N           | 1 ea |
| MDS108015 | Fold-Flat Overbed Table     | 50 lb (22.7 kg) | 29.5"-41"    | E0274 | N           | 1 ea |
| MDS104950 | Tilt-Top Overbed Table      | 50 lb (22.7 kg) | 30"-42"      | E0274 | N           | 1 ea |



MDS107015

## Trapeze Parts

- Assist with patient transfers, movement and positioning
- Free-standing and bed-mounted options

| Item No.   | Description                   | Wt. Capacity        | HCPCS | Retail Pkg. | Pkg. |
|------------|-------------------------------|---------------------|-------|-------------|------|
| MDS80615T  | Standard Trapeze Bar          | 250 lb (113.4 kg)   | E0910 | N           | 1 ea |
| MDS80616T  | Standard Trapeze Base         | 250 lb (113.4 kg)   | A9900 | N           | 1 ea |
| MDS80618T  | Economy Trapeze Bar           | 250 lb (113.4 kg)   | E0910 | N           | 1 ea |
| MDS80619T  | Economy Trapeze Base          | 250 lb (113.4 kg)   | A9900 | N           | 1 ea |
| MDS500TPZ  | Lightweight Bariatric Trapeze | 500 lb (226.8 kg)   | E0912 | N           | 1 ea |
| MDS1000TPZ | Bariatric Trapeze             | 1,000 lb (453.6 kg) | E0912 | N           | 1 ea |




MDS1000TPZ

# Patient Lifts

## Elevate resident and caregiver safety

The Association of Safe Patient Handling Professionals (ASPHP) recommends one mobile lift for every eight residents who need one. Our lifts can help you reduce risks of caregiver injuries and lift residents with exceptional safety, convenience and dependability.



## Mobile Lifts

Easy-to-use and reliable mobile lifts for all essential clinical tasks, including resident transfers and repositioning

### Manual Hydraulic Patient Lifts

- Elevating hydraulic lift for patient transportation
- Base opens with an easy hand lever
- 6-point cradle type

| Item No.    | Description                    | Wt. Cap. | Pkg. |
|-------------|--------------------------------|----------|------|
| MDS88200D   | 6-Point Cradle, Includes Sling | 400 lb   | 1 ea |
| MDS88200DLB | Low Base, Includes Sling       | 450 lb   | 1 ea |
| MDS88200HL  | 6-Point Cradle, No Sling       | 600 lb   | 1 ea |



MDS88200D

### Electric Total Mobile Lifts

- Can be lowered to floor height for lifting a resident who has fallen
- Compatible with lift scales for convenient weighing of residents
- State-of-the-art control panel displays battery charge level and more
- Removable battery pack to reduce downtime with optional off-board charging
- Warranty: three-years on all mechanical and electrical parts, one-year on the battery

| Item No.    | Description                                    | Wt. Cap. | Pkg. |
|-------------|--|----------|------|
| MDS400ELB   | Powered Lift with Manual Base                  | 400 lb   | 1 ea |
| MDS450EL    | Powered Lift with Powered Base                 | 450 lb   | 1 ea |
| MDS600ELB   | Powered Lift with Manual Base                  | 600 lb   | 1 ea |
| MDS700EL    | Powered Lift with Powered Base                 | 700 lb   | 1 ea |
| MDSMELSCALE | Electric Scale Accessory for all Total Lifts   | 700 lb   | 1 ea |
| MDS700ARMS  | Gait Training Arms Accessory for MDS700EL Only | —        | 1 ea |
| MDSSELSCALE | Scale Accessory for all Total Lifts            | 600 lb   | 1 ea |



MDS700EL  
(Shown with optional gait training arms)

### Patient Lift Clip Cradle

- Enables caregivers to position residents with one hand, providing better control and helping prevent injuries
- Helps to eliminate grabbing and pulling
- Makes the caregiver's job easier
- Compatible with Medline mobile powered patient lifts MDS450EL and MDS700EL


| Item No.   | Description                                     | Pkg. |
|------------|---|------|
| MDS700MCRD | Patient Lift Clip Cradle (use with Clip Slings) | 1 ea |



MDS700MCRD

## Stand Assist Lifts

For transferring activities with residents who have some ability to bear weight or to provide support for residents to stand and practice bearing weight.



MDSS00SA

### Electric Stand Assist Lifts

- Adjustable knee pad for finding the perfect fit
- State-of-the-art control panel displays battery charge level and more
- Removable battery pack to reduce downtime with optional off-board charging
- Warranty: three-years on all mechanical and electrical parts, one-year on the battery

| Item No.  | Description                                | Wt. Cap. | Pkg. |
|-----------|--|----------|------|
| MDS400SAB | Powered Stand Assist with Manual Base      | 400 lb   | 1 ea |
| MDS500SA  | Powered Stand Assist with Powered Base     | 500 lb   | 1 ea |
| MDS600SAB | Powered Stand Assist Lift with Manual Base | 600 lb   | 1 ea |



STA204

### Manual Standing Transfer Aids

- Non-powered transfer device allows a single caregiver to assist residents safely
- Enables residents to support themselves to a standing position
- Features folding seats, dual knee pads, foot rest and grab bar for user
- Hooks for optional support strap

| Item No. | Description                            | Wt. Cap. | Pkg. |
|----------|--|----------|------|
| STA204   | Resident Transfer Aid, Open/Close Base | 450 lb   | 1 ea |
| STA182   | Resident Transfer Aid, Fixed Base      | 400 lb   | 1 ea |



## Raizer® Mobile Lifting Chair

Now there's a better way to lift someone after a fall.

- Gentle on caregiver—requires virtually no physical effort
- Fast—under a minute, one assistant can set up and raise a large adult from the floor to a standing position
- Easy to carry—two handy bags hold seven lightweight components
- Simple setup—color-coding guides the placement of components, and a loud click confirms each connection
- Easy to clean—cleaning's a snap with the smooth surfaces on all Raizer components

| Item No.   | Description                 | Pkg. |
|------------|-----------------------------|------|
| MDSRAIZER2 | Raizer Mobile Lifting Chair | 1 ea |



### Technical Specifications

|   |  |
|---|--|
| <b>Weight capacity</b>                  | 330 lbs.   |
| <b>Lifting time</b>                     | 20–30 seconds                                    |
| <b>Product weight</b>                   | 28 lbs. total                                    |
| <b>Charger</b>                          | Two pieces: Soniel 12v, 100–240V 12V car adapter |
| <b>Number of lifts on a full charge</b> | 40 at maximum weight capacity                    |
| <b>Charging time</b>                    | Six hours  |
| <b>Warranty</b>                         | Two years from invoice date                      |

### All In 4 Easy Steps:



1 Place the seat motor under the person's knees



2 Attach components to the seat motor



3 Mechanically lift the resident using the console buttons or remote control



4 Transfer the resident to a walker or wheelchair

# Mobility Rack

Free with \$3,500 minimum purchase order of Durable Medical Equipment products

## Recommended planogram:

### Walking Aids

| Item No.      | Description                        | Pkg.    | Page |
|---------------|------------------------------------|---------|------|
| G07850M       | Flip Tray for Guardian Walkers     | 3/cs    | 26   |
| MDS86615K     | Walker Basket                      | 2/cs    | 26   |
| MDS86410KDBW  | Upright Folding Walker, 5" Wheels  | 1/ea    | 23   |
| MDS4PTCANETIP | 4-Point Cane Tips                  | 4/cs    | 11   |
| MDS86425W     | Quad Cane Tips                     | 6/cs    | 11   |
| MDS86426W     | Cane Tips                          | 4/cs    | 11   |
| MDS86222W     | Quad Cane, Small Base, Black       | 2/cs    | 10   |
| MDS86420      | Offset Cane, Black                 | 6/cs    | 9    |
| MDS86228      | Quad Cane, Large Base, Black       | 2/cs    | 10   |
| MDS86430W     | Folding T-Handle Cane              | 6/cs    | 10   |
| 4UCANE        | Folding Hybrid Cane, Offset, Black | 1 ea    | 11   |
| 4UCANEB       | Folding Hybrid Cane, Offset, Blue  | 1 ea    | 11   |
| MDS86850ESKD  | Basic Steel Rollator, Burgundy     | 1 ea    | 12   |
| MDS86615SKI   | Ski Gliders for Walkers            | 6 pr/cs | 26   |
| MDS86615T     | Walker Tips                        | 1 pr    | 26   |
| MDSWBALL      | Walker Ball Glides                 | 6 pr/cs | 26   |
| MDS86615BAG   | 2-Button Walker Bag                | 4/cs    | 26   |
| MDSRUBSKI     | Walker Glides                      | 6/cs    | 26   |
| MDS86615L     | 10" Leg Replacement                | 1 pr    | 26   |
| MDS86615W5    | Footpiece, 5" Wheels               | 1 pr    | 26   |
| MDS86615W5T   | 5" Long Extension Wheels           | 1 pr    | 26   |
| MDS808210KDR  | Lightweight Steel Transport Chair  | 1 ea    | 41   |

### Aids for Daily Living

| Item No.      | Description                             | Pkg. | Page |
|---------------|---|------|------|
| MDS86031RCE   | 31" Reacher, Blue                       | 6/cs | 58   |
| MDS6800BA     | Bed Assist Bar, Black                   | 3/cs | 58   |
| MDS6800SWIV   | Swivassist Swiveling Bed Assistance Bar | 3/cs | 58   |
| MDS6800SWIVRH | Swivassist Bar                          | 1 ea | 58   |
| MDS6800BAM    | Bed Assist Bar with M-Shaped Handle     | 1 ea | 58   |

### Bath Safety

| Item No.    | Description              | Pkg. | Page |
|-------------|--------------------------|------|------|
| G2-101KRX1  | Bath Chair with Back     | 1 ea | 49   |
| G2-201KRX1  | Bath Chair without Back  | 1 ea | 50   |
| G3-100KBRX1 | Transfer Bench with Back | 1 ea | 51   |



MDSMOBRACK

Easy mobility for re-stocking



# Our complete selection will help you every step of the way.

**Medline is your trusted DME partner**, with a wide range of premium products that enhance comfort and mobility, so patients can live their best life.



HCPCS Disclaimer: HCPCS codes and Home Health Consolidated Billing codes provided by Medline are intended as general guidelines only. Medline does not guarantee coverage or reimbursement of any products. You must address all coverage and reimbursement issues (including the correctness and accuracy of codes) with your individual payers.



**We make  
healthcare  
run better™**

**Medline Industries, LP**  
Three Lakes Drive, Northfield, IL 60093  
Medline United States | 1-800-MEDLINE (633-5463)  
medline.com | info@medline.com

**Medline Canada**  
1-800-268-2848 | medline.ca | canada@medline.com  
**Medline México**  
01-800-831-0898 | medlinemexico.com | mexico@medline.com

Follow us    

Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure. © 2022 Medline Industries, LP. All rights reserved. All trademarks, regardless of whether they appear with a trademark symbol, are the property of their respective owners unless otherwise noted. Trademark designations refer to the United States only. MKT19W3805227 / LIT474 / 7M / PP / 43