

# The Innovative Reclining Manual Wheelchair For Kids, Teens & Young Adults





**Recline  
up to 180°**



**Expertly designed to provide a whole new found level of freedom, independence and comfort by combining:**

- ✓ A Range Of Positioning Options
- ✓ Reclining Back With Infinite Adjustments Up To 180°
- ✓ Dual Axle For Seat-To-Floor Height Adjustability
- ✓ Antitipping For Maximum Support & Stability
- ✓ Colored Steel Frame, Lightweight & Foldable Design





- Colored Steel Frame
- Full back and neck support with wings
- Reclining back with infinite adjustments up to 180°
- Back stabilizer bar
- Pelvic belt
- Flip-back desk length arm rests
- Dual axle provides adjustable seat to floor height
- Quick-release mag-style rear wheels
- Swing-away height adjustable elevating legrests
- Anti-tippers - prevents tipping when in recline positions or when tilting chair to go over a curb.



### AVAILABLE COLORS



### ACCESSORY



**3**  
YEAR  
LIMITED  
WARRANTY

SPECIFICATIONS	ZREC1200	ZREC1400	ZREC1600
SEAT WIDTH	12"	14"	16"
SEAT DEPTH	12"	14"	16"
OVERALL OPEN WIDTH	20"	20"	24"
OVERALL FOLDED WIDTH	12"	12"	12"
BACK HEIGHT INCLUDING HEADREST	25.4"	25.4"	30.9"
BACK HEIGHT WITHOUT HEADREST	17.52"	17.52"	19.69"
OVERALL HEIGHT	40.75"	40.75"	49.53"
SEAT TO FLOOR HEIGHT	15"/16"	15"/16"	19.6"/20.6"
OVERALL HEIGHT INCLUDING LEG SUPPORT AND FOOTREST	38.58"	38.58"	52.76"
BACKREST RECLINES	90°-180°	90°-180°	90°-180°
NET WEIGHT	44 LBS	44 LBS	56.75 LBS
CARTON SHIPPING WEIGHT	49.6 LBS	49.6 LBS	62.35 LBS
WEIGHT CAPACITY	150 LBS	2000 LBS	300 LBS



It's time for your child to enjoy the gift of easier, safer and more comfortable movement, indoors and outdoors.

# Dream Big Roam Free

This Time, with

**ziggo**Pro



Powered by Circle Specialty's 25+ years of experience in revolutionary mobility and rehab equipment tailored exclusively to the needs of kids and teens, you have every good reason to make **ziggo**Pro a part of your family.