

## What The Healthcare Community Is Saying About The AgileLife TMS...

*"We are so grateful for the AgileLife, it enables me to keep [my husband] at home. I am anxious to share the message of how to help keep loved ones at home. Because of the AgileLife I am no longer using my strength to lift him in and out of bed and I am no longer worried that we both will fall. - Wife of Veteran*

*"This system makes my life easier. Instead of taking extra time using a mechanical lift and straining my back, with the touch of a button our patient can be moved." -CNA*

*"The AgileLife TMS has been a "game changer" for me. The mechanical lift caused extreme pain and this is so comfortable. I felt a sense of freedom when this came and can't say enough about it." - Client in Skilled Nursing Facility*

*"The AgileLife TMS allowed my client to easily transfer in and out of bed as many times per day as needed/wanted. Overall we were extremely pleased with the functionality of the AgileLife TMS." - Geriatric Care Manager*

*"My mom has been extremely blessed and well-served on every level in her home health care along with the AgileLife. The outdated mechanical lift which caused my mom pain and fear was also the most undignified experience and I would not wish it on my enemy. Mom is now walking independently thanks to the AgileLife."*

*- Client's Daughter*

## PRODUCT SPECIFICATIONS

### Bed Technical Specifications R1.0 | R2.0

|                          |                 |
|--------------------------|-----------------|
| Bed Frame Length         | 93"   91.5"     |
| Bed Frame Width          | 43"   41.5"     |
| Sleep Deck Width         | 36"   36.5"     |
| Sleep Deck Length        | 80"   80"       |
| Sleep Deck Height (low)  | 15.5"   15.5"   |
| Sleep Deck Height (high) | 28"   28"       |
| Weight Capacity          | 500lbs   500lbs |

### Wheelchair Technical Specifications

|                |             |
|----------------|-------------|
| Overall Width  | 28"   33"   |
| Overall Height | 35"   40"   |
| Overall Depth  | 32"   41.5" |
| Seat Width     | 18"   22"   |
| Seat Depth     | 16"   18"   |

### System Technical Specifications

|                     |             |
|---------------------|-------------|
| Overall Length      | 114"   122" |
| Hospital Grade Plug | Yes   Yes   |
| Volts               | 115 VAC     |


### Additional Accessories Offered

- Variety of Therapeutic Mattresses
- Variety of Therapeutic Seat Cushions
- Integrated Commode
- Leg Rests and Leg Support Pad
- High Back Extension
- Arm Trough
- Adjustable Lateral Support
- Seat Belt
- Integrated Bed Alarm

# AgileLife™

## Transfer & Mobility System

## The Only Zero-Lift Transfer & Mobility System



## NEXT HEALTH™

SIMPLY TRANSFORMING CAREGIVING

300 Wilson Avenue | Norwalk, CT 06854  
 Phone: 203-939-1153 | info@nexthealthinc.com  
 www.AgileLifeMobility.com

# AgileLife™

Transfer & Mobility System

The AgileLife System Is  
An Integration Of Several  
Assistive Technologies

## TRANSFER IN APPROXIMATELY 2 MINUTES

Every day individuals, healthcare professionals, and family caregivers face the challenges of safe patient handling and the individual's need for healthy movement and independence.

At **Next Health** we have developed breakthrough assistive technology that provides greater mobility and independence while addressing critical therapeutic and safety needs. Never before has a single solution addressed so many issues across the healthcare continuum.

### The System Includes:

- A Zero-lift transfer system
- A fully powered hospital bed
- A manual wheelchair
- In-bed repositioning
- Transfer to commode (optional)
- Bed alarm (optional)
- In-bed scale
- Therapeutic pressure redistribution mattress and seat cushion
- Alternating air mattress (optional)
- Tilt-in-space wheelchair
- System data logger to measure and record system activity

Assisted Living



Our Veterans



"Aging in Place"



Healthcare Facilities



300 Wilson Avenue | Norwalk, CT 06854  
Phone: 203-939-1153 | [info@nexthealthinc.com](mailto:info@nexthealthinc.com)  
[www.AgileLifeMobility.com](http://www.AgileLifeMobility.com)