THE AQUA WALKER

Item# F-605UW

The Aqua Walker is constructed of durable PVC tubing and stainless steel swivel casters with solid rubber wheels.

This adjustable, waterproof walker aids in independent mobility in the pool.

Armrest cradles provide extra support. The cradle can be attached and removed without the use of tools.





FEATURES:

- Adjustable height range of 32" 45" (81 cm 114 cm)
- 4" removable PVC armrest cradles
- Tool free conversion from standard walker to armrest walker
- 3" solid rubber wheels, two locking swivel casters

PROUDLY SOLD BY: