



Solutions for Daily Living

BATH TRANSFER SYSTEMS

TRANSPORTATION

TOILETING

BATHING

EXERCISER

GAIT TRAINERS



Table of Contents

TRANSPORTATION

Selecting A Car Seat	3
Spirit™ APS™ (Adjustable Positioning System™) Car Seat	4
Accessories for Spirit APS Car Seat	5
Therapedic™ IPS™ (Integrated Positioning System™) Car Seat	6
Accessories for Therapedic IPS Car Seat	7
Therapedic IPS Transport Bases	8

BATHING

Commitment to Bathroom Safety	9
Ultima™ Stainless Steel Bath Chair	10
Accessories for Ultima Stainless Steel Bath Chair	11
Contour™ Supreme™ Reclining Bath Chair	12
Contour Deluxe™ Tilt-in-Space Bath Chair	13
Accessories for Contour Bath Chairs	14
Wrap-Around Bath Supports	15
AquaLift™ Bath Lifter System	16

TOILETING

Positioning Commode	17
Optima™ TSS™ (Toilet & Shower System)	18
Accessories for Optima TSS	19
Toilet Support System	20
Accessories for Toilet Support System	21

BATHROOM TRANSFER SYSTEMS

Elite™ and Omni™ Bath/Shower Transfer System	22
EliteVerse™ and OmniVerse™ Transfer System	23
Versa™ Bath/Shower/Commode Transfer System	23
Elite, Omni, Versa Rolling Chairs and Accessories	24
Dimensions and Specifications for Elite, Omni, and Versa	25
PVC Hi-Back Shower/Commode Chair	26
PVC Hi-Back Foldable Shower/Commode Chair	27

GAIT TRAINERS AND EXERCISER

TheraGait™ PVC Gait Trainer	28
PVC Walker Chairs	29
Rickshaw™ Rehab Exerciser	30
Product Advantages	31

Symbol Legend

Look for these helpful icons throughout our catalog.



SLIDER SAFE™ Auto-Stop Safety System™
Prevents chair from falling off base



Grow & Save
Product is upgradable in size

Company Policies

Customer Service

For over 30 years, Columbia Medical has been dedicated to providing excellent customer service. Questions? No problem. Contact us at (800) 454-6612 or (562) 282-0244, FAX (310) 305-1718, or you can send an e-mail to info@columbiamedical.com. Our Customer Service Representatives are available Monday through Friday, 7am - 5pm PST/PDT to answer questions regarding any of our products. For added convenience, our website is available 24 hours a day, 7 days a week at www.columbiamedical.com.

Manufacturer's Warranty

Every Columbia Medical product comes with our one year manufacturer's warranty. This warranty ensures that when you receive our products, they will arrive free of defects in material and workmanship. We will repair or replace any part or entire product at our discretion for a period of up to one year from the date of delivery to the original purchaser. This warranty does not include cost of inconvenience, shipping, property damage, misuse, abuse, accident, or similar incidents. The warranty will not apply if the product has been modified, tampered with or repaired by unauthorized individuals. Custom orders are non-refundable and non-returnable.

Shipping

In-stock products will ship within three business days. Shipping charges are included on the invoice when applicable. Products are shipped by FedEx ground. Orders under \$100, for International and Transfer Chair are charged a shipping fee.

Proof of Delivery

Proof of Delivery confirmation is available by contacting Customer Service.

Return Policy

Columbia Medical will provide a refund for all unused products returned within 30 days from the date of purchase, with an approved Return Authorization. Returned items are subject to a reprocessing fee. Products must be new, unused condition, not tampered with, and in original packaging. To obtain a Return Authorization, please call 800-454-6612 or 562-282-0244, Monday through Friday, 7am - 5pm PST/PDT.

Custom orders are non-refundable and non-returnable.

Product information and photography included in this catalog were as accurate as possible at the time of publication. Subsequent refinements may be evident in the actual product. Specifications are subject to change without notice.

For most current specifications, please visit us on the web at www.columbiamedical.com.

COLUMBIA MEDICAL CAR SEATS ARE SPECIFICALLY DESIGNED TO PROVIDE BETTER POSITIONING FOR CHILDREN AND ADULTS WITH SPECIAL NEEDS

What do IPS™ and APS™ mean?

IPS (Integrated Positioning System™)

The Therapedic™ Integrated Positioning System (IPS) car seats have high profile sides to provide “**built-in**” **lateral positioning**. The IPS models (#2000 and #2500) have a more “closed” seat to back angle and are more upright than the APS models. The IPS models **ALWAYS** require tethering to a vehicle tether anchor for installation.

Additional lateral positioning pads may be added to the IPS for better fit and can be removed as the child grows. Height adjustable head pads or headwings come standard with the car seat. Other options include an abductor, an extensor thrust wedge and memory foam pads.

APS (Adjustable Positioning System™)

The Spirit™ Adjustable Positioning System (APS) car seats have **low profile sides** to make transfers easier. The open seating area allows for **wider shoulder and seated leg widths** than the IPS models. Height adjustable headwings provide lateral head positioning. For users who need more comprehensive positioning, the Spirit APS has **optional swing-away hip and trunk supports, and memory foam pads** similar to a wheelchair.

If the Spirit APS is installed with a lap/shoulder belt and the shoulder belt falls within the **labeled safety zone**, additional tethering to a vehicle tether anchor is not required.

Why is a special needs car seat so important?

Growth Capability

Children with special needs often require car seats for a longer period of time than typically developing children. Special needs car seats are designed for and tested to higher weight and height capacities than conventional car seats. **Our car seats can accommodate heights up to 66 inches and weights from 25 to 130 lbs.**

While most conventional car seats have a typical product life span of six years or less, **our car seats have a product life span of TEN years under normal use conditions.**

Specification	2000 IPS	2500 IPS	2400 APS
Max. height	60 inches	66 inches	66 inches
Weight range	20 - 102 lbs.	40 - 130 lbs.	25 - 130 lbs.

Higher Harness Height and Weight Capacities

Most conventional car seats with higher height and weight capacities rely upon the use of the vehicle lap/shoulder belt for occupant protection. However, crash tests by transportation safety experts have proven that **a five-point harness provides better occupant protection** than a vehicle lap/shoulder belt.

In addition, five-point harness systems help to maintain better seating posture, especially in special needs car seats when higher harness heights (measurement from seat to top of user’s shoulder) and higher weight limits must be accommodated. **The Columbia Medical car seats utilize a five-point harness** that allows for harness slot heights up to 23 inches and can support up to 130 lbs. in weight.

Recline* for Improved Upper Body Positioning

Some individuals with disabilities need to **recline to maintain good upper body posture**. Aside from infant car seats, most conventional car seats are too upright for these individuals, and thus do not provide enough support. Special needs car seats are ideal when recline is needed. **Columbia Medical car seats offer up to 45° of recline**, depending on the model and vehicle seat recline.

Specification	2000 IPS	2500 IPS	2400 APS
Seat-to-back angle	88°	85°	98°
Max. recline**	15°	15°	30°

*Recline means “Tilt” in transportation terminology

**When measured on a flat surface. Most vehicle seats add to recline angle.

Accessories for Head, Trunk, Hip and Leg Positioning

Most special needs car seats have five-point harnesses that assist with postural control to promote proper breathing and improve seating stability. **Our car seats offer built-in and/or additional accessories for improved overall body alignment** and head, trunk, hip and leg positioning.

Specification	2000 IPS	2500 IPS	2400 APS
Built-in lateral positioning	Yes	Yes	No
Adjustable Swing-Away Hip/Trunk Supports	No	No	Optional
Height Adjustable Head Supports	Yes	Yes	Yes
Positioning Pad Inserts	Yes	Yes	No
Swing-Away Abductor	Optional	Optional	Optional
Memory Foam Pads	Optional	Optional	Optional

Spirit™ APS™ (Adjustable Positioning System™) Car Seat

Patent #D636,190S

CAR SEAT SAFETY MEETS COMPLEX POSITIONING

The Spirit APS is the only car seat to offer optional Swing-Away Trunk and Hip Supports usually found only in complex seating/wheelchair systems.



Standard configuration

Features

- Accommodates users weighing 25 - 130 lbs. and up to 5'6" tall
- Vehicle shoulder belt may be used to tether in many installations
- Standard 5-point harness, unlike many other seats that utilize the automobile's lap and shoulder restraint
- Meets flammability code FMVSS 302 for automotive upholstery
- Crash tested¹

Benefits

- Can be outfitted with optional adjustable and swing-away trunk and hip supports for better positioning
- Low profile sides for easier transfer of users
- Open shell design provides more room for varying shoulder and seated leg widths
- Multiple accessories allow restraint systems to grow with the child
- Fits in most cars, vans and school buses

Includes

- Adjustable Soft-Padded Headwings
- Durable, machine-washable padded cover
- Recline bar (-30° from vertical when measured on a flat surface)
- 5-Point Safety Harness
- Latex and silicone free


Available Colors



MODEL NUMBER	2400
Product Weight (car seat shell)	27 lbs. basic; 33 lbs. with swing-away trunk and hip supports
Min. Weight Capacity	25 lbs.
Max. Weight Capacity	130 lbs.
Max. Height Capacity	66
Shoulder Harness Heights	11, 13½, 16, 18½, 20½, 23
(a) Back Height	28
(b) Seat Depth (standard)	12
Seat Depth with Extension	16
(c) Inside Seat Width at Front	17 (19 with extension)
(d) Inside Seat Width at Back	15
(e) Shoulder Width	18
Seat to Back Angle	98°

Note: Dimensions are in inches unless otherwise noted

¹MODEL #2400 meets or exceeds crash testing requirements in the US (FMVSS #213)



Accessories for Spirit™ APS™ (Adjustable Positioning System™) Car Seat



Swing-Away Abductor

Improves hip alignment and stabilizes seating posture. Available in blue or walnut colored upholstery.

MODEL #2415 - For use with #2400

MODEL #2419 - For use with 4" seat depth extender #2425



Seat Depth Extender

Increases seat depth by 4".

MODEL #2425¹



Incontinence Cover²

Leak-resistant cover helps shield seating from fluid and soiling.

MODEL #2405 Car Seat Cover²

MODEL #2406 Hip Support Pad Covers (pair)

MODEL #2409 Swing-Away Abductor Cover



Viscoelastic Foam Padding

Set of seat and back foam inserts with an added layer of memory foam to cradle the body prominences for pressure reduction and comfort.

MODEL #2410



BuckleGuard

Prevents child from releasing the five-point harness buckle. For use with all Columbia Medical car seats.

MODEL #2075



Extensor Thrust Wedge

Changes seat to back angle. Increases front seat height by 3" and rear seat height by 1".

MODEL #2455



Tether Strap (Latch Compatible)

Required accessory. Secures the top of the car seat to the vehicle seat.

MODEL #2400P-01



Swing-Away Hip Supports (sold as single units)

Provide adjustable hip positioning and stability. Can completely swing-away to allow for easy transfer of passenger.

MODEL NUMBER	2473-H	2471-H
Hip Width Ranges	10" - 14"	13" - 16"
Brackets	2" x 3", 2" x 3"	2" x 2", 2" x 2"

Hip Support Pads (sold separately as single units)

MODEL NUMBER	2481-S	2481-M	2481-L
Dimensions	3" x 6"	4" x 8"	4" x 12"



Swing-Away Trunk Supports (sold as single units)

Provide adjustable trunk positioning and stability. Can completely swing-away to allow for easy transfer of passenger.

MODEL NUMBER	2474-T	2473-T	2472-T	2471-T
Trunk Width Ranges	6" - 12"	10" - 14"	12" - 15"	14" - 16"
Brackets	2" x 4", 2" x 4"	2" x 3", 2" x 3"	2" x 2", 2" x 3"	2" x 2", 2" x 2"

Trunk Support Pads (sold separately as single units)

MODEL NUMBER	2482-S	2482-M	2482-L
Dimensions	3" x 3"	3" x 5"	4" x 6"

¹When ordering this item for an existing Spirit APS Car Seat, a new seat cover upholstery (MODEL #2402) must be purchased; please specify color Blue or Walnut.

²Covers for swing-away abductor and swing-away hip pads are sold separately.

Therapedic™ IPS™ (Integrated Positioning System™) Car Seat

THE CHOICE OF PARENTS AND THERAPISTS WHEN "BUILT-IN"
LATERAL POSITIONING IS REQUIRED



*Shown with head pads



*Shown with headwings

Includes

- Independently Adjustable Height Head Pads or Adjustable Soft-Padded Headwings
- Shoulder Harness Pads
- Crotch Strap Pad
- Lateral Positioning Pads
- 5-Point Safety Harness
- Dual Hook Tether Strap
- Easy to Follow Installation Instructions

Available Colors



Blue



Rose

Features

- Two sizes available, accommodating users weighing 20 - 130 lbs. and up to 5'6" tall
- Durable, machine washable, comfortable padded cover
- Standard 5-point harness restraint for better protection than vehicle lap/shoulder belt
- Meets flammability code FMVSS 302 for automotive upholstery
- Latex and silicone free
- Crash tested*
- Positioning options available

Benefits

- Fits in most cars, vans and school buses
- Model #2000 is certified for use in aircraft
- Lightweight, weighing much less than most comparable car seats
- Optional seat depth extender allows car seat to grow with the child



MODEL NUMBER	2000	2500
Product Weight	11 lbs.	18 lbs.
Weight Capacity	20 - 102 lbs.	40 - 130 lbs.
Height Capacity	Up to 60	54 - 66
Shoulder Harness Heights	11½, 14½, 17½, 20½	13½, 17, 20½, 24
(a) Seat to Top Height	29	33
(b) Seat Depth	10	12
Seat Depth with Extender	12½ or 15	16
Inside Seat Width at Front	15½	19½
Inside Seat Width at Back	13	17½
(c) Outside Width	19	23
Outside Depth	11	12
Seat to Back Angle	88°	85°

Note: Dimensions are in inches unless otherwise noted

*Both MODEL #2000 and #2500 meet or exceed crash testing requirements in the US (FMVSS #213)

Accessories for Therapedic™ IPS™ (Integrated Positioning System™) Car Seat



Seat Depth Extender

Increases seat depth.

For use with **MODEL #2000**

MODEL #2025¹ - Increase seat depth to 12½"

MODEL #2026¹ - Increase seat depth to 15"

For use with **MODEL #2500**

MODEL #2525¹ - Increase seat depth to 16"



Swing-Away Abductor

Separates the knees to provide hip abduction.

For use with **MODEL #2000**

MODEL #2015¹

MODEL #2018¹ - For use with MODEL #2025

MODEL #2019¹ - For use with MODEL #2026

For use with **MODEL #2500**

MODEL #2515¹

MODEL #2519¹ - For use with MODEL #2525



Extensor Thrust Wedge

Changes seat to back angle.

Increases front seat height by 2" or 3".

For use with **MODEL #2000**

MODEL #2050 - 2" front height to 1" rear

MODEL #2055 - 3" front height to 1" rear

For use with **MODEL #2500**

MODEL #2550 - 2" front height to 1" rear

MODEL #2555 - 3" front height to 1" rear



Block Style Trunk Support Pads

Additional pads for customized positioning. For use with all Therapedic IPS Car Seats. Available in four packs. Block style, two pads 2½" thick and two pads 1½" thick.

MODEL #2008 (Blue or Rose)



Independently Height Adjustable Head Pads

Individual head cushions (set of two)

MODEL #2001 (Blue or Rose)



Height Adjustable Headwings

Adjustable soft-padded headwings

MODEL #2501 (Blue or Rose)



Lateral Positioning Pads

Four additional pads for customizing lateral positioning. Inserted under the cover.

MODEL #2010 - For use with MODEL#2000

MODEL #2510 - For use with MODEL#2500



BuckleGuard

Prevents child from releasing the five-point harness buckle. For use with all Columbia Medical car seats.

MODEL #2075



Viscoelastic Foam Padding

Set of seat and back foam supports with an added layer of memory foam to cradle the body prominences for pressure reduction and comfort.

MODEL #2011 - For use with MODEL #2000

MODEL #2511 - For use with MODEL #2500

¹When ordering this item for an existing Therapedic IPS Car Seat, a new seat cover upholstery (MODEL #2002 or MODEL #2502) must be purchased; please specify color Blue or Rose.

Therapedic™ IPS™ Transport Bases

LIGHTWEIGHT, FOLDING FRAMES DESIGNED SPECIFICALLY FOR USE WITH THERAPEDIC™ INTEGRATED POSITIONING SYSTEM™ (IPS™) CAR SEATS. OUR TRANSPORT BASES AND UMBRELLA TRANSPORT BASES ARE **IDEAL ALTERNATIVES TO TRADITIONAL POSITIONING WHEELCHAIRS** AND PERFECT FOR USE WHILE SHOPPING, VISITING DOCTORS, GOING OUT, OR TRAVELING.



Therapedic IPS Transport Base

- Designed for comfort – Five angles of seat tilt: 20°, 25°, 30°, 35°, 40°
- Height adjustable footrest
- Built to go – Large wheels for a smooth ride on any surface
- Convenient – Folds down to carry, travel, or store
- Easy to use – Seat locks into the base with ease
- Durable – Lightweight, yet strong frame

Accessories: Tote Bag

Attaches to the Mobility Base

MODEL #2175 (Blue or Rose)



MODEL NUMBER	2100	2520
FOR USE WITH MODEL	2000	2500
Product Weight	25 lbs.	33 lbs.
Weight Capacity	20-102 lbs.	40-130 lbs.
Open	22½ W x 49 L x 36¼ H	25½ W x 54½ L x 37¾ H
Folded	22¼ W x 41 L x 11 H	25½ W x 47½ L x 13 H
Leg Drop Range	5 -16	6 - 13
Front Wheel Diameter	7	7
Rear Wheel Diameter	10	10

Note: Dimensions are in inches unless otherwise noted

Therapedic IPS Umbrella Transport Base

- Versatile – Use alone or with the Therapedic IPS Car Seats
- Built to go – Large wheels provide a smooth ride on any surface
- Convenient – Folds compactly to carry, travel, or store
- Easy to use – Seat easily locks into the base
- Durable – Strong, yet lightweight frame

Accessories

- Heel Loops **MODEL #2301**
- Shoulder Straps **MODEL #2303**
- Canopy **MODEL #2304**

MODEL NUMBER	2300
FOR USE WITH MODEL	2000
Product Weight	22 lbs.
Maximum Weight Capacity	170 lbs.
Open	25 W x 36½ L x 47½ H
Folded	14 W x 47½ L x 15 H
Leg Drop Range	14 - 18
Seat Width	16
Seat Depth Range	14 - 17
Seat Length	14
Seat Back Height	25
Wheel Diameter	7

Note: Dimensions are in inches unless otherwise noted



Umbrella Transport Base MODEL #2300 shown with MODEL #2000



Commitment to Bathroom Safety

COLUMBIA MEDICAL OFFERS A FULL LINE OF BATHING/ SHOWERING/COMMODE POSITIONING SUPPORT OPTIONS TO MAXIMIZE SAFETY AND PROVIDE EASE OF USE

Columbia Medical's bathing aids have been designed to make home care safer and easier. Our products fit a wide range of home environments and meet the varying needs of users and caregivers. In many instances, our bathing systems are a lower cost alternative to home modifications. Our PVC line of Blue Hub products offers unsurpassed durability and versatility. Columbia Medical offers a full line of equipment to assist in the bathroom.

Checklist for selecting bathing equipment

The User and Caregiver

- Age, weight and size of the user
- Functional independence of the user
- Caregivers and their ability to assist
- Other individuals who will also be using the same bathroom
- Postural supports needed to accomplish safe bathing
- Recline or tilt-in-space positioning needs for postural control
- Dystonia or other reflexive movements
- Access for hair washing and other personal hygiene care

Durability and Adjustability of Equipment for Growth

- Expected length of use of the equipment
- Accommodation for growth and change of the user

The Bathroom Environment

- Whether a bathtub or a roll-in shower will be used
- If the same piece of equipment will be used for both bathing and/or toileting
- The size and configuration of the space
- Space around the commode/tub/shower for equipment
- Doorway widths to allow access of equipment
- Thresholds to cross
- Carpet or flooring in bathroom and adjacent rooms
- Storage space for equipment

Other Assistive Equipment In Use

- Will the person transfer from a bed, chair or wheelchair
- Lift equipment used for transfers
- Wheelchair turning radius

Bath Chair/Support Comparison Chart

Specification	Ultima™ Bath Chairs	Supreme™ Bath Chairs	Deluxe™ Bath Chairs	Wrap Around Bath Supports
Material Construction	Stainless Steel	PVC	PVC	PVC/ABS
Positioning Features	Back and Calf Recline	Back and Calf Recline	Tilt-in-Space	Upright sitting - Hi-Back™ or Lo-Back™ support
Rolling Base	Optional	Optional	Optional	N/A

Bath Lifts

Specification	Aqualift™
Features	Transfer flap, tool-free installation
Positioning Features	Up to 22° tilt/recline. Can be ordered with or without positioning Hi-Back or Lo-Back support.

Bath Transfer Systems

Specification	Elite™ & EliteVerse™ Transfer Systems/Rolling Chairs	Omni™ & OmniVerse™ Transfer Systems/Rolling Chairs	Versa™ Transfer Systems/Rolling Chairs
Material Construction	Stainless Steel	Stainless Steel	Stainless Steel
Function	Bath/Shower/Transfer	Bath/Shower/Commode/Transfer	Bath/Shower/Commode/Transfer
Positioning Features	Back and Calf Recline	Back and Calf Recline	Upright

Shower/Commode Chairs

Specification	Optima™ TSS™ (Toilet & Shower System)	5015 Foldable	5050 Self-Propelled	5060, 5070, 5075, 5080 Hi-Back	5050 Bariatric
Material Construction	Stainless Steel	PVC	PVC	PVC	PVC
Max. Weight	250 lbs.	250 lbs.	300 lbs.	525 lbs.	700 lbs.
Positioning Features	Back tilt with high or low back support	Upright	Upright	Upright	Upright

Ultima™ Stainless Steel Bath Chair

Patent #D507,904; D508,175; D503,556; 7,328,953B2



Custom sizes available

Features

- Independently adjustable seat, back and leg angles
- Back angles: 0°, 30°, 60°, 90° (from vertical)
- Leg angles: 0°, 30°, 60° (from horizontal)
- Durable surgical grade stainless steel frame
- Bacteria and fungal resistant ComfortMesh™ seat cover
- Latex and silicone free
- Includes two positioning belts

Benefits

- Seat length options available to meet the needs of growing children
- Upgrade to the next size for a fraction of the cost of a new chair
- Full head, trunk and calf support in seated or reclined positions
- Designed to contour to the body
- Folds flat for storage

Available Mesh Colors



Blue



Rose



Accessories (See page 11 for details)

MODEL NUMBER	8520	8620	8720	8920
Product Weight	19 lbs.	20½ lbs.	21½ lbs.	23 lbs.
Weight Capacity	100 lbs.	130 lbs.	160 lbs.	180 lbs.
Height Capacity	< 40	40 - 54	54 - 68	68 - 80
Overall Length Flat	44	45	53	61
Leg Length Flat	39½	39½	41½	44½
(a) Overall Width	17	17	18	20
(b) Inside Seat Width	14	14	15	17
(c) Seat Depth	10	11	13	15
(d) Calf Length	10	10	13	15
(e) Back Height	20	22	27	31
(f) Seat Height Above Floor	4 - 15	4 - 15	4 - 15	4 - 15

Note: Dimensions are in inches unless otherwise noted

Accessories for Ultima™ Stainless Steel Bath Chair

Conversion Kits

Growth kits are available to convert your Ultima bath chair to a larger size. Each conversion kit includes a new cover with two belts, stainless steel frame and legs.

PRODUCT
Ultima Conversion Kit

MODEL #	CONVERTS TO
8629	8620

MODEL #	CONVERTS TO
8729	8720

MODEL #	CONVERTS TO
8929	8920

Lateral Support

Cushion designed for better torso positioning.

MODEL #8830
(Blue or Rose)

Detachable Armrests

MODEL #5150

Anterior Support Tray

Attaches to armrests MODEL #5150.

MODEL #5151

Padded Abductor

Attaches easily to the Ultima chair to improve hip alignment and stabilize seating posture.

MODEL #8818 (Blue Only)



Head Support

Cushion designed for better head positioning.

MODEL #8875 (Blue or Rose)

Positioning Belt

Quick release belt usable in multiple positions. Two standard positioning belts included with each bath chair.

MODEL #8835 (1 ea.) (Blue or Gray)

Detachable Rolling Base

Stainless steel base for use in roll-in showers. Includes four 4" locking casters and straps to secure chair to base.

MODEL #8890

Max. Weight Capacity: 180 lbs.

Width at Base: 27"

Adjustable Length Range: 31" – 44"

Raises Bath Chair by: 13"

Pelvic Positioning Belt (Not Shown)

Improves posture and safety.

MODEL #8834

Ultima Bases

Columbia Medical offers three types of bases for the Ultima bath chairs:

- Detachable Rolling Base (see above)
- Detachable Extension Base (see below)
- Attached Rolling Base (see below)



Attached Rolling Base

Aluminum base is secured to Ultima chair for added stability. Ideal for use in roll-in showers. Includes four 4" locking casters.

MODEL #8895

Max. Weight Capacity: 180 lbs.

Dimension: 22" x 22"

Seat Height Range: 23" – 27½"



Detachable Extension Base

Stainless steel extension legs for use in bathtubs. Includes straps to secure chair to base.

MODEL #8880

Max. Weight Capacity: 180 lbs.

Width at Base: 18"

Adjustable Length Range: 31" – 38"

Raises Bath Chair by: 9½"

Accessories for Attached Rolling Base



Footrests

Provide foot support.

MODEL #5125



Seat Angle Adjustment Brackets

Offer different seat angle positions. Raises front seat height 2½", 3", or 3½".

MODEL #5126



Rear Anti-Tippers

For added safety and stability.

MODEL #5127

Contour™ Supreme™ Reclining Bath Chair

Patent #D606,389; D606,387; D606,388; D592,492; 3,795,218

INDUSTRY LEADING **BLUE HUB** CONSTRUCTION BATH CHAIR
FEATURING INDEPENDENTLY ADJUSTABLE SEAT AND BACK



Features

- Proprietary "Blue" PVC hub provides greater strength and durability for use in clinical applications
- Seat and back sections adjust independently for multiple positions
- Back angles: 0°, 15°, 30°, 45°, 60° (from vertical)
- Seat angles: 0°, 15°, 30° (from horizontal)
- Simple button press for positioning
- Bacteria and fungal resistant ComfortMesh™ seat cover
- Latex and silicone free
- Includes two positioning belts


Benefits

- Full head, trunk and calf support in seated or reclined positions
- Folds flat for storage
- Seat back cover slides down for easy shampooing on **MODEL #8510** and **MODEL #8610** only
- **MODEL #8710** has two adjustable bars for additional support

Accessories (See page 14 for details)

- Lateral Support **MODEL #8830** (Blue or Rose)
- Positioning Belt **MODEL #8835** (Blue or Gray)
- Head Support **MODEL #8875** (Blue or Rose)
- Padded Abductor: (Blue only)

Available Mesh Colors



MODEL NUMBER	8510	8610	8710*
Product Weight	14 lbs.	15 lbs.	21 lbs.
Weight Capacity	100 lbs.	130 lbs.	180 lbs.
Height Capacity	< 40	40 - 54	54 - 68
Overall Length	35	45	53
(a) Seat Depth	13	17	20
(b) Back Height	21	30	34
(c) Inside Seat Width	12	12	17
(d) Outside Width	17	17	24
(e) Base Length	27½	27½	27½
(f) Seat Height Above Floor	9¼	9¼	10¼

Note: Dimensions are in inches unless otherwise noted

*See photo on page 14

Contour™ Deluxe™ Tilt-in-Space Bath Chair

Patent #D606,389; D606,387; D606,388; D592,492; 3,795,218

PVC BATH CHAIR WITH THREE FIXED TILT POSITIONS



Features

- Proprietary "Blue" PVC hub provides greater strength and durability for use in clinical applications
- Seat and back sections adjust together
- Tilt angles: 30°, 50°, 70°
- Simple button press for tilt angle changes
- Bacteria and fungal resistant ComfortMesh™ seat cover
- Latex and silicone free
- Includes two positioning belts

Benefits

- Head and trunk support
- Seat back cover slides down for easy shampooing on **MODEL #8500 and MODEL #8600** only
- **MODEL #8700** has two adjustable bars for additional support

Accessories (See page 14 for details)

- Lateral Support **MODEL #8830** (Blue or Rose)
- Positioning Belt **MODEL #8835** (Blue or Gray)
- Head Support **MODEL #8875** (Blue or Rose)
- Padded Abductor: (Blue only)



Available Mesh Colors



MODEL NUMBER	8500	8600	8700*
Product Weight	7 lbs.	10 lbs.	12 lbs.
Weight Capacity	100 lbs.	130 lbs.	180 lbs.
Height Capacity	< 40	40 – 54	54 - 68
Overall Length	29	40	45
(a) Seat Depth	10	16	19
(b) Back Height	21	27	33
(c) Inside Seat Width	12	12	17
(d) Outside Width	17	17	20
(e) Base Length	28	28	28
(f) Seat Height Above Floor	2 - 4	2 - 4	2 - 4

Note: Dimensions are in inches unless otherwise noted

*Design of MODEL #8700 is slightly different. Contact Customer Service for details.

Accessories for Contour™ Bath Chairs

Patent #D606,389; D606,387; D606,388; D592,492; 3,795,218



MODEL #8710 Bath Chair shown with optional accessories

Head Support

Cushion designed for better head positioning.

MODEL #8875 (Blue or Rose)

Positioning Belt

Quick release strap usable in multiple positions. Two standard positioning belts included with each bath chair.

MODEL #8835 (1 ea.) (Blue or Gray)

Lateral Support

Cushion designed for better torso positioning.

MODEL #8830 (Blue or Rose)

Padded Abductor

Attaches easily to bath chair to improve hip alignment and stabilize seating posture. (Blue Only)

MODEL #8815 - Fits MODEL #8500, 8510

MODEL #8816 - Fits MODEL #8600, 8610

MODEL #8817 - Fits MODEL #8700, 8710

PVC Rolling Shower Base

Frame fits Contour Supreme Reclining and Deluxe Tilt-in-Space models. Ideal for use in roll-in showers. Includes 4" EZ Clean locking casters and two height positions.

MODEL #8870

Product Weight: 16 lbs.

Width: 17"

Length: 28"

Raises Bath Chair Height: 16" or 24"

EZ Clean™ Casters

Four locking rustproof wheels with rubber surface for added safety. Removable axle for easy cleaning. Fits Supreme Reclining Bath Chair, Deluxe Tilt-in-Space Bath Chair, and Extension Legs.

Height: 4"

MODEL #5505

Extension Legs

PVC legs raise bath chair height. Fits all Contour bath chairs.

MODEL #8850

Product Weight: 8 lbs.

Width: 17"

Length: 28"

Raises Bath Chair Height: 11"


MODEL #8860

Product Weight: 11 lbs.

Width: 17"

Length: 28"

Raises Bath Chair Height: 15"



Wheelbase Extension Kits

Fits all PVC Extension Legs and Rolling Base.

MODEL #8872 Side Wheelbase Extension Kit

MODEL #8873 Front and Rear Wheelbase Extension Kit

Wrap-Around Bath Supports



Lo-Back shown without Padding



Hi-Back shown with Padding



PROVIDE EXTRA SUPPORT AND IMPROVED SEATING POSTURE WHILE SITTING

Features

- Hi-Back™ and Lo-Back™ options available (Hi-Back model includes four-point harness for added trunk support)
- Orthopedic contoured back with support belts
- Padded back models available for pressure reduction
- PVC frame
- Bacteria and fungal resistant ComfortMesh™ seat
- Latex and silicone free

Benefits

- Orthopedic back height adjusts vertically
- May be used as bath support or floor sitter
- Wide-stance stable base fits into any tub
- Easily adjustable back heights
- Compact design for easy storage or travel
- Rustproof frame and tear resistant mesh fabric for extra durability
- Mounting Base allows you to upgrade or replace the support back to the next size for a fraction of the cost of a new support

Available Mesh Colors*



Blue




Rose

*Please note padded backs are Blue.

WRAP-AROUND BATH SUPPORTS (PADDED BACK IS BLUE)	LO-BACK SUPPORTS WITH PADDING			HI-BACK SUPPORTS WITH PADDING		
MODEL NUMBER	7210	7310	7410	7610	7710	7910
Product Weight	5½ lbs.	6½ lbs.	7¼ lbs.	9½ lbs.	9¾ lbs.	12 lbs.
Overall Depth	29	29	29	29	29	29
(a) Seat Depth	14½	14½	14½	14½	14½	14½
(b) Seat Width Overall	18	18	18	18	18	18
(c) Seat Width Inside	15	15	15	15	15	15
(d) Back Width at Top Inside	7	11	13	7	8	11
(e) Back Width at Bottom Inside	7	11	13	12	12	14
(f) Back Height Range	9 - 14	10 - 14	11 - 16	17 - 22	19 - 24	24 - 28
OTHER OPTIONS						
MODEL NUMBER (w/o Padding) ¹	7200	7300	7400	7600	7700	7900
MODEL NUMBER (Pads only)	7205	7305	7405	7605	7705	7905

Note: Dimensions are in inches unless otherwise noted

¹ Inside back width increases by 1" for Wrap-Around supports without padding



Accessories for Wrap-Around Bath Supports

Padded Abductor

Attaches to improve hip alignment and stabilize seating position.

MODEL #7815 (Blue)

Pelvic Support Belt

Adds extra security when attached to base. Maintains correct pelvic angle.

MODEL #7830 (Blue or Gray)

AquaLift™ Bath Lifter System

FULLY PORTABLE, EASY TO STORE,
LIGHTWEIGHT, RECLINING BACK BATH
LIFT. IDEAL FOR HOME USE. DESIGNED
FOR PEOPLE WHO HAVE DIFFICULTY
GETTING IN AND OUT OF THE TUB.

Features

- Low voltage 14.4V "Smart" battery lifts up to 300 lbs. This safety mechanism prevents the chair from lowering when there is not enough power to fully raise after use
- Suction cups keep lifter from sliding
- Back gradually reclines from 8° off of vertical at highest position to 22° at bottom position (seat remains flat)
- Side flap for easy lateral transfer
- Easy to use soft-touch control
- Battery charger
- Latex and silicone free

Benefits

- Easy to operate
- No tools required for installation
- Fits most bathtubs
- Provides versatile positioning
- Easy to transport
- Easily disassembles for storage or travel

Accessories

- Lateral Support **MODEL #8100P-13**
(Blue or Rose, see pg. 14)
- Positioning Belt **MODEL #8835** (Blue or Gray)
- Head Support **MODEL #8100P-12** (Blue or Rose)
- Lo and Hi-Back Wrap-Around Supports, with or without padding (see charts below for specs)



Lo-Back™ Support with Padding



Hi-Back™ Support with Padding




Head Support Model #8100P-12

MODEL NUMBER	8100
Product Weight	21.25 lbs.
Weight Capacity	300 lbs.
(a) Seat Width	14½
(b) Seat Depth	15
(c) Seat Height Range	2 - 16
(d) Back Height	30
Folded Dimensions	15 W x 28 L x 7.8 H
Battery Specification	14.4V DC nominal, 7A max
Lifts Per Full Charge	- Approximately 6 lifts with an average of 165 lbs. - Approximately 3 lifts when lifting the max. weight of 300 lbs.
Reclines	- Reclines 82° when lift is at its highest position and 57° when at its lowest position (seat does not tilt)

WRAP-AROUND SUPPORT DIMENSIONS				
MODEL NUMBER	8113	8114	8116	8117
Product Weight	5 lbs.	6 lbs.	7 lbs.	8 lbs.
Support Height	8	10	16	18
Support Inside Back Width at Top	11	13	7	8
Support Inside Back Width at Bottom	11	13	12	12
Complementary padded back (Model #)	5405	5405	5605	5805

Note: Dimensions are in inches unless otherwise noted



AquaLift Only

POSITIONING COMMODES WITH ORTHOPEDIC CONTOURED BACKS CAN BE USED FREESTANDING OR CAN FIT OVER ANY TOILET



Hi-Back™ Support without Padding



Lo-Back™ Support without Padding



Lo-Back Support with Padding

Features

- Orthopedic contoured back in Hi-Back and Lo-Back options (Hi-Back models include harness for added support)
- Padded back models available for pressure reduction
- Includes a commode pail, footrest, padded reducer ring 7" W x 8" L hole (5" W x 6" L hole can be requested) with removable abductor/splashguard and pelvic belt
- Adjustable seat depth, back height and footrest
- Lightweight steel frame
- Latex and silicone free



Benefits

- Versatile design for use as a bedside commode or over-the-toilet commode
- Easily adjustable back heights
- Easy push-button pail release
- U-Shape Mounting Base allows you to upgrade or replace the support back to the next size for a fraction of the cost
- Not for shower use

Accessory

- Casters **MODEL #5505**
(set of four 4" EZ Clean™ locking casters)



	LO-BACK SUPPORTS WITH PADDING			HI-BACK SUPPORTS WITH PADDING			
MODEL NUMBER	5210 ² (RSD)	5310	5410	5610 ² (RSD)	5710 ² (RSD)	5810	5910
Product Weight	27 lbs.	28 lbs.	29 lbs.	30 lbs.	30½ lbs.	30½ lbs.	33 lbs.
Weight Capacity	175 lbs.	175 lbs.	175 lbs.	175 lbs.	175 lbs.	175 lbs.	175 lbs.
(a) Seat Depth Range	14 - 16	13 - 19	13 - 19	14 - 16	14 - 16	13 - 19	13 - 19
(b) Inside Back Width at Top	7	11	13	7	8	8	11
(c) Inside Back Width at Bottom	7	11	13	12	12	12	14
(d) Back Height Positioning Range	8 - 17	8 - 18	11 - 19	17 - 25	18 - 26	18 - 26	23 - 30
(e) Fixed Seat Heights	17 7/8 19 3/8 20 7/8 22 3/8	17 7/8 19 3/8 20 7/8 22 3/8	17 7/8 19 3/8 20 7/8 22 3/8	17 7/8 19 3/8 20 7/8 22 3/8	17 7/8 19 3/8 20 7/8 22 3/8	17 7/8 19 3/8 20 7/8 22 3/8	17 7/8 19 3/8 20 7/8 22 3/8
(f) Back Height at Center	6 3/4	7	10 1/2	16	18	18	22 1/2
	LO-BACK SUPPORTS			HI-BACK SUPPORTS			
MODEL NUMBER (Back support w/o Padding) ¹	5200 ² (RSD)	5300	5400	5600 ² (RSD)	5700 ² (RSD)	5800	5900
MODEL NUMBER (Pads only)	5205	5305	5405	5605	5705	5805	5905

¹Inside back width increases by 1" for Positioning Commodes without pads

²Reduced seat depth (RSD): sits 4" forward to lessen seat depth when attached to U-Shape Mounting Base

Note: Dimensions are in inches unless otherwise noted



Hi-Back Support with Padding

Optima™ TSS™ (Toilet & Shower System)

MULTI-USE SYSTEM FOR COMPLEX POSITIONING NEEDS

Features


- Durable, adjustable height frame
- Unique Hi-Back™ and Lo-Back™ wrap around supports, available in three sizes
- Height, angle, and depth adjustable back support
- Soft-padded cushion for pressure reduction, available in three sizes
- 5° seat angle for positioning
- Built-in seat depth for growth
- Width adjustable, swing-away armrests, available in two heights
- Height adjustable, swing-away footrests
- EZ Clean™ Casters with removable axle for easy cleaning
- Latex and silicone free

Benefits

- Multi-use system; functions as a stand-alone commode, over-the-toilet system, or roll-in shower chair
- System can grow in seat depth and accepts multiple back sizes
- Configuration and accessory options can accommodate complex positioning needs
- Easily adjusts to varying toilet configurations
- Small footprint and tight turning radius fits in small spaces
- Armrests provide support and stability
- Armrests and footrests swing away for easier transfer
- Anterior support tray available
- Easily disassemble into a compact package for transport or storage



Shown with optional accessories



MODEL NUMBER	5500 (ROLLING BASE FRAME ONLY)	
Product Weight	Frame only: 23.5 lbs.	Frame with Medium Seat & Back Supports: 35 lbs.
Weight Capacity	250 lbs.	
(a) Overall Frame Depth (at maximum seat depth)	Without Footrest: 20½	With Footrest: 26¾
(b) Overall Frame Width	Without Casters: 18½	With Casters: 20
(c) Inside Frame Width	16	
(d) Clearance to Floor (for over toilet)	13¾, 14¾, 15¾, 16¾, 17¾, 18¾	
(e) Armrest Width Range	Inside Width: 16½, 18½, 20½, 22½	Outside Width For Clearance: 20, 22, 24, 26
(f) Armrest Height Above Seat	Low Armrest: 6½ (S/M seat), 6 (L seat)	High Armrest: 8½ (S/M seat), 8 (L seat)
(g) Footrest to Floor	7-20 (Footrest stems can be flipped up-side-down for maximum ranges)	

COMMODE CUSHION	5562	5563	5564
Seat Depth	12 - 14	14 - 16	17 - 19
Seat Width	13	14	16
Seat Height Above Floor (in 1" increments)	20½ - 24½	20½ - 24½	21¼ - 25¼
Seat to Tray (w/low armrest, w/high armrest)	6½, 8½	6½, 8½	6¼, 8¼
Seat to Footrest	5 - 18	5 - 18	4½ - 17½

PADDED BACK SUPPORT	LOW-BACK OPTIONS			HIGH-BACK OPTIONS		
	5542	5543	5544	5546	5548	5549
(h) Back Height Range from Seat	6¼ - 11¼	7 - 11½	10¼ - 14¾	16 - 20½	18 - 22½	22½ - 27
(i) Back Angle Range	-30° To +45°	-30° To +45°	-30° To +45°	-30° To +45°	-30° To +45°	-30° To +45°
Inside Width at Top, Bottom	7, 7	11, 11	13, 13	7, 12	8, 12	11, 14

Note: Dimensions are in inches unless otherwise noted


Accessories for Optima™ TSS™ (Toilet & Shower System)

Swing-Away Armrests

Available in two heights.

MODEL #5550-L Low, 8"

MODEL #5550-H High, 10½"



Cushion (Shown With Abductor/Splashguard)

Coming Soon!

Commode Pail

MODEL #5500P-02

Flip-up, Swing-Away Footrests

Range 5" - 18" seat to footrest.
Swings out of the way for transfers.

MODEL #5525

Pelvic Positioning Belt (Not Shown)

Improves posture and safety.

MODEL #5203

Anterior Support Tray

Fits 16½" and 18½" armrest widths. For use with #5550 Armrests only.

MODEL #5151

Rear Anti-Tippers

For added safety and stability.

MODEL #5527

EZ Clean™ Casters (set of 4)

Locking casters with removable axle for easy cleaning.

MODEL #5505

Adjustment Options



Base Options for Optima TSS



Transfer System (OptiVerse)



21" Wide Base
(17" Inside Width)

Compact for Easy Transport and Storage



Toilet Support System



Hi-Back Toilet Support without Padding

TOILET SUPPORTS FEATURE AN ORTHOPEDIC CONTOURED BACK WITH SUPPORT BELTS

Features

- Orthopedic contoured backs in Hi-Back™ and Lo-Back™ options (Hi-Back models include harness for added support)
- Padded back models for pressure reduction
- Easy to clean fabrics and materials
- Latex and silicone free

Benefits

- Easy to attach and remove from toilet
- Padded backs provide pressure reduction through seamless contour design
- Increased privacy, security, independence and dignity
- Easily adjustable back height
- U-Shape Mounting Base allows you to upgrade or replace the support back to the next size for a fraction of the cost
- Additional U-Shape Mounting Base MODEL #4001 can be purchased so support can be used on other toilets

Accessories

(See page 19 for details)



Lo-Back Toilet Support with Padding



U-Shape Mounting Base



LO-BACK TOILET SUPPORTS WITH PADDING				HI-BACK TOILET SUPPORTS WITH PADDING			
MODEL NUMBER	4210 ^{2,3} (RSD)	4310	4410	4610 ^{2,3} (RSD)	4710 ^{2,3} (RSD)	4810	4910
Product Weight	5½ lbs.	6½ lbs.	7¾ lbs.	9½ lbs.	9 ¾ lbs.	9¾ lbs.	12 lbs.
Inside Back Width at Top	7	11	13	7	8	8	11
Inside Back Width at Bottom	7	11	13	12	12	12	14
Back Height Range	7 - 14	10 - 14	10 - 17	17 - 22	19 - 24	19 - 24	24 - 28
U-Shape Mounting Base Dimension	18½ W x 3 D x 10 H			18½ W x 3 D x 10 H			
OTHER OPTIONS							
MODEL NUMBER (w/o Padding) ¹	4200 ^{2,3} (RSD)	4300	4400	4600 ^{2,3} (RSD)	4700 ^{2,3} (RSD)	4800	4900
MODEL NUMBER (Pads only)	4205	4305	4405	4605	4705	4805	4905

¹Inside back width increase by 1" for Positioning commodes without pads

²Reduced seat depth (RSD): sits 4" forward to lessen seat depth when attached to U-Shape Mounting Base

³Includes #4950 plastic reducer ring

Note: Dimensions are in inches unless otherwise noted

Accessories for Toilet Support System

Armrest and Footrest Options

MODEL #4921 - Armrest only, does not include U-Shape Mounting Base

MODEL #4001 - U-Shape Mounting Base, does not include armrests
18½" W x 3" D x 10" H

MODEL #4922 - Footrest only, does not include U-Shape Mounting Base

MODEL #4925 - Armrest and Footrest with brackets, does not include U-Shape Mounting Base

MODEL #4930 - Armrest with U-Shape Mounting Base

Armrest Width Outside: 19"

Armrest Width Inside: 15½"

Height Adjustment: 20" - 31"

Length Adjustment: 19" - 22"

Footrest Dimensions: 17" W x 7" H x 7" D

Seat to Footrest: 17" to floor



Hi-Back™ Support Option



Plastic Reducer Ring

Fits on any toilet seat without tools or modification. Seamless and easy to clean. Built-in splashguard.

MODEL #4950

Hole: 6" W x 7" L

Outside:

10" W x 13" L

Abductor:

3½" W x 2" H



Padded Reducer Ring

For those who require a reduced opening. Includes replaceable abductor splashguard. Seamless and easy to clean. Fits on any toilet seat without tools or modification, 1½" thick.

MODEL #4965

Hole: 5" W x 6" L

Outside:

14" W x 15½" L

Abductor:

3½" W x 3½" H

MODEL #4960

Hole: 7" W x 8" L

Outside:

14" W x 15½" L

Abductor:

3½" W x 4½" H



Soft-Flex™ SplashGuard

For adults and larger children. Slides off easily after use. Fits open or closed front toilet seat. Cleans easily.

MODEL #4975

5" W x 5½" L

MODEL #4976

Replacement Straps
2 ea. for splashguard

Elite™ Bath/Shower Transfer System

US Patent #D508,175; D503,556; D611,715; D507,904; 7,328,953



MAXIMIZE SAFETY AND CARE WITH OUR BATH/SHOWER/ COMMODE TRANSFER SYSTEMS

Features

- Multi-functional transfer systems for use in home or facility without any bathroom modifications
- Folding legs for easy storage
- Custom tub lengths available
- Four recline positions. Back angle: 0°, 30°, 60°, 90°
- Chair locks into position once transfer is complete
- Soft, bacterial and fungal resistant ComfortMesh™ seat cover
- Includes two trunk positioning belts
- Latex and silicone free

Benefits

- Provide safer bathroom transfers that fit the individual's bathing and postural needs
- Eliminate lifting of individual within the bathroom
- Adaptable to multiple bathroom configurations
- Built with sturdy aluminum frames to maximize safety and life span
- Bacterial and fungal resistant ComfortMesh™ seat cover
- Fully adjustable chair that can be retrofitted to the next size for a fraction of the cost of a new chair



Accessories

(See page 24 for details)

Omni™ Bath/Shower/Commode Transfer System

US Patent #D508,175; D503,556; D611,715; D507,904; 7,328,953



SLIDER SAFE™ Auto-Stop Safety System™



As a stand-alone commode chair



Over-the-toilet

EliteVerse™ and OmniVerse™ Transfer Systems

US Patent #D508,175; D503,556; D611,715; D507,904; 7,328,953

INTRODUCING **NEW VERSATILE TRANSFER SYSTEMS IDEAL FOR SMALL BATHROOMS.** THE SYSTEM INCORPORATES **ELITE OR OMNI SEATING** ON A ROLLING BASE AND INTERFACES WITH **A COMPACT TRANSFER FRAME.**



EliteVerse
Bath/
Shower
System

Features

- Compact, 18" x 18", tub base transfer frame fits inside bathtub
- Adjustable or various bridge length options
- Locking system for added safety



OmniVerse
Bath/
Shower/
Commode
System

Benefits

- Reduce the need for multiple pieces of equipment
- Accommodate small bathroom spaces
- Designed for independent use, or with assistance
- Shower door or curtain can be drawn when bridge is removed

Versa™ (Fixed Back) Bath/Shower/Commode Transfer System

Features

- Compact 18" wide tub base transfer frame
- Locking system for added safety
- Height adjustments to fit over most toilets and tubs
- Removable 13" bridge¹
- Lightweight aluminum construction
- Flip-up armrests and swing-away footrests
- Latex and silicone free
- Includes commode pail, locking casters, and footrests

Benefits

- Reduces the need for multiple pieces of equipment
- Designed for independent use, or with assistance
- Shower door or curtain can be drawn when bridge is removed

Accessories

- Padded Reducer Ring 7" W x 8" L Hole **MODEL #4960**
- Padded Reducer Ring 5" W x 6" L Hole **MODEL #4965**
- Footrests - Additional set for tub base **MODEL #5125**
- Seat Lap Belt **MODEL #5135**
- 7", 9" or 11" bridge lengths available
- Adjustable Bridge - adds additional 13.5" to 17.75" to bridge length **MODEL #5100-P04**
- Adjustable Bridge - adds additional 17.5" to 21.75" to bridge length **MODEL #5100-P05**



¹ Shown with optional adjustable bridge

Elite™, Omni™, and Versa™ Rolling Chairs Shown with Accessories for Rolling Chairs and Transfer Systems

Conversion Kits

Growth kits are available to convert Elite and Omni bath chairs to a larger size. The kit includes a new cover and two positioning belts.

PRODUCT	MODEL #	CONVERTS TO	MODEL #	CONVERTS TO
Elite Conversion Kit	8139	8130 (Medium)	8149	8140 (Large)
Omni Conversion Kit	5139	5130 (Medium)	5149	5140 (Large)



Head Support*

Cushion designed for better head positioning.

MODEL #8875

Lateral Support*

Cushion designed for better torso positioning.

MODEL #8830

Detachable Armrests

Attach and remove easily. Four height positions.

MODEL #5150

Versa Swing-Away Armrests (Not Shown)

MODEL #5105

Rear Anti-Tippers

For added safety and stability. Prevent chair from tipping.

MODEL #5127

Pelvic Positioning Belt*

MODEL #8834



Positioning Belt*

Quick-release strap that can be threaded through loops in multiple positions. Two belts included with each chair.

MODEL #8835

Anterior Support Tray*

(Must have #5150 armrests to use)

MODEL #5151

Padded Abductor/SplashGuard* (Not Shown)

Attaches easily to bath chairs to improve hip alignment and stabilize seating posture. For Elite only.

MODEL #8818

Padded Reducer Ring* (Not Shown)

Available for use with MODEL #5040 and #5140.

MODEL #8818

Footrests

Foot Support.

MODEL #5125

Seat Angle Bracket*

Raises front seat height 2½", 3", or 3½" to assist with pelvic stability and postural control.

MODEL #5126


* Not available on the Versa

ROLLING CHAIR	ELITE			OMNI			VERSA
MODEL NUMBER	8625	8725	8925	5020	5030	5040	5110
Weight	33	34	35	35	36	37½	40
Weight Capacity (lbs.)	130	175	200	130	175	200	300
Dimension (W x D)	21 x 19½			21 x 19½			21 x 19½
Back Angle Adjustment	0°, 30°, 60°, 90°			0°, 30°, 60°, 90°			N/A
Calf Angle Adjustment	0°, -30°, -60°			0°, -30°, -60°			N/A
a. Seat Depth	10	13	15	12	14	17	18
b. Inside Seat Width	14	15	17	13	14	16	19
c. Back Height	22	27	31	21	26	30	15
d. Calf Length	10	13	15	9	13½	15½	N/A
e. Seat Hole Dimension (W x L)	N/A	N/A	N/A	6 x 7	7 x 8	8 x 9	9 x 14½
f. Seat Height Range	23½ - 27½			23½ - 30½			25½ - 29½
g. Rolling Base - inside width	17			17			17
h. Rolling Base - height clearance	17 - 21			17 - 21			17 - 21
i. Floor to bottom of bucket	N/A			10½ - 14½			10½ - 14½


Note: Dimensions are in inches unless otherwise noted




Dimensions and Specifications for Elite™, Omni™, and Versa™




	ELITE			OMNI		
COMPLETE SYSTEM MODEL NUMBER	8120	8130	8140	5120	5130	5140
SEATING MODEL NUMBER	8122	8132	8142	5122	5132	5142
Weight Capacity	130 lbs.	175 lbs.	200 lbs.	130 lbs.	175 lbs.	200 lbs.
Back Angle Adjustment	0°, 30°, 60°, 90°			0°, 30°, 60°, 90°		
Calf Angle Adjustment	0°, -30°, -60°			0°, -30°, -60°		
(a) Seat Depth	10	13	15	12	14	17
(b) Inside Seat Width	14	15	17	13	14	16
(c) Back Height	22	27	31	21	26	30
(d) Calf Length	10	13	15	9	13½	15½
(e) Seat Angle (Model #5126)	Up to 10°					



FOLDABLE TUB BASE TRANSFER FRAME - MODEL NUMBER 5113						
Chair + Rolling Base Weigh	31 lbs.	34 lbs.	36 lbs.	36 lbs.	38 lbs.	40 lbs.
Tub Base Weight	13½ lbs.			13½ lbs.		
Rolling Base Dimension	21½ W x 20¼ D			21½ W x 20¼ D		
Seat Height Range (1" increment)	23½ - 27½			26-30		
(g) Rolling Base - inside width	16			16		
(h) Rolling Base - height clearance	15½ - 19½			15½ - 19½		
(i) Tub Base Width at Top	31			31		
(j) Tub Base Footprint Depth Range	30½ - 33			30½ - 33		
(k) Tub Base Clearance Height Range*	13¾ - 20			13¾ - 20		
(l) Tub Base Footprint Width Range	28 -30			28 -30		
(m) Floor to bottom of Bucket	N/A			9¾ - 13¾		
Overall Length (Transfer System)	52			52		



COMPACT 18"X18" TUB BASE TRANSFER FRAME - MODEL NUMBER 5103						
Chair + Rolling Base Weight	41 lbs.	44 lbs.	46 lbs.	46 lbs.	48 lbs.	50 lbs.
Tub Base + Bridge Weight	20 lbs.			20 lbs.		
Rolling Base Dimension	19½ W x 17 D			19½ W x 17 D		
Seat Height Range (¾" increment)	23½ - 26¼			26 - 29		
(g) Rolling Base - inside width	17½			17½		
(h) Rolling Base - height clearance	17½, 18¼, 19, 19¾, 21½			17½, 18¼, 19, 19¾, 21½		
(i) Tub Base Width	18			18		
(j) Tub Base Footprint Depth	17			17		
(k) Tub Base Clearance Height Range	14 - 18			14 - 18		
(l) Tub Base Footprint Width	17½			17½		
(m) Floor to bottom of Bucket	N/A			12 - 16		
(n) Overall Length (Transfer System)	53½			53½		



VERSA TRANSFER SYSTEM - MODEL NUMBER 5100	
Weight Capacity	300 lbs.
Chair + Rolling Base Weight	40 lbs.
Tub Base + Bridge Weight	20 lbs.
Rolling Base Dimension	19½ W x 17 D
Tub Base Dimension	18 W x 17 D
Inside Seat Depth	18
Seat Hole	9 W x 14½ L
(a) Seat Width	19
(b) Back Width	16½
(c) Back Height	15
(d) Seat Height Adjustment Range (¾" increment)	24½, 25¼, 26, 26¾, 27½
(e) Overall length (Transfer System)	53
(f) Maximum Clearance Height (Toilet, Tub)	20

Note: Dimensions are in inches unless otherwise noted

* Custom heights available. Contact Customer Service for details.

PVC Hi-Back Shower/Commode Chair



MODEL #5060 shown with factory installed Footrest **MODEL #5000-P04**



MODEL #5075

VERSATILITY, CONVENIENCE, SAFETY, AND VALUE FOR LARGER USERS TOILETING AND BATHING NEEDS

Features

- Open front padded seat
- Tall backrest
- Adjustable height
- 4" Deluxe dual wheel casters, locking rear casters

Available Mesh Colors*



Blue

Benefits

- Durable and safe construction, easy to clean
- Comfortable and convenient to use
- Can use with support pads and belts for extra support
- Fits a range of user sizes
- Casters lock for safety and will not mark or rust

MODEL NUMBER	5060	5070	5075	5080
Product Weight	30 lbs.	31 lbs.	33 lbs.	47 lbs.
Weight Capacity	300 lbs.	350 lbs.	400 lbs.	525 lbs.
Overall Height Range	53 - 56½	53 - 56½	53 - 56½	53 - 56½
Seating Area Depth	18	19	19½	20
Back Height from Seat	33	33	33	33
Inside Width	17	20½	24	27
Outside Width	21	24½	28¼	30
Overall Depth	25	26¼	26¼	26¼
Floor to Top of Seat¹	20	20	20	20
Floor to Bottom Clearance	16¾	18	N/A	N/A
Seat Hole	12 W x 8¼ L	12 W x 8¼ L	12 W x 8¼ L	12 W x 8¼ L
ACCESSORY OPTIONS	5060	5070	5075²	5080³
7 qt. Commode Pail	5000-P01	5000-P01	N/A	N/A
10 qt. Commode Pail	5000-P02	5000-P02	N/A	N/A
Footrest – Factory Installed	5000-P04	5000-P04	5000-P04	5000-P05

Note: Dimensions are in inches unless otherwise noted

PVC Hi-Back Bariatric Shower/Commode Chair

SAFETY, COMFORT, AND VALUE IN TOILETING AND BATHING FOR BARIATRIC INDIVIDUALS



MODEL #5085 shown with optional Footrests **MODEL #5000-P06**

Features

- Extra wide, bariatric design supports up to 700 lbs.
- Open front padded seat
- Tall backrest
- Adjustable height
- 4" Deluxe dual wheel casters, locking rear casters

Benefits

- Safe and comfortable use for bariatric users
- Durable and safe construction, easy to clean
- Comfortable and convenient to use
- Can use with support pads and belts for extra support
- Fits a range of user sizes
- Casters lock for safety and will not mark or rust

Accessories

- Footrests **MODEL #5000-P06**

MODEL NUMBER	5085³
Product Weight	74 lbs.
Weight Capacity	700 lbs.
Overall Height	51 - 54½
Seating Area Depth	22½
Back Height from Seat	30
Inside Width	31½
Outside Width	35¼
Overall Depth	29½
Floor to Top of Seat¹	21
Floor to Bottom Clearance	N/A
Seat Hole	12 W x 8¼ L

Note: Dimensions are in inches unless otherwise noted

¹Based on lowest position

²10 qt. Commode pail included

³16 qt. Commode pail included

*Other colors available as special option. Contact Customer Service for details.

Available Mesh Colors*



Blue

PVC Hi-Back Foldable Shower/Commode Chair

INNOVATIVE FOLDING DESIGN OFFERS EASE OF USE AND STORAGE

Features

- Compact folding design
- Open front padded seat
- Tall backrest
- 4" Deluxe dual wheel casters, locking rear casters

Benefits

- Easy to store when not in use
- Durable and safe construction, easy to clean
- Easy, comfortable and convenient to use
- Can use with support pads and belts for extra support
- Casters lock for safety and will not mark or rust

Available Mesh Colors*



Blue

MODEL NUMBER	5015
Product Weight	24 lbs.
Weight Capacity	250 lbs.
Overall Height	54
Seating Area Depth	18
Back Height from Seat	31¾
Inside Width	17
Outside Width	21
Overall Depth	25¼
Floor to Top of Seat¹	22¼
Floor to Bottom Clearance	17½
Seat Hole	12 W x 8¼ L

Note: Dimensions are in inches unless otherwise noted



Chair folds for storage



PVC Self-Propelled Shower/Commode Chair

ALLOWS USERS MORE INDEPENDENCE DURING TOILETING AND BATHING

Features

- 24" plastic rear wheels with locking brakes
- Open front padded seat
- Anti-tip wheels in rear
- Fold-away footrest

Benefits

- Durable and safe construction, easy to clean
- Allows user to self-propel chair for independent mobility; will not rust or mark
- Comfortable and convenient to use
- Anti-tip wheels avoid possibility of tipping backwards
- Footrest lowers for comfort, folds up when not in use

Available Mesh Colors*



Blue

MODEL NUMBER	5050
Product Weight	37 lbs.
Weight Capacity	300 lbs.
Overall Height	38
Seating Area Depth	18
Back Height from Seat	17½
Inside Width	17
Outside Width	29
Overall Depth	36½
Floor to Top of Seat¹	20½
Floor to Bottom Clearance	17¼
Seat Hole	12 W x 8¼ L
Footrest	14 W x 6½ L
Footrest Height Above Floor	5½

Note: Dimensions are in inches unless otherwise noted

¹Based on lowest position



*Other colors available as special option. Contact Customer Service for details.

TheraGait™ PVC Gait Trainer

VERSATILE, ADJUSTABLE, AND ECONOMICAL GAIT TRAINER

Features

- Adjustable sling seat and padded armrests
- Rear safety gates
- Adjustable height; ergonomic angle handles
- 4" Deluxe dual wheel casters, locking rear casters

Benefits

- Durable construction, easy to clean
- Fit to individual use, easy to adjust
- Improves stability and increases safety
- Fits a range of user sizes; provides maximum comfort
- Casters lock for safety and will not mark or rust

Available Colors



Blue



Includes support harness



MODEL NUMBER	2760
Product Weight	41 lbs.
Maximum Weight Capacity	250 lbs.
Top Frame Inside Width	19½
Top Frame Outside Width	29½
Bottom Frame Inside Width	22
Bottom Frame Outside Width	26
Armrest Width	5
Armrest Height Range (2" Increment)	39¼ - 59¼
Overall Length	48
Top Frame Length (Includes pad)	32½

Note: Dimensions are in inches unless otherwise noted

TheraGait™ PVC Compact Gait Trainer

COMPACT SIZE FOR USE AT HOME OR IN FACILITIES WITH LIMITED SPACE

Features

- Compact size
- Adjustable height and sling seat
- Front safety gate
- 4" Deluxe dual wheel casters, locking rear casters

Benefits

- Durable construction, easy to clean
- Easy to maneuver in home and tight spaces
- Easy to adjust and fit to individual
- Improves stability and increases safety
- Casters lock for safety and will not mark or rust

Available Colors



Blue



MODEL NUMBER	2730
Product Weight	25 lbs.
Weight Capacity	200 lbs.
Overall Length	28½
Top Frame Inside Width	17
Top Frame Outside Width	21
Bottom Frame Inside Width	20½
Bottom Frame Outside Width	24½
Hand Rail Range	36½ - 39½
Seat Width	15½ W x 7 L
Overall Height to Pole/Swing Gate Entrance	43½

Note: Dimensions are in inches unless otherwise noted

ADJUSTABLE WALKER CHAIR OFFERING PERFORMANCE, STABILITY AND VALUE

Features

- Adjustable height
- Padded seats
- Front safety gate
- 4" Deluxe dual wheel locking rear casters
- Sliding footrest (**MODEL #2710** only)
- Wide base (**MODEL #2720** only)

Benefits

- Durable construction, easy to clean
- Custom fit to individual use, easy to adjust
- Offers a comfortable and stable sitting platform
- Improves stability and increases safety
- Casters lock for safety and will not mark or rust
- Comfort footrest allows for dual use as a transport chair (**MODEL 2710** only)
- Minimizes possibility of tipping; increases safety

Accessory

- Retractable Side Stabilizers

MODEL #2710-P01¹

Available Colors

 Blue

MODEL NUMBER	2710	2720
Product Weight	31 lbs.	34 lbs.
Weight Capacity	250 lbs.	250 lbs.
Hand Rail Range (1½" Increment)	40½ - 43½	38½ - 46½
Top Frame Inside Width	17	17¼
Top Frame Outside Width	21	21¼
Bottom Frame Inside Width	23	17¼
Bottom Frame Outside Width	27	31
Seating Area Depth	12	10
Overall Length	37	39¼
Seat Height Range	22¼ - 25¼	19 - 27
Overall Width	27 ¹	31

¹Overall width is 35" with optional Retractable Side Stabilizer MODEL #2710-P01
Note: Dimensions are in inches unless otherwise noted



MODEL #2710 shown with optional factory installed Retractable Side Stabilizers **MODEL #2710-P01**. Sliding footrest standard



MODEL #2720 with standard wide extended base

Rickshaw™ Rehab Exerciser

(aka PowerFit™)

**EXERCISE EQUIPMENT DESIGNED
FOR WHEELCHAIR USERS
TO STRENGTHEN ARM AND
SHOULDER MUSCLES NEEDED
FOR PROPULSION, TRANSFER, AND
PRESSURE REDUCTION LIFTS**




Weights not included

Features

- Innovative design with removable crosspiece allows for independent arm workout
- Fits around the width of any standard wheelchair
- Heavy-duty steel construction
- Latex and silicone free

Benefits

- Ideal for use in home or facilities
- Strengthens the arm and shoulder muscles needed for propulsion, transfer and pressure reduction lifts
- Can be positioned away from or facing the user



MODEL NUMBER	3000
Product Weight	37½ lbs.
Weight Capacity Per Side	100 lbs.
(a) Outside Width	36
(b) Inside Width	30
(c) Overall Length	50
(d) Arm Length	26
Handle Heights	28 and 32

Note: Dimensions are in inches unless otherwise noted

Columbia Medical “Blue” PVC Hub Story

U.S. Patents #3,795,218; D606,389S; D606,387S; D606,388S; D592,492

Why are the connecting hubs blue?

Connecting hubs are critical when building PVC devices. To improve durability and increase the strength in this key component, Columbia Medical formulated the proprietary rehab **Columbia Medical “Blue” PVC Hub**. Compared to conventional hubs, the **Columbia Medical “Blue” PVC Hub** provides greater strength and durability for use in clinical applications.

PVC Hubs and Tubes

Features

- Proprietary rehab grade Blue Connector
- Industry leading durability
- Rounded bent corners instead of joint connectors
- The only Blue and White medical supports in the rehab industry

Benefits

- Less breakage
- Stronger than standard PVC connectors and tubes for increased safety and durability
- Material designed specifically for use in clinical applications
- Easy to clean and keep sanitary

INDUSTRY LEADING DURABILITY



EZ Clean™ Casters*

Features

- Removable center axle for cleaning



ComfortMesh™ Material

Featured on most Columbia Medical bath products

Features

- Bacterial and fungal resistant
- Latex and silicone free

Benefits

- Stronger and softer than standard mesh covers to last longer and improve user comfort
- Easy to clean
- Does not contain common allergenic materials



*EZ Clean Casters is a trademark of Care Products, Inc.



Solutions for Daily Living

Optima™ TSS™ (Toilet & Shower System)

MULTI-USE SYSTEM FOR COMPLEX POSITIONING AND SUPPORT NEEDS



MULTIPLE ADJUSTMENTS



- Versatile system accommodates complex positioning needs
- Utilizes unique wrap-around back supports
- Offers maximum adjustable options
- System can grow with the individual
- Small footprint and tight turning radius fits in small spaces
- Multi-use system, functions as a stand-alone commode, over-the-toilet system, and roll-in shower chair, lowering overall equipment cost

QUALITY • INNOVATION • SAFETY

- ENHANCING THE QUALITY OF LIFE FOR INDIVIDUALS WITH DISABILITIES
- MEETING THE UNIQUE NEEDS OF OUR CLIENTS
- PROVIDING SAFE SOLUTIONS FOR DAILY LIVING

For over 30 years, Columbia Medical has established a reputation for developing high quality, innovative products that enhance the lives of individuals with disabilities. Our commitment to quality, innovation and safety are unsurpassed.

Car Seat Safety Meets Complex Positioning

Spirit™ APS™ car seat accommodates users up to 130 lbs.

FEATURING INNOVATIVE POSITIONING SYSTEMS



Spirit™ APS™ (Adjustable Positioning System™) Car Seat

is the only car seat to feature optional swing-away hip and trunk supports normally found in complex seating and wheelchair systems.

What's more, the Spirit provides:

- Adjustable padded headwings
- 5-point harness
- Optional swing-away abductor
- Optional seat depth extender
- Optional pressure reducing viscoelastic foam