

Instructions for Classic Whirlpool Disinfectant Cleaner

For broad spectrum disinfection of gram negative and gram positive bacteria, including *Pseudomonas aeruginosa* and *Staphylococcus aureus*, add two ounces Classic Whirlpool Disinfectant Cleaner to one gallon of water. Remove heavy soil or gross filth from the surface to be disinfected, then apply properly diluted solution with a mop, cloth, sponge, or hand pump trigger sprayer so as to wet the surface thoroughly. Allow to remain wet for 10 minutes and then let air dry. If higher detergency is desired, increase dilution to 4 to 6 ounces per gallon of water. Prepare a fresh solution for each use or when solution becomes visibly dirty. To disinfect hard, non-porous, inanimate surfaces (such as fiberglass and stainless steel tubs and chair surfaces, chrome plated intakes and lifts, etc.), apply properly diluted Classic Whirlpool Disinfectant Cleaner so as to wet all surfaces thoroughly. For routine disinfection, proper dilution is 1:64 (2 ounces of product per gallon of water). Allow to remain wet for 10 minutes, then let air dry. For heavily soiled areas, a pre-cleaning step is required. Prepare a fresh solution of reach use as above.

FOR CLEANING BATH AND THERAPY EQUIPMENT: After using the whirlpool unit, drain the water and refill with fresh water to just cover the intake valve. Add 2 ounces of Classic Whirlpool Disinfectant Cleaner for each gallon of water in the unit at this point. Briefly start the pump to circulate the solutions. Turn off pump. Wash down the unit sides, seat of the chairlift and any/all related equipment with a clean swab or sponge. After the unit has been thoroughly cleaned, drain solutions from the unit and rinse any/all clean surfaces with fresh water. The unit is ready for re-use.

SANITIZING: To sanitize non-porous, non-food contact surfaces, add 2 ounces of Classic Whirlpool Disinfectant Cleaner per gallon of water. Remove heavy soil or gross filth from the surface to be disinfected, then apply properly diluted solution with a mop, cloth, sponge, or hand pump trigger sprayer so as to wet the surface thoroughly. Allow to remain wet for 30 seconds and let air dry.