

### 3-Position Recliner, Bariatric

HCPCS: E1031



FIGURE A



FIGURE B

- 3 positions: upright (Figure A), deep recline and elevated leg rest/footrest (Figure B)
- Retractable lock bar secures the chair in the desired position
- Self-adjusting suspension in the seat prevents bottoming out
- Side panels "pop off" for easy cleaning and maintenance
- Large, blow-molded tray locks in 5 positions and can be stored on the side of the chair
- 5" casters (2 with locks) make it easy to move and secure the recliner
- Wide, comfortable armrests are angled to provide support in any reclined position
- Neutral-colored, sturdy steel frame is easy to clean and maintain
- Moisture barrier on seat prevents seepage
- Comfortable built-in headrest
- Foam and vinyl meet California Technical Bulletin 117

#### PRODUCTS

#### Description

D574EW-BR	3 Position Recliner,ExtraWide,1/cs
D574EW-J	3 Position Recliner,ExtraWide,Jade,1/cs
D574EW-R	3 Position Recliner,ExtraWide,1/cs

#### SPECIFICATIONS

#### Description

Arm Height From Floor	26.5"
Arm Height From Seat	9.5"
Carton Shipping Weight	127 lbs.
Overall Dimensions	47.5" (H) x 32" (W) x 41.5" (D)
Seat Dimensions	21" (H) x 25" (W) x 20" (D)
Warranty	Limited Lifetime on Frame, 1 Year on Upholstery and Moving Parts
Weight Capacity	500 lbs
Width Between Arms	24.5"

#### PARTS

#### Description

D574P-1064	Slide Tube,Table Tray,D574,1/pair, EA
D574P-1036	Table Fitting,D574,1/ea, EA
D574P-1027	Front Casters W/Hardware, D574, 1pr/bx, BX
D574P-1020	Cut Side Panel,D574,Left&Rt,1pr, BX
D574P-1025B	Rear Locking Caster,D574,1pr/bx, BX
D574P-1037	Washer & Screw,D574tray,1/ea, EA
D574P-1061	Bolt Set,Footrest Sup.Bar,1/ea, EA
D574EW-1003	Push Bar for D574EW,1/ea, EA

PARTS	Description
D574EWBACK-BR	Back,D574EW-BRBl.Ridge,1/bx, BX
D574EWBACK-J	Back,D574EW-Jade,1/ea, EA
D574EWSEAT-BR	Seat W/Legrest,D574EW-BL,Bl.Ridge,1/bx, BX
D574EWSEAT-J	Seat W/Legrest, D574EW-Jade, 1/ea, EA
D574EWSEAT-R	Seat W/Legrest,D574EW-R,Rosewood,1/bx, BX
D574EW-1010BR	Footrest,D574EW,BlueRidge,1/bx, EA
D574EW-1011	Footrest Support Bar,D574EW,1/ea, EA
D574EW-1035	Recliner Mech.D574EW,RT&LFT,1/ea, EA
D577P-1012	Protecting Panel,D577,1/ea, EA
D574EW-1065	Tray Table Assm.D574EW,1/cs, EA
D574EW-1017N	Frame, D574EW, New, 1/ea, EA
D574EW-1070BR	Armrest,D574EW,BlueRidge,1/ea, EA