

## Clinical Care Recliner



HCPCS:



FIGURE A



FIGURE B



FIGURE C

- Ideal for use in dialysis, blood collection, respiratory and cardiac care
- Chair can be positioned in 4 positions upright (Figure A), TV (Figure B), Full Recline (Figure C) and Trendelenburg
- Deep recline and Trendelenburg position can only be attained by an attendant
- Gas cylinder controls deep recline and Trendelenburg position
- Side panels "pop-off" for easy cleaning and maintenance
- Large, blow-molded side trays includes recess for cup
- 5" Casters (2 with locks) make it easy-to-move and secure recliner
- Gas cylinder release handle for quick back release
- Release handle can be operated from either side of the recliner
- Wide, comfortable armrests are angled to provide support in any reclined position
- Upholstered arm supports provide extra comfort in recline position
- Neutral-colored, sturdy steel frame is easy to clean and maintain
- Comfortable built-in headrest
- Foam and vinyl meet California Technical Bulletin 117 Standard

### PRODUCTS

#### Description

D577-BR	Clinical Care Recliner, Blue Ridge
D577-J	Clinical Care Recliner, Jade, 1/cs
D577-R	Clinical Care Recliner, Rosewood, 1/cs

### SPECIFICATIONS

#### Description

Carton Shipping Weight	135 lbs.
Floor to Seat Height	26.5"
Overall Dimensions	47.5" (H) x 30.5" (W) x 38" (D)
Product Weight	115 lbs.
Seat Dimensions	21.5" (H) X 21.5" (W) X 20" (D)
Warranty	Limited Lifetime on Frame, 1 Year on Upholstery and Moving Parts
Weight Capacity	250 lbs
Width Between Arms	22"

### PARTS

#### Description

D574P-1048-R	HeadrestAdjBracket,D574/D574EW,Rght,1ea, EA
D577P-1004BR	Side Wing Pad,D577 BlueRidge,1ea, EA
D577P-1004J	Side Wing Pad,D577 Jade,1ea, EA
D577P-1004R	Side Wing Pad,D577 Rosewood,1ea, EA
D577P-1001	Recliner Cable,D577,1/ea, EA
D577P-1005	Recliner Release Handle,D577,1/ea, EA

PARTS	Description
D577P-1060	Back Push Bar,D577,1/ea, EA
D577BACK-J	BackOnly,D577-J,Jade,1/ea, EA
D577P-1010J	Footrest,D577,Jade,1/bx, EA
D577P-1075	Formed Side Panels,Lft& Rt,1/bx, BX
D577P-1026	Front Casters,D577,1pr/bx, BX
D577P-1018BR	Armrest,BlueRidge,D577,1/ea, EA
D577P-1018J	Armrest,Jade,D577,1/ea, EA
D577P-1065	Side Table Assembly,D577,1/ea, EA
D577P-1022	Reclining Piston Release Mount,D577,1ea, EA
D577P-1050	Reclining Piston ,D577,1ea, EA
D577P-1058	Hardware Pack,D577,1set/bx, BX
D577P-1012	Protecting Panel,D577,1/ea, EA
D574P-1048	HeadrestAdjBracket,D574/D574EW,Left,1ea, EA
D574P-1035	Recliner Mechanism,D574,Rt.&Lf,1/bx, BX