

Deluxe, Tilt-Top Overbed Table

HCPCS: E0274



FIGURE A



FIGURE B

- Easily assembled
- Walnut, wood grain low-pressure laminate top
- Top can be tilted 33° in either direction (Figure A)
- Top can be raised or lowered in infinite settings between 29.5" - 46"
- Swivel casters allow for easy maneuverability
- Chrome-plated steel "H" base provides security and stability
- Tabletop locks securely when height adjustment handle is released It can be raised with the slightest upward pressure
- For proper assembly, do not remove bottom screw (Figure B)

PRODUCTS

Description

13008	Over Bed Table,Tilt,1/cs
-------	--------------------------

SPECIFICATIONS

Description

Base Dimensions	15.5" (D) x 26.5" (W)
Carton Shipping Weight	23 lbs.
Height Adjustable	29.5" - 46"
Tabletop Surface Dimensions	15" (D) x 30" (W)
Warranty	Limited Lifetime
Weight	21 lbs.
Weight Capacity	50 lbs

PARTS

Description

13008CP	Caster Plug, 13008,4/bx
13008H	HardwareSet/13008,1/ea
13007CL	Caster for 13007 w/Lock,1/ea
13008AM	Height Adjust.Mast,13008,1/ea
13007C	Caster,Overbed table.1/ea
13008B	Bushing,Plastic,for 13008,1/ea
13008K	Angle Adj.Knob,13008/1ea
13008AB	Bushing, AngleBolt,13008,1/ea
13008EL	Elevating Lever, 13008,1/ea
13008AAB	Adjustment Bolt,13008,1/ea
13008-T	Replacement Top,13008,1/ea
13008TB	Table Top Bracket Only, 13008, 1/ea