



Enhance Your Life

TRIAXE Tour

<ul style="list-style-type: none"> • Lightweight Aluminum Construction
<ul style="list-style-type: none"> • Rear Wheel Drive
<ul style="list-style-type: none"> • Easy Paddle Acceleration
<ul style="list-style-type: none"> • Durable - 375lb Weight Capacity
<ul style="list-style-type: none"> • Electromagnetic Brakes
<ul style="list-style-type: none"> • Folds Flat for Transport & Storage



Model	TRIAXE Tour T3050
Overall Length	38"
Overall Width	19"
Tires	8" Solid Front Tire 9.5 Pneumatic Rear Tires with Anti-Flat
Maximum Speed	6mph
Range	Up to 18 miles
Turning Radius	27"
Weight (without battery pack)	58lbs
Color	Fiery Red, Reflex Blue
Brake	Electric Brake Systems
Drive System	Rear Wheel Drive, Brushless Motor
Handlebar	Aluminum with Rubber Hand Grips
Maximum Weight Capacity	Up to 375lbs
Seat Width	13" standard, 18" optional
Seat Height	23"
Battery	24V 12AH
Motor	24v 180w
Battery Charger	24V 2AH