

Comfort - Functionality - Durability

Features:

- The Lumex FR577RG Recliner features three positions for patient comfort. The fourth position is a caregiver-activated Trendelenburg (or Shock) Position
- Improved wider and softer armrest for patient comfort
- Top section of headrest can be moved forward for enhanced patient comfort
- Wide range of fabrics / colors
- Quick Ship available
- Easy-to-activate gas cylinder mechanism for quick Trendelenburg activation
- Durable, high-impact side tables sit flush with recliner arm to conveniently hold prep materials. Side table surface is smooth and features tapered edges for easy cleaning

- Easily removed side panels enable access to recliner for cleaning and maintenance
- Four Tente® swivel casters facilitate transport and maneuverability in confined areas. Three total locking casters with one directional locking caster
- 290 lb maximum weight capacity, evenly distributed

Options:

- Heat and Massage
- Fully Upholstered Armrests
- IV Pole Holder / IV Pole
- Headrest Cover
- Universal Pillow

Dimensions

Overall Height.....	47" / 119 cm	Seat Depth.....	21" / 53 cm
Overall Width.....	31" / 79 cm	Depth Fully Reclined	
Overall Depth.....	35" / 89 cm	In Trendelenburg	74" / 188 cm
Width Between Arms	22" / 56 cm	Width With Both Side Tables	
Arm Height From Seat.....	8" / 20 cm	in Upright Position	51" / 130 cm
Seat Height.....	22" / 56 cm	Maximum Weight Capacity	
Seat Width.....	21.5" / 55 cm	Evenly Distributed.....	290 lb / 132 kg



Trendelenburg Lever
Located on Push Handle



Articulating Headrest
in Forward Position



Smooth Tables with
No-Seam Tapered
Edge



Multiple Positions for Enhanced Patient Comfort