



Multi-Positioning Seat (MPS) 5 year warranty

Versatile Seating with Built-in Growth

NEW Cushion sizes available for both Small and Large MPS!
Cushions grow the seating system from 9 months to young adult!

- Adjustable 5-point harness
- Hip-Flex® Seat Cushion
- Height adjustable Trunk and Head Cushions



Use the Special Tomato® MPS system everywhere!



At home, secures to most chairs by using attachment straps



At school, secure the Special Tomato® in a classroom chair, or use it with the Mobile Base as your all-in-one classroom chair



In cars, use the crash-tested seat with LATCH and Tether system



On the school bus, transport the Special Tomato® in Push Chair (79200000, see next page)



In an airplane, both sizes approved by FAA for air travel for seats purchased after April 1, 2006



Features & Accessories



1 Mobile Floor Base

79070100 Small
79070300 Large

2 Headrest with Laterals

MHL shown in Lilac.
See page 7.

3 Standard Headrest

MHS shown in Cherry.
See page 7.

4 Optional Cushions

MS5 shown in Aqua.
Four sizes of Seat Cushions
and two Back Cushions for
even more growth!

5 MPS Tray

79060159 Dark Gray, Small
79060359 Dark Gray, Large
Depth adjustable. 60° tilt adjustment.
Tray easily attaches to back of Mobile
Floor Base for storage.

6 Push Chair Sun Shade

79210000
Provides protection from the sun and
a convenient window for monitoring.

7 MPS Push Chair Base

79200000 Umbrella style folding
stroller for easy transport. Folds
down to 15"W × 11"H × 48"L. Supports
150 lbs. and weighs only 20 lbs. 30°
fixed seat tilt. Designed for use with
standard tie-down systems for bus or
van transport and meets ANSI/
RESNA WC19 regulations.

8 Push Chair Travel Bag Accessory

79220000

Customize Your Size



Small Shell (MS)

Shells

Choose a Shell size to determine your options



Large Shell (ML)

Options

Headrests Cushions for both Small (MS) and Large Shells (ML)

MHS (Standard)



7.5"W x 6.5"H

MHL (w/ Laterals)




7"W x 6"H

MBH (Basic)



7.5"W x 6.5"H
1.5"D


MBL (Basic w/ Laterals)



7"W x 6"H
2.5"D

Back Cushions for Small Shell (MS)

MS2



8"W x 8.5"H

MS4



9"W x 11"H

ML2



10"W x 11"H

Back Cushions for Large Shell (ML)

ML2



10"W x 11"H

ML4



11"W x 13"H

ML6



13"W x 13"H

Seat Cushions for Small Shell (MS)

MS1



7"W x 6"D

MS3



9"W x 8"D

MS5



10"W x 10"D


MS7



10"W x 12"D

Seat Cushions for Large Shell (ML)

ML1



11"W x 11"D

ML3



12"W x 13"D

ML5



13"W x 15"D

ML7



13"W x 17"D

Customize Your Size

Fill in your custom sizes here for ease when ordering

MS

Example

MS

Shell

MHS

Head

MS2

Back

MS1

Seat

	Small MPS Seat Cushions				Large MPS Seat Cushions			
MPS Seat Specifications	MS1	MS3	MS5	MS7	ML1	ML3	ML5	ML7
Target User Age	9 mos - 3 yrs	2 - 5 yrs	4 - 8 yrs	5 - 9 yrs	6 - 9 yrs	8 - 11 yrs	10 - 14 yrs	12 yrs - Adult
Target User Height	28" - 37"	32" - 43"	40" - 50"	43" - 52"	45" - 58"	48" - 60"	48" - 62"	58" - Adult
Maximum User Wt.	80 lbs.				130 lbs.			
Seat Width	7"	9"	10"	10"	11"	12"	13"	13"
Seat Depth	6"	8"	10"	12"	11"	13"	15"	17"
Seat to Footrest	6" - 16"				8" - 18"			
Hip Flex Angle	-15° to +15°				-15° to +15°			
Seat Tilt Angle Forward	5°				5°			
Seat Tilt Angle Back	-30°				-30°			
Floor to Seat Ht.	14" - 28.5"				17" - 28.5"			

	Small MPS Back Cushions		Large MPS Back Cushions		
MPS Back Specifications	MS2	MS4	ML2	ML4	ML6
Target User Age	9 mos - 3 yrs	2 - 9 yrs	6 - 9 yrs	8 yrs - Adult	10yrs - Adult
Target User Height	28" - 37"	32" - 52"	45" - 58"	48" - Adult	48" - Adult
Maximum User Wt.	80 lbs.		130 lbs.		
Back Width	8"	9"	10"	11"	13"
Back Height	8.5"	11"	11"	13"	13"
Top of Back Cushion	11" - 16"	15" - 20"	15" - 20"	19" - 22"	19" - 22"
Top of Head Cushion	18" - 23"	21" - 26"	24" - 28"	26" - 30"	26" - 30"

MPS Tray	Small (79060159)	Large (79060359)
Dimensions	24"W x 17"D x 2"H	
Weight	8 lbs.	

MPS Hi-Low Seating System

79300109 Dark Gray, Small
79300309 Dark Gray, Large
Includes Seat Shell with
Footrest, 8-point Harness,
Hi-Low Base with Casters
(Must Choose Custom Seating
Cushions Separately)



Seat Cushions	Cherry	Lilac	Lime	Aqua	Dark Gray	Light Gray
MS1	79050551	79050552	79050553	79050554	79050559	79050560
MS3	79050151	79050152	79050153	79050154	79050159	79050160
MS5	79050251	79050252	79050253	79050254	79050259	79050260
MS7	79050651	79050652	79050653	79050654	79050659	79050660
ML1	79050751	79050752	79050753	79050754	79050759	79050760
ML3	79050351	79050352	79050353	79050354	79050359	79050360
ML5	79050451	79050452	79050453	79050454	79050459	79050460
ML7	79050851	79050852	79050853	79050854	79050859	79050860

Back Cushions	Cherry	Lilac	Lime	Aqua	Dark Gray	Light Gray
MS2	79120551	79120552	79120553	79120554	79120559	79120560
MS4	79120151	79120152	79120153	79120154	79120159	79120160
ML2*	79120751	79120752	79120753	79120754	79120759	79120760
ML4	79120351	79120352	79120353	79120354	79120359	79120360
ML6	79120451	79120452	79120453	79120454	79120459	79120760

Headrests	Cherry	Lilac	Lime	Aqua	Dark Gray	Light Gray
MHS (Standard)	79100051	79100052	79100053	79100054	79100059	79100060
MHL (w/Laterals)	79110051	79110052	79110053	79110054	79110059	79110060
MBH (Basic)	79130051	79130052	79130053	79130054	79130059	79130060
MBL (Basic w/ Laterals)	79140051	79140052	79140053	79140054	79140059	79140060

*ML2 designed to fit both the Small and Large system, providing a wider Back Cushion option for Small size systems.
Warranty: 5-Year Manufacturer's Warranty on cushions, and seat shell; 2-Year Manufacturer's Warranty on Hi-Low base, and strapping.