



MULTICHAIR 5200

The MULTICHAIR 5200 is designed for independent wheelchair users who cannot access the tub because the toilet is in the way, or for able-bodied walkers not comfortable rotating their hips as required to use a traditional stationary transfer bench. The connector bridge is removable for shower curtain closure or improved caregiver access, and the shuttle rides on tough, corrosion-proof Delrin® wheels.

Features

- One tub section, one removable connector bridge and one commode section
- Waterproof nylon cloth back
- Cantilever fold-back / locking padded chair arms
- Height adjustable cantilever fold-back / locking padded arms
- Seamless padded self-skinning urethane back cushion and 20" x 20" x 2" thick naugahyde seat cushion with 1/2" waterfall-style foam ring around the commode opening and slot help protect the ischial and coccyx areas to reduce risk of pressure sores
- Seat height adjusts from 20-1/2" to 24-1/2" (easily fits over the highest toilets required by the ADA)
- Won't rust or corrode—constructed of aircraft-grade aluminum, brass, stainless steel and durable plastics

Specifications

- Approximate weight: 28 pounds
- Weight capacity: 300 pounds
- 17" overall width at the outside of leg crutch tips for each section
- Passed ANSI / BIFMA X5.1 test performed by Gaynes Labs, Inc. (a national testing laboratory). Weights of 250 lbs. and 300 lbs. were dropped on the chair, with no structural breakage or loss of serviceability

Dimensions and Notes

- A = 35-1/2" to 39-1/2" (adjustable)
- B = 15-1/2" (top edge of seat cushion to top of back cushion)
- C = 16-1/2" (inside minimum clearance)
- D = 20-1/2" to 24-1/2" (adjustable, floor to top of seat cushion)

