

INNOVATIVE PRODUCTS FOR SENIORS AND PEOPLE WITH DISABILITIES



MULTICHAIR 6000Tx
(carrying case included)



A truly cost-effect alternative to installing a roll-in shower, the 6000Tx eliminates risky bathroom transfers. The seat adjusts from 20" to 24" high for access over the highest toilets required by the ADA, and the exclusive self-aligning connector bridge easily removes for shower curtain closure or improved caregiver access. Light, strong, and portable - it travels easily in its standard-equipment carrying case.

Features

- Swing-away / removable push handles
- Adjustable waterproof nylon cloth back
- Cantilever fold-back / locking padded chair arms
- Self-skinning padded urethane back and seat cushion with 1/2" waterfall-style foam ring around the commode opening and slot help protect the ischial and coccyx areas to reduce r risk of pressure sores
- Seat height adjusts from 20-1/2" to 24-1/2" (easily fits over the highest toilets required by the ADA)
- Adjustable swing-away footrests
- Four 5" casters with total-lock brakes
- Different length bridges available to fit nearly any bathroom.
- Black nylon carrying case (29" x 20" x 10"-13.5" expandable)
- Won't rust or corrode—constructed of aircraft-grade aluminum, brass, stainless steel and durable plastics

Specifications

- Approximate weight: 42 lbs.
- Weight capacity: 300 lbs.
- Passed ANSI / BIFMA X5.1 test performed by Gaynes Labs, Inc. (a national testing laboratory). Weights of 250 lbs. and 300 lbs. were dropped on the chair, with no structural break age or loss of serviceability

