



# Jazzy<sup>®</sup>

## SELECT ELITE

The Jazzy Select<sup>®</sup> Elite delivers a reliable blend of power, performance and style. In-line, front-wheel drive technology gives the Elite excellent stability and manoeuvrability for solid performance indoors and out.



**Pride<sup>®</sup>**  
Mobility Products  
Australia Pty. Ltd.  
— LIVE YOUR BEST —

## SPECIFICATIONS

<b>Weight capacity</b>	Up to 136 kg (300 lbs.)
<b>Maximum speed<sup>1</sup></b>	Up to 6.44 km/h (4 mph)
<b>Ground clearance</b>	3.81 cm (1.5") (at anti-tips and motors)
<b>Turning radius<sup>2</sup></b>	61 cm (24")
<b>Overall length<sup>2,6</sup></b>	87.3 cm (37.375") without front riggings
<b>Overall width<sup>2</sup></b>	58.4 cm (23")
<b>Drive wheels</b>	22.86 cm (9") solid
<b>Anti-tips</b>	7.62 cm (3") solid, front-mounted
<b>Casters</b>	15.24 cm (6") solid, rear-mounted
<b>Suspension</b>	Limited
<b>Drivetrain</b>	Two-motor, in-line, front-wheel drive
<b>Braking system</b>	Intelligent braking (electronic, regenerative disc brakes)
<b>Standard electronics</b>	40 amp, PG GC 3 controller
<b>Standard seat</b>	Comfort high-back
<b>Specialty controls</b>	N/A
<b>Battery charger</b>	Off-board, 3 amp
<b>Per-charge range<sup>1,4</sup></b>	Up to 21.97 km (13.65 miles)
<b>Battery requirements</b>	(2) 12 volt deep-cycle
<b>Battery size (largest)</b>	(U-1) 30-34 AH
<b>Battery weight<sup>3</sup></b>	11.11 kg (24.5 lbs.) each
<b>Base weight</b>	51 kg (112.5 lbs.)
<b>Standard seat weight</b>	17.7 kg (39 lbs.) (Comfort high-back)
<b>Max. seating dimensions</b>	50.8 cm x 50.8 cm (20" x 20")
<b>Warranty</b>	2-year limited

The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.

- 1) Range and speed vary with user weight, terrain type, battery charge, battery condition and tyre pressure. This specification can be subject to a variance of (+ 10% or - 5%) .
- 2) Due to manufacturing tolerances and continued product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 3) Battery weight may vary based on manufacturer.
- 4) Tested in accordance with ANSI/RESNA, WC Vol. 2, Section 4, standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Testing conducted at maximum weight capacities. This specification can be subject to a variance of +10%, -5%.
- 5) AGM or gel-cell type required
- 6) Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

## FEATURES

- In-line, front-wheel drive technology
- 40 amp, PG GC 3 controller
- Easy access to batteries from rear
- Simple main frame design for easy serviceability
- Comfort high-back seat with headrest

## OPTIONS

- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder

## COLOURS



Jazzy Red



Jazzy Blue