

Model Number	3W: SC610 4W: SC710
Weight Capacity	400 lbs.
Maximum Speed ¹	Up to 5.2 mph
Ground Clearance ³	3W: 2.2" at motor, 1.26" at anti-tip 3.4" at side deck 4W: 2" at motor, 1.5" at anti-tip 3.1" at side deck
Turning Radius ³	3W: 42" 4W: 54"
Overall Length ³	3W: 45.25" 4W: 47"
Overall Width ³	22.25"
Seat-to-Ground Height	20.5"-22"
Seat-to-Deck Height	17.5"-19"
Front Tires	3.5" x 10" solid
Rear Tires	3.5" x 10" solid
Suspension	N/A
Range Per Charge ^{1,2}	Up to 16 miles
Total Weight Without Batteries ⁵	3W: 133 lbs. 4W: 137 lbs.
Heaviest Piece When Disassembled ³	3W: 54 lbs. (front section) 4W: 61 lbs. (front section)
Standard Seating	Type: Compact, foldable Weight: 28 lbs. Material: Black Vinyl Dimensions: 18" width, 17" depth
Drive System	Rear-wheel drive sealed transaxle with 24-volt DC motor
Dual Braking System	Regenerative and electromechanical
Battery Requirements ^{4,6}	Size: (2) U1 standard, Weight: 24.5 lbs. each Size: (2) 40 AH (optional), Weight: 32.5 lbs. each
Battery Charger	3.5-amp, off-board (U1), 5-amp, off-board (40 AH)
Warranty	Frame: lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Battery: 6-months

HCPCS Code: K0801 (Group 1, Heavy Duty)

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- 2 Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Includes standard seat weight.
- 6 AGM or gel cell type required.

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

(US) 800-800-8586 • (Canada) 888-570-1113

pridemobility.com |    

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

VIC10_061219

VICTORY[®] 10

**ENHANCE
YOUR ACTIVE
LIFESTYLE
WITH ADDED
COMFORT AND
PERFORMANCE**



LIVE YOUR BEST[®]

pridemobility.com |    

STYLE, COMFORT, PERFORMANCE

18" x 17"
Pinchless hinge,
molded plastic,
black vinyl seat
with 8-point
swivel and
sliders for
increased
comfort



LED battery meter

Standard
front basket

Easy-to-grip
tiller adjustment
knob

Long-lasting,
maximum-
intensity LED
headlight
offers optimal
pathway
illumination

Exclusive stylish, lightweight, non-scuffing, black, low-profile tires



Tiller console



Feather-touch disassembly

The Victory® 10 series offers added style, comfort and performance for your active lifestyle. Standard features include LED pathway light, LED battery meter, feather touch disassembly, wraparound delta tiller, solid tires and a front basket. A 400 lb. weight capacity, a per charge range of up to 16 miles and a maximum speed up to 5.2 mph compliment the versatility of the Victory 10. Available in Viper Blue and Candy Apple Red.

- Exclusive one-hand feather touch disassembly
- Charger port located on the tiller lets you conveniently charge your scooter
- Delta tiller with ergonomic wraparound handles lets you operate the scooter with one hand and rest your wrist
- Easy-to-read LED battery meter
- Rearview mirror
- Long-lasting, maximum-intensity LED headlight offers optimal pathway illumination

Available Colors



Viper
Blue



Candy
Apple Red

Optional Seating

Limited recline, black vinyl, high back with height adjustable headrest, flip-up height and width adjustable armrests

- 18" x 18" (sliders)
- 20" x 20" • 22" x 18"

Deluxe contour high back black vinyl with height adjustable headrest, flip-up height and width adjustable armrests (no sliders)

- 18" x 18" • 18" x 20"
- 20" x 20" • 22" x 20"

Accessories

- Cane clip attachment
- Cane/crutch holder (single, double)
- Cell phone holder
- Cup holder
- Forearm crutch holder
- Oxygen tank holder
- Quad cane holder
- Rear basket
- Additional rearview mirror
- Saddlebag
- Safety flag
- Walker holder
- Weather cover
- XLR USB charger
- Lap belts: 50", 60" and 70"

Military Patches

