

# CARE CLINER SERIES

---

6530 (STANDARD W/ NYLON CASTERS) / 6531 (STANDARD W/ STEEL CASTERS)  
6540 (X-LARGE W/ NYLON CASTERS) / 6541 (X-LARGE W/ STEEL CASTERS)



## KEY FEATURES

---

- Standard headrest cover
- Padded armrests
- Dual fold-down side tables
- Removable side panels
- Trendelenburg positioning

## OUTSTANDING VALUE AND QUALITY

Lower price does not mean lower quality when you purchase a Winco Care Cliner.

Our Care Cliner comes with all the superior features you might find on higher-priced clinical chairs

ENHANCING THE PATIENT & CAREGIVER EXPERIENCE



# ANATOMY OF A CARE CLINER



## SPECIFICATIONS

MODELS	6530 / 6531	6540 / 6541
DESCRIPTIONS	STANDARD SIZE	XL SIZE
Seat Height	21" (53.3 cm)	21" (53.3 cm)
Seat Width	21" (53.3 cm)	25" (63.5 cm)
Seat Depth	20" (51 cm)	20" (51 cm)
Overall Height	49" (124.5 cm)	49" (124.5 cm)
Overall Width*	32.25" (81.9 cm)	36.25" (90.1 cm)
Back Height above Seat	31.5" (80 cm)	31.5" (80 cm)
Overall Length (Upright)	39" (99 cm)	44" (111.8 cm)
Overall Length (Reclined)	76.5" (194.3 cm)	76.5" (194.3 cm)
Armrest Height	27" (68.6 cm)	27" (68.6 cm)
Weight	94 lbs (42.6 kg)	111 lbs (50.4 kg)
Weight Capacity	275 lbs (124.7 kg)	2450 lbs (204 kg)

\* Add 9.25" to Overall Chair Width for each raised side table.

## OPTIONS



TB-133 (TB)



Leg Extensions (RE)



Side Cushions (SC)



TB-133 (TB)



Heat Only (H)



Heat and Massage (HM)

**NOTE** Winco reserves the right to change specifications without notice.

**WARRANTY:** Contact Winco or visit our website for warranty details.

**VINYL:** Available in over 50 vinyl colors and with optional TB133 protection. Visit our website or contact us for swatches.



## CONTACT US

To arrange a demonstration or to speak with one of our associates, please feel free to contact us at:

**800-237-3377 | WincoMfg.com | 5516 SW 1st Lane, Ocala, FL 34474 | Info@WincoMfg.com**