

Caremor Cliner 5351/5361



Our most comfortable and versatile recliner

Ease of operation, outstanding durability, and exceptional value have made this one of our bestselling recliners. The Caremor Cliner provides patient comfort and support as well as user friendly features.



5351 w/ tray
5361 w/out tray

KEY FEATURES

- Infinite positioning back
- Padded armrests w/ arm guard
- Webbed elastic suspension back and seat
- Removable side panels
- Continuous seat & legrest cushion

Proudly Made in the USA



1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

www.wincomfg.com



PRODUCT SPECIFICATIONS

CAREMOR CLINER 5351/5361

- 1 Push handle and back adjustment handle on back for infinite positioning.
- 2 Adjustment handle can be installed on side of chair for patient use.
- 3 Comfortable, padded armrests with guards to keep upholstery free from damage.
- 4 Recessed hand-grips assist in entering and exiting the chair.
- 5 Removable side panels for easy cleaning and maintenance.
- 6 5" Total lock casters (in front) allow for better handling and control.



5351/5361

Seat Height	.20.5" (52cm)
Seat Width	.20" (50.8cm)
Seat Depth	.18" (45.7cm)
Overall Height	.47" (119.4cm)
Overall Width*	.26.25" (66.7cm)*
Back Height Above Seat	.31.5" (80cm)
Overall Length (Upright)	.37" (94cm)
Overall Length (Reclined)	.62" (157.5cm)
Arm Height	.27" (68.6cm)
Weight	5351: - 81 lbs. (36.7kg) 5361: - 76 lbs. (34.5kg)
Weight Capacity	.300 lbs. (136kg)


*ADD 1.75" TO WIDTH WHEN TRAY IS IN STORED POSITION

SPECIFICATIONS

OPTIONS



STANDARD COLORS



Printed colors shown may vary slightly from actual color of product.

OPTIONAL COLOR UPGRADES

(ADDITIONAL CHARGES APPLY)



SILVERTEX VINYL

(ADDITIONAL CHARGES APPLY)



COLORS